CA-MI

HI-4 KIT

REF RE 300350

HI-4

REF RE 300350/01

IT Manuale d'uso

EN Instruction manual

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DE Handbuch

ES Manual de instrucciones

**COMPOSIZIONE DEL DISPOSITIVO/ COMPOSITION OF THE DEVICE
 COMPOSITION DU DISPOSITIF / ZUSAMMENSETZUNG DES GERÄTS /
 COMPOSICIÓN DEL DISPOSITIVO**

HI-4 KIT **REF** RE 300350

1
 Forcella nasale
 Nosepiece
 Embout nasal
 Nasenstück
 Horquilla nasal

2
 Boccheruola
 Mouthpiece
 Embout buccal
 Mundstück
 Boquilla

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HI-4 KIT is an Medical Device (Risk Class IIa - MMD 93/42/EEC and subsequent changes) highly efficient in medical therapies administered by means of aerosol nebulization. This device is designed for the nebulization of for atomising antibiotics and bronchodilator drugs (distribution of drugs with aerosol therapy).

**Please read the instructions carefully before use.
Drug Administration must be under medical control.**

WARNINGS - GENERAL SAFETY REGULATIONS

**Please read the instructions carefully before use.
Keep the instructions in a safe place.**

- Always observe instructions and safety warnings as indicated in the instructions manuals of devices for aerosoltherapy. The device should be used only original accessories shown in this manual.
- The device has small components which might be removed and easily swallowed. Use by minors and disabled people require presence of an adult with his faculties. Don't leave the device unattended in places easily accessible by minors and disabled people.
- Dispose of device ads provided by standards in force.
- This medical device must be destined exclusively for the use for which it has been designed ad described in this manual.
- Device not suitable for anaesthesia and lung ventilation.
- In presence of pathologies with microbial contamination and infection hazard a single-patient use of accessories and nebulizer is recommended (always consult your doctor).
- The device should be used only with compressor for aerosoltherapy that comply with standards in force.
- Use the "nose piece" accessory only if expressly indicated by your doctor and paying attention **NEVER** to introduce inside the nose the nasal bifurcation, but only bring it as close as possible.
- The materials used to contain the drugs are made with highly stable thermoplastic polymers that are resistant against chemicals. Such materials were tested with commonly used drugs (Salbutamol, Beclametasone dipropionate, Acetylcysteine, Budesonide, Ambroxol) and no interaction phenomenon was observed. Interactions cannot however be excluded given the variety and the continuous evolution of the drugs that are used. Remember to:
 - To consume the drugs as quickly as possible after opening its package;
 - To avoid keeping the drug in the tray-like container for too long and to clean it immediately after every application;
 - If the tray-like container presents any abnormal situation (such as softening or cracks), do not introduce any solution and do not proceed with the inhalation. Contact the

technical service and describe the methods and type of drugs used.

WARRANTY





- HI-4 KIT has a 2-year warranty from the date of purchase against any material or manufacturer defects;
- The warranty replaces/repairs the defective parts;
- The device must only be repaired by an authorised CA-MI Srl technical service centre;
- The warranty shall be null and void if the device is tampered with or the defect is caused by improper use;
- The warranty is valid from the date of purchase of the product on confirmation of the tax receipt or purchase invoice.

TECHNICAL CHARACTERISTICS

Maximum Pressure	300kPa (3,0 Bar)			
Maximum Flow	16 l/min			
Operating Pressure	95 kPa (0,95 Bar)			
Operating Air Flow	8,0 l/min a 95kPa (0,95 Bar)			
Max. capacity	8 ml			
Min. capacity	2ml			
Nebulization (**)	POS.A	POS.B	POS.C	POS.D
	0,40 ml	0,60ml	0,70ml	0,80ml
MMAD (EN 13544-1) *	POS.A	POS.B	POS.C	POS.D
	3.32µm	4.07µm	4.23µm	4.18µm
GSD (EN 13544-1) *	POS.A	POS.B	POS.C	POS.D
	4.12	2.64	2.74	2.49
Connection air tube	ø 5 x 11 (mm)			
Working Condition	Room temperature: 5 ÷ 40 °C Room humidity percentage: 10 ÷ 93% RH Atmospheric Pressure: 700 ÷ 1060 hPa			
Conservation condition and Transport	Room temperature: -25 ÷ 70°C Room humidity percentage: 0 ÷ 93% RH Atmospheric Pressure: 500 ÷ 1060 hPa			

(*) Data determined with Cascade Impactor 290 series, compliant with EN 13544-1, by spraying 2 ml of NaF 1.0%.

(**) Data from free nebulisation of 2 ml of NaCl 0.9% (average nebulisation value per minute).


SYMBOLOLOGY

CE 0123	CE marking in conformity with EC directive 93/42/EEC and subsequent changes
	Consult the instruction manual
	General warnings and/or specifications
	Manufacturer
LOT	Lot Number
REF	Identification device

USE OF THE DEVICE

- Using your fingers, press the two side tabs of the nebuliser cup that hold the top of the nebuliser cup against the bottom of the same component;
- Pour the medicine prescribed by the doctor into the nebulizer;
- Re-close the nebulizer, re-screwing the lid;
- Connect the air tube to the air outlet nozzle of the device
- Connect the desired accessory to the nebulizer: child mask or adult mask, mouth-piece or nosepiece;

Inhalation Valve: Three vents on top are mechanically opened by patient inhalation force, entraining ambient air. The airflow generated by the compressor along with the airstream entrained by patient enhance the nebulisation activity inside the chamber, thus increasing the volume of aerosol delivered to the patient.

Exhalation Valve: While remaining closed during inhalation to prevent aerosol to be wasted in the air, the exhalation valve opens when patient starts exhaling, directing the exhaled air in the ambient and preventing mixing with medication in the chamber.

Speed Selector: By rotating and stopping the speed selector on one of the 4 available positions (I-II-III-IV), speed of nebulization can be doubled while maintaining a stable and optimal particle size.

Always use the nebulizer facing upwards so that substances and / or medicines cannot escape from the nebulizer during the normal use.

NEVER INHALE IN HORIZONTAL POSITION NEVER BEND THE NEBULIZER OVER 60°

CLEANING THE DEVICE

Before using and/or after cleaning, pay special attention to ensure that all the accessories supplied with the device are intact. Switch off the device before cleaning it and disconnect the mains cable from the electrical socket.

For each individual patient it's recommended to use the nebulizer for 6 months or for a maximum of 120 treatments. The nebulizer must be replaced after a long period of inactivity, if it is deformed or broken, or if the nebulizer nozzle is blocked by dry medicine, dust, ecc.. **Only use the original nebulizer supplied by CA-MI with the device**

In the presence of infection or microbial contamination prone pathologies, we recommend using your personal accessories and nebulizer (always consult your doctor).

CLEANING

Before and after each use proceed with cleaning all of the components of the nebulizer (with the exception of the air tube) according to one of the two methods described below.

Method 1: Thoroughly clean the components for 5 minutes, using warm drinking tap water (about 40°C) and/or mild soap.

Method 2: Clean the components (except for the air tube) by immersing them in a solution with 60% water and 40% white vinegar. When finished, thoroughly rinse with warm drinking water (approx. 40°C).

After cleaning, rinse thoroughly by removing the excess water and allow to air dry in a clean place.

**DO NOT BOIL OR AUTOCLAVE THE AIR TUBE
AND MASKS
DO NOT WASH ACCESSORIES IN A DISHWASHER**

WASHING

If there are pathologies with risks of infection and microbial contamination, it is the end user's responsibility to proceed with suitable washing. The washing procedure can only be carried out if the components to be treated have undergone specific cleaning (see chapter on cleaning).

Proceed as follows for the washing procedure:

- Fill a container, of a suitable size to contain all the individual components, with a solution of drinking water and disinfectant (hypochlorite-based solution readily available in a pharmacy) by following the proportions indicated on the packaging of the disinfectant itself;
- The period of time for which is to be immersed in this solution is indicated on the packaging of the hypochlorite solution in accordance with the chosen concentration for preparing the solution;
- Rinse thoroughly with lukewarm drinking water to remove all traces of the solution. Dry and store in a dry, dust-free environment.

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- El dispositivo tiene que ser reparado solamente por el servicio técnico autorizado de CA-MI Srl;