GIMA

CAVI BIPOLARI
BIPOLAR CABLES
CABLES BIPOLAIRES
ZWEIPOLIGE KABEL
CABLES BIPOLARES
FIOS ELÉTRICOS BIPOLARES
ΔΙΠΟΛΙΚΑ ΚΑΛΩΔΙΑ
أسلاك ثنائيّة القطب

MANUALE D'USO E MANUTENZIONE
USE AND MAINTENANCE BOOK
INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN
BETRIEBS UND WARTUNGSANWEISUNGEN
MANUAL DE USO Y MANTENIMIENTO
MANUAL DE USO E MANUTENÇÃO
ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ
دليل الإستعمال والرعاية

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FEATURES

The GIMA bipolar cables are devices used with bipolar electro-surgery devices. The cables are made of light and flexible silicone, they are 3 meters long and can be sterilized in an autoclave up to 137°. They have a universal connection that makes them easy to use with other GIMA devices and with most of the electric scalpels produced in Europe.

INFORMATION

The product is designed to be used only by qualified doctors, specialized in electro-surgery techniques and the maintenance of this type of product.

The incorrect use or poor or inappropriate maintenance can make the product deteriorate quickly.

These products can cause injuries if not used correctly!

Do not put in contact with devices that are not insulated!

Do not use electro-surgical devices on patients with Pacemakers or near inflammable liquids or anesthetics!

The electric scalpels used together with these cables will destroy tissues and are extremely dangerous if used incorrectly. Carefully read the user instructions supplied by the manufacturers of the electric scalpels.

The product must be sterilized before the first use: we recommend sterilization in the autoclave up to 137°.

UNPACKING

Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled.

They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

USE

Before using the device make sure the cables are fully compatible with the electro-surgical devices.

Connect the connectors to the chosen surgical device and then the feeding device and verify that the connections are correct and that no metallic parts are left exposed.

Do not bending the cable or fix it in any way with clips or similar devices that could damage its functionality.

Wrong connections and/or improper use could cause shocks: do not leave the cable attached to the feeding device without first positioning the surgical device.

Before replacing the surgical instrument disconnect the cable from the feeding device.

MAINTENANCE

It is very important to examine the product to identify any breaks, defects or worn parts. Pay special attention to the connectors. The cables must be kept in a dry and clean place.

If the product is defective or too worn, do not use or try to repair it.

After use the product must be cleaned thoroughly with a proper solution. Remove blood and any other impurities using a soft cloth and rinse properly.

Before sterilizing the device, run water out of the connectors and dry the cable.

Do not use ultrasound or mechanical cleaning methods because they could damage the product or shorten its technical life.

Be sure the product does not come into contact with acids or other disinfectants which could corrode it.

Sterilization with ethylene oxide

Roll up the cable without pressing or twisting and put it in the ethylene oxide container (ETO).

Sterilization in the autoclave

Make sure the cable does not come into contact with the metal parts of the sterilizing unit. Do not exceed the temperature of 137°C.

Sterilization with hot air

Hot air sterilization or other systems that may damage the cable are not recommended.

GUARANTEE

Thank you for having acquired our product.

This product meets the most stringent requirements regarding the selection of quality manufacturing materials and also the final control. The product has a 12-month warranty, valid from the date it is delivered by GIMA.


During the guarantee period free repair and/or replacement of any defective parts due to faulty manufacture will be given, labour, postal, transport and packaging charges etc. are not included. The guarantee therefore excludes components subject to wear and tear such as parts in rubber or PVC, etc. No compensation will be given for loss of use of the product.

Furthermore replacement or repairs effected during the guarantee period do not lengthen the duration of the guarantee.

This guarantee is void in the event of: repairs effected by unauthorised persons or with spare parts not approved by GIMA, damage or defects caused by negligence, blows, abnormal use of the appliance or faulty installation.

The guarantee will be invalidated if the registration number has been removed, cancelled or altered. Faulty appliances must be returned only to the retailer where the item was bought. Any item dispatched directly to us will be refused.

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