SINCE 1966
MADE IN ITALY

PulmoWaves²

by **FLAEM**[®]
PROLINE

IT> ISTRUZIONI PER L'USO
EN> INSTRUCTIONS FOR USE
FR> INSTRUCTIONS POUR L'UTILISATION
ES> INSTRUCCIONES PARA EL USO
DE> GEBRAUCHSANWEISUNG

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Device for respiratory rehabilitation physiotherapy

These operating instructions are provided for device model P0112EM F1000 and aerosol therapy device accessory model P0112EM.

INTENDED USE. Medical device for the removal of tracheobronchial secretions and for respiratory physiotherapy, the treatment of which must be prescribed by a doctor.

INDICATIONS FOR USE. • Treatment of conditions of the respiratory system, in particular for the removal of tracheobronchial secretions (COPD, cystic fibrosis, bronchiectasis and/or as indicated by the treating physician).

• Aerosol therapy can be combined with prescription.

⚠ CONTRAINDICATIONS. • The medical device should NOT be used for patients who are unable to breathe on their own or who are unconscious. • Do not use the device in anaesthetic or assisted ventilation circuits. • The device must not be used on infants. • The device must not be used on children under 3 years of age.

• The device should not be used on patients with cardiovascular instability.

INTENDED USERS. The devices are intended for use by legally authorised medical personnel/health workers (doctors, nurses, therapists, etc.). The device can be used directly by the patient for home therapy, but under the guidance of medical personnel or the legally authorised healthcare professional.

⚠ TARGET GROUP OF PATIENTS. Adults, children over 3 years of age. Before using the device, carefully read the user manual and ensure the presence of an adult responsible for safety if it is intended for use by children or people with limited capabilities (for example: physical, mental or sensory). It is up to the medical personnel to assess the patient's condition and capabilities in order to determine, when prescribing therapy, whether the patient is able to operate the device safely on their own or whether the therapy should be administered by a responsible person. Please refer to medical personnel for the evaluation of the use of the device on particular types of patients such as pregnant women, lactating women, incapacitated persons or persons with limited physical capabilities.

OPERATING ENVIRONMENT. This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., or even at home.

⚠ WARNINGS ABOUT ANY MALFUNCTIONS. • Should your device fail to perform, please contact the authorised service centre for clarifications. • The manufacturer should be contacted to report problems and/or unexpected events relating to operation and if necessary for clarification of use and/or maintenance/hygienic preparation. • Please also refer to the case history of faults and their resolution.

WARNINGS. • Only use the device for its intended purpose. This medical device is not intended as a life-saving device. Any other use is considered improper and can be dangerous. The manufacturer is not liable for any improper use. - Always consult your doctor when determining the dosage of therapy. • Follow your doctor's or therapist's instructions regarding the possible use of types of medicine, dosage and indications for aerosol treatment. • Pay attention to the indications provided with the drug and avoid using the devices with substances and dilutions other than those recommended. • The apparatus must be transported inside the supplied bag and enclosed in its protective pouch. • The patient is also considered to be the operator.

• If you should experience allergic reactions or other problems while using the device, stop using it immediately and consult your doctor. • Keep this manual safe for further reference. • If the packaging is damaged or opened, contact the distributor or service centre. • Do not expose the device to particularly extreme temperatures. • The time required to switch from storage to operating conditions is about 2 hours.

• It is forbidden to access the device opening in any way. • Repairs may only be carried out by personnel authorised by the manufacturer. • Unauthorised repairs invalidate the warranty and may pose a danger to the user.

• **Choking hazard:** Some parts of the device are small enough to be swallowed by children; keep the device out of the reach of children.

• **Strangulation hazard:** Do not use the supplied connecting hose and cables outside of their intended use, they could cause a strangulation hazard, take special care for children and people with special difficulties, often these people are not able to correctly assess the dangers.

• **Fire hazard:** - This device is not suitable for use in the presence of an anaesthetic mixture that is flammable with air, or with oxygen or nitrous oxide. - Never obstruct the ventilation slots located on both sides of the device.

• **Electrocution hazard:** - Before using the product for the first time, and periodically during its service life, check the integrity of the device structure and of the power cable to make sure there is no damage. In the event of damage, do not plug in the cable and immediately take the product to an authorised service centre or to your trusted dealer.

- Keep the power cable away from animals (e.g. rodents), otherwise such animals could damage the insulation of the power cable.

- Always keep the power cable away from hot surfaces. - Do not handle the plug and power cord with wet hands. - Do not handle the apparatus with wet hands. - Do not use the device in damp environments (for example, while taking a bath or shower). - Do not immerse the appliance in water; if this happens, pull the plug immediately. - Do not remove or touch the device immersed in water; unplug the power cable first. Immediately take the device to an authorised FLAEM service centre or to your trusted dealer.

- No maintenance/cleaning must be carried out during use.

• **Risk of ineffectiveness of therapy:** - Only use original Flaem spare parts, no liability is accepted if non-original parts are used - Observe the hygienic preparation of the spray chamber as described under Hygienic Preparation.

- Only use the device in a dust-free environment, otherwise therapy may be impaired.

- Do not obstruct or insert objects into the filter and its housing in the unit. - If aerosol therapy is used with substances that are too dense, it may be necessary to dilute them with a suitable physiological solution, as prescribed by your doctor.

- When using aerosol therapy, use the nebulizer unit in the correct position, preferably vertical. Do not tilt the nebuliser beyond an angle of 30 degrees in any direction to avoid the drug flowing into the mouth or being excessively dispersed, reducing the effectiveness of the treatment.

• **Risk of infection:** - Follow hygienic preparation before and after each use. Make sure that the parts are not stored near accessories or devices for other therapies (e.g. infusion devices).

- The supplied components are single-patient use. - Once the therapy is finished, do not leave the medicine inside the ampoule and proceed with the hygienic preparation. - If the nebuliser is used for several types of medicine, the residues must be completely removed. Perform hygienic preparation after each inhalation, also to achieve the highest degree of hygiene and to optimise the service life and operation of the device.

• **Risk of injury:** - Do not place the device on a soft supporting surface such as a sofa, a bed or a tablecloth. - Always operate it on a hard surface clear of any objects.

WARNINGS ON INTERFERENCE RISKS DURING USE IN DIAGNOSTIC INVESTIGATIONS

This device is designed to meet the current requirements for electromagnetic compatibility. As far as EMC requirements are concerned, electromedical devices require special care when being installed and used, and are therefore required to be installed and/or used in accordance with the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. The device may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnosis or treatment. For more information visit www.flamnuova.it.

TROUBLESHOOTING

Before performing any operation, switch off the apparatus and disconnect it from the mains.

Problem	Cause	Remedy
The device does not work.	The power cable (A12) has not been plugged properly into the device power socket (A11) or into the mains socket.	Plug the power cable correctly into the socket (A11) and into the relative mains socket.
The nebuliser (B2) does not nebulise and/or dispenses the drug correctly.	The switch (B1) of the aerosol unit (B) is on "OFF". The cable connector (B4) has not been plugged into the 12VDC socket (A5). The cable (B4) has not been plugged into the aerosol unit (B).	Set the switch (A1) on the main unit (A) to "ON". Insert the cable connector (B4) into the socket (A5). Insert the cable (B4) into the aerosol unit (B).
	The medication or saline solution has not been inserted into the nebuliser (B2) (bowl).	Pour the correct amount of medication or saline solution into the nebuliser (B2).
	The nebuliser (B2) has not been properly fitted.	Disassemble the nebuliser (B2) and reassemble it correctly as shown in the diagram on the cover.
	The nebuliser nozzle (B2) is clogged	Dismantle the atomiser (B2), remove the upper part and the valve support (B2.2+B2.3) and clean. Failure to remove medication deposits from the nebuliser (B2) affects its efficiency and proper operation. Strictly follow the instructions in the Hygienic Preparation chapter.
The digital pressure gauge (C) does not light up the LEDs and/or does not work properly.	The nebuliser (B2) does not nebulise and/or dispenses the drug correctly.	Insert the power cable correctly into the socket (A11) and into the relative mains socket.
	The air tubes (A9) (A10) are not properly connected to the main unit (A) and/or the PEP connector (C3).	Check that the tubes are correctly connected to the relative sockets of the main unit (A). See the connection diagram on the cover.
	Air tubes (A9) (A10) are bent and/or damaged and/or kinked. The ring nut (C1) has not been positioned correctly.	Uncoil the tubes (A9) (A10) and check them for crushing, punctures or kinks; if they are, replace them. Correctly adjust the PEP ring nut (C1) by matching the holes.
	The flow regulator filter (A3) is dirty.	Replace the filter (A3).
	Filter (A8) clogged or dirty, bent tube (A9) and/or faulty circuit.	Restore the tube (A9), replace the filter (A8) and/or contact an authorised service centre to have the electronic circuit replaced.
	The electronic circuit does not work	Contact a FLAEM service centre to replace the electronic circuit.

When I inhale, the compressor does not activate and no pulsed air comes out of the mouthpiece.	The power cable (A12) has not been inserted properly into the device power socket (A11) or into the mains socket.	Insert the power cable (A12) correctly into the socket (A11) and into the relative mains socket.
	The switch (A1) is set to "OFF". The pressure switch is faulty.	Set switch (A1) to "ON" or contact a FLAEM service centre to replace the pressure switch.
	The air tubes (A9) (A10) are not properly connected to the unit (A) and/or the PEP connector (C3)	Check that the tubes are correctly connected to the relative sockets of the device (A). See the connection diagram on the cover.
	The flow regulator filter (A3) is dirty.	Replace the filter (A3).
The device (A) is very noisy.	The pressure regulator (A3) is not correctly inserted in its housing	Properly insert the pressure regulator (A3) in its housing.
If, after checking the conditions described above, the device still is not working properly, we recommend that you contact your trusted dealer or an authorised FLAEM service centre nearest to you.		

DISPOSAL.

 In accordance with Directive 2012/19/EC, the symbol on the appliance indicates that the apparatus to be disposed of (excluding accompanying components) is considered as waste and must therefore be subject to "separate collection". Therefore, the user must deliver this waste (or have it delivered) to the separate collection centres set up by the local authorities, or hand it over to the retailer when purchasing a new appliance of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment from recycled materials and limit the negative environmental and health effects caused by improper waste management. Unauthorised disposal of the product by the user entails the application of the administrative sanctions provided for in the laws transposing Directive 2012/19/EC of the member state or country where the product is disposed of.

Packaging

 Product box	 Tube packaging bag	 Heat-shrink film	 Sack packaging bag
--	--	--	--

REPORTING SERIOUS EVENTS. Please report any serious incidents occurring in connection with this device to the manufacturer and to the competent authority of the Member State where you reside. An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

SYMBOLS ON DEVICE OR PACKAGING

 Unique device identifier	 Serial number of the device	 Class II device	 Manufacturer
 Type BF applied part	 Production date	 Batch Code	 Alternating current
 Attention	 Phthalate- and bisphenol-free	 Keep dry	 Model number
 Medical device	 Temperature limits	 Humidity limits	 Atmospheric pressure limits
 "ON"	When the device is switched off, the on/off switch stops compressor operation in only one phase of the two-phase power supply.	 Before use: Caution check instructions for use	
 "OFF"		 Single patient	
IP21 Protection rating of the envelope: IP21. (Protected against solid bodies larger than 12 mm. Protected against access with a finger			
CE0051 Medical CE marking ref. regulation 2017/745 EU as amended			
OUT 12V Output Power for Nebuliser			

INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES. • Interactions: The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely excluded. • It aerosol therapy is utilized, use the medicine as soon as possible once it has been opened and avoid leaving it in the nebuliser; once therapy has ended, do not leave the medicine inside the nebuliser and proceed with hygienic preparation.

TECHNICAL CHARACTERISTICS MAIN UNIT (A)

Power supply/power	230 V ~ / 140 VA
Frequency	50Hz
Fuse	1x T2A 250V
Noise (at 1 m)	55 dB (A) (approx.)

TECHNICAL CHARACTERISTICS AEROSOL UNIT

Power supply	12 V
Nebulisation ml/min ⁽¹⁾	0,32 ml/min'
MMAD (2)	4,18 µm
Breathable fraction <5 µm ⁽²⁾	58,9 % (approx.)
Medicine tray capacity	8 ml

⁽¹⁾ Data collected according to Flaem's internal procedure. ⁽²⁾In vitro characterisation performed by TÜV Rheinland Italia S.r.l. in cooperation with the University of Parma. More details are available on request.

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ENVIRONMENTAL CONDITIONS

Operating conditions:

Ambient temperature: Between +5°C and +40°C
Relative air humidity: Between 10% and 95%
Atmospheric pressure: Between 69 KPa and 106 KPa

Storage and transport conditions:

Ambient temperature: Between -25°C and +70°C
Relative air humidity: Between 10% and 95%
Atmospheric pressure: Between 69 KPa and 106 KPa

APPLIED PARTS

Type BF applied parts are: patient accessories (C2)

DURATION

Main unit (A) excluding pipes (A7, A9, A10)	Service life 1000 hours.
Aerosol unit (B) excluding nebuliser (B2)	Service life 400 hours
For components (C1, C2, C4, C5)	The average life expectancy is 1 year, however, it is advisable to replace them every 6 months in intensive use to ensure maximum therapeutic efficacy at all times.
The component (A8 is disposable)	

APPLIANCE AND MATERIAL INFORMATION

A Main unit A1 Switch A2 Digital pressure gauge A3 Flow regulator with filter A4 Pulsed air outlet A5 12Vdc output power socket for aerosol unit A6 Pressure switch inlet A7 Filter tube (single-patient) A8 Filter (disposable). It is not a FLAEM-manufactured device A9 Pressure switch tube A10 Pulsed air tube A11 Power supply socket A12 Power cable A13 Tube connector clip A14 Protective gasket	B Aerosol unit accessory (compressor + nebuliser) B1 Switch B2 Nebulizer B2.1 Lower part B2.2 Upper part B2.3 Support with valve B3 Aerosol unit power socket B4 Power cord for aerosol therapy unit B5 Cable plug B4 C Accompanying component (single-patient) C1 PEP adjustment ring nut C2 Mouthpiece C3 PEP fitting C4 Nose clip C5 Valve for use without nebuliser device D Comfortable, spacious carry bag
--	---

Polypropylene

INSTRUCTIONS FOR USE

Before each use, wash your hands thoroughly and clean your appliance as described in section "HYGIENIC PREPARATION". The parts (A7-A9-A10-C2-C3) are for single-patient use only to avoid potential risks of transmitting infection. It is recommended to replace them every six months. The A8 filter is disposable.

Insert power supply cord (A12) into socket (A11) and then connect it to a mains outlet adequate to the device voltage. The position of the socket must be such that the device can be easily unplugged from the mains network.

Tube assembly (Fig. 1)

- Connect tube A10 to the pulsed air outlet A4 and tube A7 to the pressure switch inlet.
- Insert the plug (B5) of the cable (B4) into the socket (A5).
- Screw the ring nut by turning it clockwise.
- Connect the filter (A8) to the tube (A7).

ATTENTION. The protection cap of the 15Vdc socket must be removed to insert the plug and must be replaced after use to protect the socket.

Mouthpiece (Fig. 2)

- Connect mouthpiece C2 to PEP connector C3.
- Connect the two tubes A9-A10 coming out of A4-A8 to PEP connector C3.

Device connection (Fig. 3)

- Insert the coaxial cable plug (B4) coming out of (A5) into socket (B3) of the nebuliser device.
- Insert the protective gasket (A14) into the power supply cable plug (A12).
- Place the plug on power supply cable (A12) into socket (A11) on the device **A**.

Nebuliser device (Fig. 4)

- Open the upper part (B2.2) of the nebuliser (B2) counterclockwise and pour the saline solution into (B.2.1).
- Apply the valve support (B.2.3) to the lower part (B.2.2).
- Close the nebuliser (B2) of Nebuliser Device **B** and connect the fitting (C3) and mouthpiece (C2) as shown in the connection diagram.

Use without nebuliser device (Fig. 5)

It is possible to use the Pulmowaves device without the nebuliser, by closing the PEP fitting (C3) and with the valve (C5) as shown in the picture.

Activating the device (Fig. 6)

- Turn on the device by pressing power switch A1, then press button B1 of the nebuliser.
- The compressor will start only with the patient's inhalation through mouthpiece C2.

Flow regulator (Fig. 7)

- Adjust ring nut **C1** of PEP fitting **C3** and flow regulator **A3** of the device for the desired respiratory physiotherapy according to the values reported in the following table.

Adjus. ring fitt. PEP	Flow regulator set on MAX		Flow regulator set on MIN	
	Val. Flow	Val. Press	Val. Flow	Val. Press
Hole diameter				
Ø 2	7.5 l/min	7 mbar	1.2 l/min	2 mbar
Ø 3	6.4 l/min	6 mbar	1.1 l/min	1 mbar
Ø 4	3.9 l/min	3 mbar	1.0 l/min	0 mbar
Ø 5	3.3 l/min	2 mbar	0.8 l/min	0 mbar
Ø 6	3.2 l/min	1 mbar	0.6 l/min	0 mbar

Digital manometer (Fig. 8)

- During the patient's exhalation phase, the pressure gauge (A2) displays the positive pressure generated according to the hole selected on the ring nut (C1) of the PEP connector (C3).
- The pressure gauge is only designed to indicate the correct patient exhalation. It is not an indicator of the precise value of exhalation flow.

Use keys "+" and "-" to set the therapy duration. The therapy may last 5, 10, 15 or 20 minutes. At the end of the set therapy, the 4 LEDs flash and an acoustic warning beeps 5 times. (Fig. 9)

Pipe connection clip (Fig. 10)

Attach tubes A9-A10 and cable B4 to the proper hooks of the Clip A13 as shown in the detail.

Nose plug (Fig.11)

For the correct use of the PEP connector, the patient must put on nose clip C4, as shown in the figure.

HYGIENIC PREPARATION

Switch off the device before any cleaning procedures and always unplug the power cable from the socket.

Appliance (A) and aerosol unit (B) and outside tube (B4)

Use only a cloth moistened with antibacterial detergent (non-abrasive and free of solvents of any kind).

Supplied components and accessories

Accessories A8-C2-C3-C4 and connection tubes A7-A9-A10 are for personal, single-patient use to avoid any risks of transmitting infection.

Open the nebuliser (B2) by turning the upper part anticlockwise, remove the upper part as shown in fig.12. NOTE: For sanitisation, disinfection and sterilisation operations the PEP C3 connector must not be disassembled.

Filter

ATTENTION

- The filter (A8) is for single use and must be replaced as indicated on the package label.
- It is washable and can be disinfected or sterilised and is not recyclable.

CLEANING AT HOME - SANITISATION AND DISINFECTION

Sanitisation

Before and after each use, sanitise the accessories and the nebuliser by choosing one of the methods described below.

Method A: Sanitise the accessories under warm drinking water (approx. 40°C) with mild dishwashing detergent (non-abrasive).

Method B: Sanitise the accessories in the dishwasher with the hot cycle (70°C).

Method C: Sanitise the accessories by soaking in a solution of 50% water and 50% white vinegar, then rinse thoroughly with warm drinking water (approx. 40°C).

If you also wish to perform cleaning by DISINFECTION, go to the DISINFECTION paragraph.

Disinfection

After sanitising the accessories, disinfect them choosing one of the methods described below.

method A: Obtain a disinfectant of the electrolytic chloride type (active ingredient: sodium hypochlorite), specifically for disinfection, available in all pharmacies.

Execution:

- Fill a container big enough to hold all of the parts being disinfected with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution according to the instructions provided by the disinfectant manufacturer.

method B: Disinfect accessories by boiling them in water for 10 minutes; use demineralised or distilled water to avoid hard water deposits.

method C: Disinfect the accessories in a hot bottle steriliser (not microwave). Carry out the process strictly according to the hot bottle steriliser's instructions. For disinfection to be effective, choose a steriliser with an operating cycle of at least 6 minutes.

After the CLEANING AND/OR DISINFECTION operations dry the area of the transducer, without exerting unnecessary pressure on the transducer itself, dry the accessories and the outside of the appliance with a paper towel or even better with a jet of hot air (e.g. hair dryer). Finally, check that the "power plug" area of the device is completely dry.

CLEANING IN A CLINICAL OR HOSPITAL SETTING - DISINFECTION AND STERILISATION

Before disinfection or sterilisation, sanitise the accessories, choosing one of the methods described below.

method A: Sanitise the accessories under warm drinking water (approx. 40°C) with mild dishwashing detergent (non-abrasive).

method B: Sanitise the accessories in the dishwasher with the hot cycle.

DISINFECTION

Obtain a disinfectant of the electrolytic chloride type (active ingredient: sodium hypochlorite), specifically for disinfection, available in all pharmacies.

Execution:

- Fill a container big enough to hold all of the parts being disinfected with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution according to the instructions provided by the disinfectant manufacturer.

If you also wish to perform STERILISATION, jump to the STERILISATION paragraph.

After the CLEANING AND / OR DISINFECTION operations dry the area of the transducer, without exerting unnecessary pressure on the transducer itself, dry the accessories and the outside of the appliance with a paper towel or even better with a jet of hot air (e.g. hair dryer). Finally, check that the "power plug" area of the device is completely dry.

STERILISATION

Device: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

Execution: Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Place the packaged components in the steam steriliser.

Run the sterilisation cycle according to the operating instructions of the device by selecting a temperature of 134°C and an initial time of 10 minutes.

Storage: Store sterilised components according to the instructions for use of the system or sterile barrier packaging chosen.

The sterilisation procedure is validated in its conformity to ISO 17665-1.

Table of expected methods													
Attachment		C1	C2	C3	C4	C5	B2	A7	A9	A10	B2.1	B2.2	B2.3
Sanitisation	A	✓	✓	✓	✓	✓	-	✓	✓	✓	✓	✓	✓
	B	✓	✓	✓	✓	✓	-	-	-	-	-	-	-
	C	✓	✓	✓	✓	✓	-	300	300	300	300	300	300
Disinfection	A	300	300	300	300	300	-	300	300	300	300	300	300
	B	300	300	300	300	300	-	300	300	300	300	300	300
	C	300	300	300	300	300	-	300	300	300	300	300	300
Sterilisation		300	300	300	300	300	300	300	300	300	300 NEBULISER ASSEMBLED (B2)		
✓: expected 300: MAX 300 TIMES													

AIR FILTERING

The device is fitted with a suction filter, located inside the flow regulator (A3), which must be replaced when it is dirty or discoloured. To replace it, lift the flow regulator (A3) and extract the Filter.

Do not wash or reuse the same Filter. Regular replacement is necessary to assure proper compressor performance. The Filter must be regularly checked. Contact your dealer or authorised service centre for replacement filters.

SPARE PARTS

Description of parts	KIT code				
	ACO662P Six-month pack consumables	ACO663P Patient pack for single use recomm.	ACO664P Outpatient pack recomm.	ACO665P Outpatient pack daily	ACO001C Single- patient
Set of tubes shq L.800, L.690, L.90	1	1	-	3	-
Anti-bacterial filter	24	1	1	-	-
Mouthpiece	1	1	1	-	-
PEP dispenser	1	1	1	-	-
Nose clip	1	1	1	-	-
Clip for external tubes	5	5	-	15	-
Manual	1	1	1	1	-
Porex filter for compressor	1	-	-	-	-
RF8 Plus	-	-	-	-	1

Only use original Flaem accessories.

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© 2025 FLAEM NUOVA®
FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.
FLAEM NUOVA S.p.A.