
IT > MANUALE ISTRUZIONI
EN> USER MANUAL
DE> BETRIEBSANLEITUNG
FR > MODE D'EMPLOI
NL > GEBRUIKSAANWIJZING
PL > INSTRUKCJA OBSŁUGI
SV > BRUKSANVISNING
ES > MANUAL DEL USUARIO
EL > ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ ΧΡΗΣΗΣ

Schema di collegamento - Connection diagram - Anschlussschema - Schéma de
branchement - Verbindungsschema - Schemat połączeń - Monteringschema -
Esquema de conexión - Διαγράμμα σύνδεσης

Apparecchio per aerosolterapia

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Aerosol therapy device

These user instructions are provided for FLAEM JUNIOR F2000 and RF7-2 model devices. The FLAEM Aerosol Therapy device consists of a compressor unit (A), a nebuliser and some accessories (B).

INTENDED USE. Medical device for the administration of medication through inhalation, with inhalation therapy and medication prescribed by a Doctor.

INDICATIONS FOR USE. Treatment of respiratory diseases. Medications must be prescribed by a doctor who has assessed the patient's general condition.

⚠ CONTRAINDICATIONS. • The medical device should NOT be used for patients who are unable to breathe on their own or who are unconscious. • Do not use the device in anaesthetic or assisted ventilation circuits.

INTENDED USERS. The devices are intended for use by legally authorised medical personnel/health workers (doctors, nurses, therapists, etc.). The device can be used directly by the patient.

⚠ TARGET GROUP OF PATIENTS. Adults, children of all ages, infants. Before using the device, the user manual must be read carefully and an adult responsible for safety must be present if the device is used by infants, children of any age or persons with limited abilities (e.g. physical, mental or sensory). It is up to the medical personnel to assess the patient's condition and capabilities in order to determine, when prescribing therapy, whether the patient is able to operate the aerosol safely on their own or whether the therapy should be administered by a responsible person.

Please refer to medical personnel for the evaluation of the use of the device on particular types of patients such as pregnant women, lactating women, infants, incapacitated persons or persons with limited physical capabilities.

OPERATING ENVIRONMENT. This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., or even at home.

⚠ WARNINGS CONCERNING POSSIBLE MALFUNCTIONS. • Should your appliance fail to perform, please contact the authorised service centre for clarification. • The manufacturer should be contacted to report problems and/or unexpected events relating to operation and if necessary for clarification of use and/or maintenance/hygienic preparation. • Please also refer to the case history of faults and their resolution.

WARNINGS. • Use the device only as a therapeutic inhaler. This medical device is not intended as a life-saving device. Any other use is considered improper and can be dangerous. The manufacturer is not liable for any improper use. • Always consult your general practitioner for identification of treatment. • Follow the instructions of your doctor or respiratory rehabilitation therapist regarding the type of medicine, dosage and treatment indications. • If you should experience allergic reactions or other problems while using the device, stop using it immediately and consult your doctor. • Keep this manual safe for further reference. • If the packaging is damaged or opened, contact the distributor or service centre. • Do not expose the device to particularly extreme temperatures. • Do not place the device near heat sources, in sunlight or in an excessively hot environment. • The time required to switch from storage to operating conditions is about 2 hours. • It is forbidden to access the compressor unit opening in any way. Repairs may only be carried out by personnel authorised by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a danger to the user.

Suffocation risk: Some components of the appliance are small enough to be swallowed by children, so keep the appliance out of the reach of children.

Strangulation risk: Do not use the supplied connecting hose and cables outside of their intended use, they could cause a strangulation hazard, take special care for children and people with special difficulties, often these people are not able to correctly assess the dangers.

Fire risk: This device is not suitable for use in the presence of an anaesthetic mixture that is flammable with air, or with oxygen or nitrous oxide.

Risk of electrocution: - Before first use, and periodically during the life of the product, check the integrity of the structure of the appliance and the power cord to ensure that there is no damage; if it is damaged, do not plug it in and immediately take the product to an authorised service centre or your dealer. - Keep the power cable away from animals (e.g. rodents), otherwise such animals could damage the insulation of the power cable. - Always keep the power cable away from hot surfaces.

- Never obstruct the ventilation slots located on both sides of the compressor unit.

- Do not handle the compressor unit with wet hands. Do not use the compressor unit in humid environments (e.g. while bathing or showering). Do not immerse the compressor unit in water; if this happens, pull the plug immediately. Do not pull out or touch the compressor unit while immersed in water, un-

plug it first. Take it immediately to an authorised FLAEM service centre or your dealer.

Risk of ineffectiveness of therapy: - Performance may vary with particular types of drugs (e.g. those with high viscosity or in suspension). For further information, please refer to the package leaflet provided by the drug manufacturer. - Use the nebuliser in the correct position, as upright as possible; do not tilt the nebuliser beyond an angle of 30 degrees, in any direction, to prevent the drug from spilling into the mouth or being overdispersed, reducing the effectiveness of the treatment. - Pay attention to the indications provided with the drug and avoid using the devices with substances and dilutions other than those recommended. - Only use the device in a dust-free environment, otherwise therapy may be impaired. - Do not obstruct or insert objects into the filter and its housing in the unit.

- In the case of substances that are too dense, dilution with suitable saline solution may be necessary, as prescribed by a doctor. - Only use original Flaem accessories or spare parts, no liability is accepted if non-original parts or accessories are used.

Risk of infection: - We recommend personal use of the accessories to avoid any risk of infection.

- Follow hygienic preparation before each use. Ensure that tubing and accessories are not stored in the vicinity of other accessories or devices for different therapies (e.g. infusions). - At the end of therapy, do not leave the medicine inside the nebuliser. - If the nebuliser is used for several types of medicine, the residues must be completely removed. Therefore, perform hygienic preparation after each inhalation, also to achieve the highest degree of hygiene and to optimise the service life and operation of the device.

Risk of injury: - Do not place the device on a soft supporting surface such as a sofa, a bed or a tablecloth. - Always operate it on a hard surface clear of any objects.

WARNINGS ON INTERFERENCE RISKS DURING USE IN DIAGNOSTIC INVESTIGATIONS. This device is designed to meet the current requirements for electromagnetic compatibility. As far as EMC requirements are concerned, electromedical devices require special care when being installed and used, and are therefore required to be installed and/or used in accordance with the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. The device may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnosis or treatment. For more information visit www.flaemnuova.it.

CASE HISTORY OF FAULTS AND THEIR RESOLUTION

Before performing any operation, switch off the appliance and disconnect it from the mains.

Problem	Cause	Remedy
The device is not working.	Power cable not correctly plugged into the appliance socket or mains power socket.	Plug the power cable correctly into the sockets.
The device does not nebulise or nebulises poorly.	The medicine was not inserted into the nebuliser.	Pour the right amount of medicine into the nebuliser.
	The nebuliser has not been assembled correctly.	Disassemble and reassemble the nebuliser correctly according to the connection diagram on the cover.
	The nebuliser is obstructed.	Hygienic preparation of the nebuliser. Deposits of medicine due to lack of hygienic preparation of the nebuliser impair its efficiency and function. Strictly follow the instructions in the HYGIENIC PREPARATION chapter.
	The accessories are not properly connected to the appliance.	Check the correct connection between the air inlet of the appliance and the accessories (see connection diagram on cover).
	The pipe is bent or damaged or twisted.	Unwind the pipe and check it for crushing or punctures. Replace it if necessary.
	The air filter is dirty.	Replace the filter
The device is louder than usual.	The filter is inserted.	Insert the filter properly into the housing.

If, after checking the conditions described above, the device still is not working properly, we recommend that you contact your trusted dealer or an authorised FLAEM service centre nearest to you. You can find a list of all Service Centres at <http://www.flaemnuova.it/it/info/assistenza>.

DISPOSAL

Compressor unit

In accordance with Directive 2012/19/EC, the symbol on the equipment indicates that the equipment to be disposed of (excluding accessories) is considered as waste and must therefore be subject to 'separate collection'. Therefore, the user must deliver this waste (or have it delivered) to the separate collection centres set up by the local authorities, or hand it over to the retailer when purchasing a new appliance of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment from recycled materials and limit the negative environmental and health effects caused by improper waste management. Unauthorised disposal of the product by the user entails the application of the administrative sanctions provided for in the laws transposing Directive 2012/19/EC of the member state or country where the product is disposed of.

Nebuliser and accessories

To be disposed of as general waste after a sanitisation cycle.

Packaging

Product box

Heat-shrinkable nebuliser film and accessories

Product packaging bag and tube

Sack packaging bag

NOTIFICATION OF SERIOUS EVENTS. Please report any serious incidents in connection with this device to the manufacturer and the competent authority of the Member State in which you reside. An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES

• Interactions: The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely excluded. • Use the medicine as soon as possible once it has been opened and avoid leaving it in the nebuliser; once therapy has ended, do not leave the medicine inside the nebuliser and proceed with hygienic preparation.

SYMBOLS ON DEVICE OR PACKAGING

 Unique device Identifier	 Serial number of the device	 Class II device	 Manufacturer
 Type BF applied part	 Production date	 Batch Code	 Alternating current
 Attention	 Phthalate- and bisphenol-free	 See instructions for use	 Model number
 Medical device	 Temperature limits	 Moisture limits	 Atmospheric pressure limits
 Function switch off	 Before use: Caution check instructions for use		
 Function switch on	 Medical CE marking ref. regulation 2017/745 EU and subsequent updates		
IP21 Protection rating of the envelope: IP21. (Protected against solid bodies larger than 12 mm. Protected against access with a finger; Protected against vertically falling drops of water).			

COMPRESSOR UNIT SPECIFICATIONS. Model: FLAEM JUNIOR F2000

Power supply:	230V~ 50Hz 210VA	Noise level (at 1 m):	57 dB (A) approx
Max. pressure:	3.5 ± 0.5 bar	Operation:	Continued
Air flow to compressor:	14 l/min approx	Dimension:	20(L) x 18(W) x 11(H) cm

NEBULIZER TECHNICAL SPECIFICATIONS. Model: RF7-2

RF7 Dual Speed Plus Nebuliser

Minimum drug capacity: 2 ml

Maximum drug capacity: 8 ml

TECHNICAL SPECIFICATIONS OF DEVICE. Model: FLAEM JUNIOR F2000 combined RF7-2

Operating pressure (with neb.):	1.30 bar.	
Selector Key Position	With valve (Fig. 1)	Max (Fig. 2)
Delivery rate: [ml/min]	0.29 ⁽¹⁾	0.65 ⁽¹⁾
(Aerosol Output) AO: [ml]	0.430 ⁽²⁾	0.417 ⁽³⁾
Percentage of fill volume emitted per minute: [%]	4.64 ⁽²⁾	6.01 ⁽³⁾
(Aerosol output rate) AOR: [ml/min]	0.093 ⁽²⁾	0.120 ⁽³⁾
(Mass Median Aerodynamic diameter) MMAD: [µm]	2.61 ⁽²⁾	2.81 ⁽³⁾
(Geometrical Standard Deviation) GSD:	2.30 ⁽²⁾	2.29 ⁽³⁾
Respirable Fraction (< 5 µm): [%]	78.2 ⁽²⁾	75.7 ⁽³⁾
Fraction of aerosol particles < 2 µm: [%]	37.6 ⁽²⁾	34.5 ⁽³⁾
Fraction of aerosol particles: >2 µm ; < 5 µm: [%]	40.6 ⁽²⁾	41.2 ⁽³⁾
Fraction of aerosol particles > 5 µm: [%]	21.8 ⁽²⁾	24.3 ⁽³⁾
Residual volume: [ml]	0.90 ⁽²⁾	0.94 ⁽³⁾

(1) The delivery rate was measured with a 0,9% Na Cl saline solution at 23°C by weight loss according to Flaem internal procedure.

(2)(3) In vitro characterization according to EN ISO 27427 with 4ml, 0.1% Salbutamol in 0.9% NaCl solution with adult ventilatory pattern. Performance are likely to be different for paediatric or infant populations.

(2) Characterization performed by TÜV Rheinland Italia S.r.l. in collaboration with the University of Parma

(3) Characterization performed by University of Parma

Further details are available on request.

APPLIED PARTS

Type BF applied parts are: patient accessories (B3, B3.2, B4, B5)
Weight: 1.800 kg

ENVIRONMENTAL CONDITIONS

Operating conditions:

Ambient temperature: Between +10°C and +40°C

Relative air humidity: Between 10% and 95%

Atmospheric pressure: Between 69 KPa and 106 KPa

Storage and transport conditions:

Ambient temperature: Between -25°C and +70°C

Relative air humidity: Between 10% and 95%

Atmospheric pressure: Between 69 KPa and 106 KPa

DURATION

Model: JUNIOR F2000

(Compressor unit)

Service life 2000 hours.

Model: RF7-2

(Nebuliser and accessories)

The expected average life is 1 year, however it is advisable to replace the nebuliser every 6 months during intensive use (or sooner if the nebuliser is clogged) to ensure maximum therapeutic effectiveness.

EQUIPMENT AND MATERIAL INFORMATION		
The equipment includes:		Information on materials
A -	Compressor unit - Model: FLAEM JUNIOR F2000 A1 - Switch, A2 - Air intake, A3 - Air filter A4 - Nebuliser holder, A5 - Carrying handle A6 - Power cable, A7 - Cable compartment	
	Nebuliser and accessories - Model: RF7-2 B1 -Connecting tube (compressor/nebuliser unit) B2 - RF7 Dual Speed Plus Nebuliser B2.1 - Lower part, B2.2 - Nozzle B2.3 - Upper part, B2.4 - Speed selector with valve	
B -	B3 - Mouthpiece with valve B3.1 - Expiratory valve, B3.2 - Non-invasive nosepiece	Polypropylene
	B4 - Adult SoftTouch mask	Polypropylene +
	B5 - Paediatric SoftTouch mask	Thermoplastic Elastomers
	B6 - Manual nebulisation control	
C -	Convenient, roomy carrying bag	
IMPORTANT NOTE: There is an identification label on the packaging, remove it and affix it in the spaces provided on page 2. When replacing the Nebuliser and accessories, perform the same procedure.		

OPERATING INSTRUCTIONS

Before each use, wash your hands thoroughly and clean your appliance as described in the section "HYGIENE PREPARATION". During application, it is recommended to adequately protect oneself from drips. This appliance is suitable for the administration of medicinal substances (solutions and suspensions), for which aerosol administration is envisaged; such substances must in any case be prescribed by a doctor. In the case of substances that are too dense, dilution with suitable saline solution may be necessary, as prescribed by a doctor.

1. Insert the plug of the power cable (A6) into a mains socket corresponding to the voltage of the appliance. It must be positioned in such a way that disconnection from the power grid is not difficult.
2. Insert the nozzle (B2.2) into the upper part (B2.3) by pressing as indicated by the 2 arrows in the 'Connection diagram' in point B2. Insert the Speed selector with valve (B2.4) into the upper part (B2.3) as shown in the 'Connection diagram' in point B2. Pour the medicine prescribed by the doctor into the lower part (B2.1). Close the nebuliser by turning the upper part (B2.3) clockwise.
3. Connect the accessories as shown in the 'Connection Diagram' on the cover.
4. Sit comfortably while holding the nebuliser in your hand, place the mouthpiece to your mouth or use a nasal piece (if provided) or mask. If you use the mask accessory, place it on your face as shown in the (Fig. 1) (with or without the use of the elastic band).
5. Switch on the device by pressing the switch (A1) and inhale and exhale deeply; after inhaling, it is recommended to hold your breath for a moment so that the inhaled aerosol droplets can settle. Then

exhale slowly.

6. When the application is finished, switch off the device and disconnect it from the mains.

ATTENTION: If, after the therapy session, an obvious deposit of moisture forms inside the tube (B1), detach the tube from the nebuliser and dry it with the ventilation of the compressor itself; this action prevents possible mould growth inside the tube.

To facilitate the connection of the connection tube (B1) to the compressor unit, act on the end of it by simultaneously rotating and inserting it, and for removal by rotating and extracting it (Fig 2).


WAYS OF USING THE 'RF7 DUAL SPEED PLUS' NEBULISER WITH SPEED SELECTOR AND VALVE SYSTEM. It is professional, fast and suitable for the administration of all types of drugs, including the most expensive ones, even in patients with chronic diseases. Thanks to the geometries of the internal ducts of the RF7 Dual Speed Plus nebuliser, a suitable and active particle size was obtained for treatment down to the lower respiratory tract.

To speed up inhalation therapy, position the speed selector with valve (B2.4) by pressing the MAX marking with your finger (Fig. 4).



To make the inhalation therapy more effective, position the speed selector with valve (B2.4) by pressing with a finger on the side opposite the Max marking, in this case you will obtain an optimal administration of the drug while minimising its dispersion in the surrounding environment, thanks to the valve system with which the ampoule, mouthpiece and mask are equipped (Fig. 3).

USE OF NEBULIZATION MANUAL CONTROL. To achieve continuous nebulisation action you should not use the manual nebulisation control (B6), especially in the case of children or persons with reduced physical, sensory, or mental capabilities. The manual nebulisation control is useful for limiting dispersion of the medication in the surrounding environment. To start nebulizing close with a finger the hole of the nebulizer manual control (B6) and breathe in gently; we recommend to hold your breath for a moment so that the inhaled aerosol droplets can be deposited (Fig. 5). Meanwhile, to disable nebulizing, remove your finger from the hole of the nebulizer manual control to avoid the waste of drug, optimizing its acquisition. Then exhale slowly (Fig. 6).

SoftTouch masks

SoftTouch masks have an outer edge made of **soft, biocompatible material** that ensures an optimal fit on the face, and are also equipped with the **innovative Dispersion Limiter**. These characteristic features allow for greater sedimentation of the drug in the patient, and again limit its dispersion.

HYGIENIC PREPARATION.

Switch off the appliance before each hygienic preparation operation and disconnect it from the mains.

Compressor unit (A) and pipe exterior (B1)

Use only a cloth moistened with antibacterial detergent (non-abrasive and free of solvents of any kind).

Nebuliser and accessories. Open the nebuliser by turning the upper part (B2.3) anticlockwise, detach the nozzle (B2.2) and the speed selector with valve (B2.4) from the upper part (B2.3) as shown in the "Connection diagram" in B2. Then proceed according to the instructions below.

Sanitisation. Before and after each use, sanitise the nebuliser and accessories by choosing one of the methods provided in the table and described below.

method A: Sanitise the accessories under warm drinking water (approx. 40°C) with mild dishwashing detergent (non-abrasive).

method B: Sanitise the accessories in the dishwasher with the hot cycle (70°C).

method C: Sanitise the accessories by soaking in a solution of 50% water and 50% white vinegar, then rinse thoroughly with warm drinking water (approx. 40°C).

After sanitising the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Disinfection. After sanitising the nebuliser and accessories, disinfect them using one of the methods provided in the table and described below. Each method can be performed for a limited number of times (see figure in table).

method A: Obtain a disinfectant of the electrolytic chloride type (active ingredient: sodium hypochlorite), specifically for disinfection, available in all pharmacies.

Execution: - Fill a container of a suitable size to hold all the individual components to be disinfected with a solution of drinking water and disinfectant, observing the proportions indicated on the packaging of the disinfectant. - Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution. - Recover the disinfected components and rinse them thoroughly with lukewarm drinking water. - Dispose of the solution according to the disinfectant manufacturer's instructions.

method B: Disinfect accessories by boiling them in water for 10 minutes; use demineralised or distilled water to avoid hard water deposits.

method C: Disinfect the accessories in a hot bottle steamer (not microwave). Carry out the process strictly according to the hot bottle steamer's instructions. For disinfection to be effective, choose a steriliser with an operating cycle of at least 6 minutes.

If you also want to carry out sterilisation, skip to the **Sterilisation** section.

After disinfecting the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

At the end of each use, store the device complete with its accessories in a dry and dust-free place.

Sterilisation Device: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060. **Execution:** Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Place the packaged components in the steam steriliser. Perform the sterilisation cycle in accordance with the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes first.

Storage: Store sterilised components according to the instructions for use of the system or sterile barrier packaging chosen.

The sterilisation procedure was validated in accordance with ISO 17665-1.

At the end of each use, store the device complete with its accessories in a dry and dust-free place.

Table of planned methods / patient accessories									
Method	B2.1	B2.2	B2.3	B2.4	B3	B3.2	B4	B5	B6
HYGIENIC PREPARATION IN THE HOME									
Sanitisation									
method A	✓	✓	✓	✓	✓	✓	✓	✓	✓
method B	✓	✓	✓	✓	✓	✓	✓	✓	✓
method C	✓	✓	✓	✓	✓	✓	✓	✓	✓
Disinfection									
method A	300	300	300	300	300	300	300	300	300
method B	300	300	300	300	300	300	300	300	300
method C	300	300	300	300	300	300	300	300	300
HYGIENIC PREPARATION IN A CLINIC OR HOSPITAL									
Disinfection									
method A	300	300	300	300	300	300	300	300	300
method B	300	300	300	300	300	300	300	300	300
Sterilisation									
	300 Assembled nebuliser				300	300	300	300	300
✓: planned ✗: not planned 300: ✓ MAX 300 TIMES									

AIR FILTRATION. The appliance is equipped with an extraction filter (A3) to be replaced with a replacement filter (D1) when it is dirty or changes colour. Do not wash or re-use the same filter. Regular replacement of the filter is necessary to help ensure proper compressor performance. The filter must be checked on a regular basis. Contact your dealer or authorised service centre for replacement filters.

To replace the filter, pull it out as shown in the (Fig. 7)

The filter is designed so that it is always fitted in its housing. Do not replace the filter during use. **Only use original Flaem accessories or spare parts, no liability is accepted if non-original parts or accessories are used.**

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FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.