
EASYRED

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VENDUTO A _____

VIA _____ LOCALITÀ _____

OPERATION AND MAINTENANCE INSTRUCTION MANUAL

CONTENTS:

| | |
|------------------------------------|----|
| INTRODUCTION | 16 |
| WARNINGS & CAUTIONS | 17 |
| PRINCIPLES OF OPERATION | 18 |
| GENERAL INSTRUCTIONS FOR OPERATION | 20 |
| CLEANING AND STERILIZATION | 20 |
| TESTING THE RESUSCITATOR | 21 |
| MATERIAL | 23 |
| PERFORMANCE CHARACTERISTICS | 24 |
| WARRANTY | 26 |

INTRODUCTIONS

Thank you for having chosen the manual resuscitator the EASYRED line by MORETTI. Its design and qualities guarantee the maximum comfort, the maximum safety and reliability. The manual resuscitator the EASYRED line by MORETTI have been designed and made to satisfy all your demands. This handbook gives you some little suggestions for a correct use of the device you have chosen and some useful advice for your safety. We suggest you to read carefully the whole handbook before using your kit. In case of any doubt, please contact your retailer, who will be able to help you and give you some correct advice.

The is designed for use as an adjunct to artificial respiration and cardio-pulmonary resuscitation. The resuscitator can be used to ventilate the aphonic patient and to augment ventilation and / or oxygen delivery to the spontaneously breathing patient.

When connected to an oxygen source, the manual resuscitator can also deliver supplemental oxygen with the aid of an oxygen reservoir. The oxygen concentration delivered by the manual resuscitator delivered by the manual resuscitator depends on the oxygen supply flow rate, tidal volume, and ventilation frequency and operator technique.

Adult manual resuscitators are supplied with over pressure are supplied with over pressure valve 60cmH₂O. The Infant and child models of the manual resuscitator are fitted with pressure relief valves which vent the delivered breath to the atmosphere if the pressure exceeds 40cm H₂O. The pressure relief valve can be overridden by the operator should a higher delivery pressure be required.

The manual resuscitator is intended for use by qualified medical and emergency personnel trained in pulmonary ventilation and advanced cardiac life support techniques. Users should read and understand the contents of this manual and demonstrate proficiency in the assembly, disassembly and use of this device prior to use.

WARNINGS AND CAUTIONS

Warning: do not use the manual resuscitator in toxic atmosphere.

Warning: remove the oxygen reservoir and reservoir valve if supplemental oxygen is not being administered. Failure to do so will affect the refill rate and maximum frequency capabilities.

Warning: do not administer supplemental oxygen in the presence of open flame.

Warning: do not use oil, grease or other hydrocarbon-based substance on any part of the manual resuscitator. Supplemental oxygen, supplied under pressure, can combine with hydrocarbons and cause explosions.

Warning: the device is intended for use by qualified medical and emergency personnel trained in pulmonary ventilation and advanced cardiac life support techniques.

Warning: always test the device in accordance with this manual after leaning and sterilization or replacement of parts.

Warning: do not use the pressure relief valve over-rid mechanism for the child and infant resuscitator.

Caution: use only replacement parts intended for use on manual resuscitators.

Caution: do not attempt to disassemble the pressure relief valve assembly. Disassembly will damage the component.

PRINCIPLES OF OPERATION 1

The manual resuscitator consists of four major components (figure 1): non-rebreathing valve assembly (A), silicone bag (B), reservoir valve (C) and oxygen reservoir (D). **Parts C and D should be removed if supplemental oxygen is not to be supplied from an external gas source.**


Figure 1. The Manual Resuscitator

Gas is delivered to the patient by squeezing the silicone bag (see figure 2). Positive pressure within the bag caused by its compression closes the intake valve (E) located at the base of the bag. This causes the 2 duckbill valve (F) to close off the expiration ports (G) of the non-rebreathing valve assembly. Further compression of the bag forces gas through the duckbill valve to the patient. If supplemental gas is used, oxygen is delivered to the reservoir (D) during compression of the silicone bag.


Figure 2. Inspiration

PRINCIPLES OF OPERATION 2

Exhalation begins when the patient exerts positive pressure (during passive exhalation) on the patient side of the duckbill valve (F) or when the operator releases pressure from the silicone bag (B). The valve lifts directing the patient's expiratory gases through the expiration ports (G) of the non-rebreathing valve (see figure 3).


Figure 3. Exhalation

The silicone bag refills for the next breath during exhalation. Negative pressure within the bag (caused by the expansion of the compressed bag) opens the intake valve (E), allowing gas to enter the bag either from the atmosphere or from the oxygen reservoir (D). The reservoir should be used whenever supplemental oxygen is delivered. (Supplemental oxygen may be administered without using the reservoir but the maximum available oxygen concentration will be reduced.) Excess oxygen vents to the atmosphere through the safety valve (H) on the reservoir (C) should the reservoir fill before the next delivered breath. If the volume of gas in the reservoir is inadequate to fill the silicone bag, room air may be drawn in through the reservoir valve safety inlet (I). The concentration of the oxygen-enriched gas entering from the reservoir will depend on factors such as oxygen flow rate, tidal volume, ventilation rate and operator technique.

PRINCIPLES OF OPERATION 3

The infant and Child models of the Manual Resuscitator are fitted with pressure release valves (see figure 4). These valves open when the pressure is in excess of 40+5 cmH₂O.


Figure 4. Pressure Relief Valve

Should higher inspiratory pressure be required, the pressure relief valve may be overridden by placing the thumb over the valve as shown in figure 5. It is also possible to lock the pressure relief valve in this position with air valve bottom. Lock Type Non-rebreathing Valve. Just twist the valve a quarter turn and it will stay in the overridden position. Arrow on the label forward to the NRV direction.


Figure 5. Overriding the Pressure Valve

GENERAL INSTRUCTIONS FOR OPERATION

1. Select the appropriate size of Manual Resuscitator (infant, child or adult). Select the appropriate size of face mask if ventilation by mask is to be performed. Connect mask to outlet on non-rebreathing valve housing.
2. If supplemental oxygen is to be used, connect supply tubing (not supplied) between the regulated gas source and the oxygen reservoir inlet. Adjust gas flow so that the reservoir expands completely during the delivery of the breath and nearly collapses as the silicone bag refills. If supplemental oxygen is not to be used, remove the reservoir and reservoir valve. Failure to do so will affect the refill rate and maximum frequency capabilities.
3. Should the non-rebreathing valve become contaminated with vomit, blood or secretions during ventilation, disconnect the Resuscitator from the patient and clear the non-rebreathing valve as follows:

Squeeze the silicone bag to deliver several sharp breaths through the non-rebreathing valve to expel the contamination. If the contamination does not clear, continue ventilation with another resuscitator or use mouth to mouth/mask techniques. Disassemble the non rebreathing valve and rinse in water. Reassemble and sterilize.
4. Once the procedure is complete, clean, sterilize and test the Resuscitator as described under Cleaning & Sterilization and Testing the Resuscitator.

CLEANING AND STERILIZATION

The Manual Resuscitator should be cleaned and sterilized as follows:

- when first using the new Resuscitator
- between patients
- whenever the Resuscitator becomes contaminated
- every 24 hours of use with the same patient.

Carry out cleaning and sterilization as follows:

1. Disassemble the Resuscitator, and, if necessary, the face mask.

See figure 6. Caution: Do not attempt to disassemble the pressure relief valve assembly as it can be cleaned without disassembly. Disassembly will damage the component.

gent. Ensure that the detergent is compatible with the component materials (see Specifications).

3. Rinse the components thoroughly in warm water to remove all traces of detergent.
4. Sterilize the components using one of the following methods:
 - **all parts can be autoclaved (max: 132°C) except the oxygen reservoir.**
 - all parts can be sterilized using ethylene oxide.
 - all parts can be sterilized using most common disinfecting solutions which are used for equipment coming into contact with the patient. Rinse thoroughly with water after sterilization.
5. Dry all components thoroughly after sterilization.
6. Inspect all components for wear. Replace as necessary.
7. Reassemble the Resuscitator as shown in figure 6.
8. Before use, test the Resuscitator as described in Testing the Resuscitator.

TESTING THE RESUSCITATOR

The manual resuscitator should be tested as follows:

- when first using the new resuscitator
- after cleaning and sterilizing.
- after any new parts have been fitted
- monthly, if the resuscitator is not frequently used

TESTING 1 - TESTING THE SILICONE BAG ASSEMBLY

1. Remove the non-rebreathing valve and the oxygen reservoir and valve (if fitted).
2. Compress the silicone bag and occlude (block) the non-rebreathing valve outlet.
3. Release the bag. The bag should expand immediately and refill. If not, check that the intake valve at the base of the silicone bag is correctly assembled.
4. While keeping the non-rebreathing valve outlet blocked, compress the bag again. The bag should not compress easily. If this occurs, check that you are blocking the valve sufficiently, and that the intake valve at the base of the silicone bag is correctly assembled.

NON-REBREATHING VALVE ASSEMBLY

1. Connect the non-rebreathing valve to the silicone bag. Connect the test lung to the outlet on the non-rebreathing valve.

2. Compress and hold the bag. The non-rebreathing (duckbill) valve inside the non-rebreathing valve should open and the test lung should fill. If not, check the connection between the Resuscitator and the test lung, and check that the non-rebreathing valve is correctly assembled.

3. Release the bag. The non-rebreathing (duckbill) valve should close and as the test lung deflates, gas should flow through the expiratory ports in the non-rebreathing valve. If not, check that the non-rebreathing valve is correctly assembled.

4. Ventilate the test lung for a minimum of ten cycles to ensure that the Resuscitator is functioning correctly. Inspiration must occur when the silicone bag is compressed and exhalation when the bag is released. If not, check that the non-rebreathing valve is correctly assembled. To check the function of the pressure relief valve (Infant and Child Resuscitators).

TESTING 2 OXYGEN RESERVOIR/RESERVOIR VALVE ASSEMBLIES

1. Attach the reservoir to the reservoir valve assembly. Attach the silicone bag.

2. Inflate the reservoir and block the reservoir port.

3. Compress the reservoir bag. Gas should escape through the safety outlet valve on the reservoir valve. If not, check that the reservoir valve is correctly assembled.

4. Connect the reservoir and reservoir valve assembly to a Resuscitator.

5. Cycle the Resuscitator through several ventilations. The safety inlet valve on the reservoir valve should open during each refill to allow room air to enter the silicone bag. If not, check that the reservoir valve is correctly assembled.

Note: If supplemental oxygen is not connected, the silicone bag will refill more slowly if the reservoir is still attached.

TESTING 3 OVERALL RESUSCITATOR FUNCTION

1. Fully assemble the Resuscitator (non-rebreathing valve, silicone bag, reservoir valve and oxygen reservoir). Connect the Resuscitator to a supplemental gas source and connect a test lung to the patient outlet on the non-rebreathing valve.

2. Set the supplemental gas flow to 15 litres per minute for the adult and child models; and to 10 l/min for the infant model.

3. Cycle the Resuscitator through several ventilations. The test lung

should inflate during inspiration and deflate during exhalation. Check for leakage at all joints and connections. Ensure that the Resuscitator refills promptly and properly and that all valves are operating correctly. If not, repeat the tests above to find where the problem lies.

Storage temp.

-40°C/10°F to 60°C/140°F

Operating temp.

-18°C/0°F to 50°C/122°F

MATERIAL

Silicone rubber

adult mask bladder, duckbill valve, flapper valve, mask grommet.

Mask retainers, relief valve seal, silicone bag.

Polysulphone

adult mask shell, bag intake valve housing, non-rebreathing valve, housing, pressure relief valve housing, pressure relief valve, stem, reservoir bag connector, reservoir valve housing, oxygen reservoir bag, pressure relief valve spring.

Polyvinyl chloride

Stainless steel

Pressure Relief

40+- 5 cmH₂O (infant & child)

Reservoir volume

60+- 5 cmH₂O (adult)
2500ml adult and child
600ml infant model

Bag volume

(stroke volume)

1500/1350ml (adult)
550/350ml (child)
280/100ml (infant)

Maximum BPM

45 breaths/min. (adult model)
105 breaths/min. (child model)
98 breaths/min. (infant model)

Oxygen Concentration

with reservoir

99%

without reservoir

45% (adult and child model)
90% (infant model)

Range of delivery pressure

adult : unrestricted of

60+-10 cmH₂O

infant & child

40+- 5 cmH₂O

(can be overridden by operator)

Stroke volume range (typical values)

| | | |
|--------|-----------------|--------|
| Adult | using one hand | 800ml |
| Adult | using two hands | 1350ml |
| Child | | 350ml |
| Infant | | 100ml |

PERFORMANCE CHARACTERISTICS

The performance characteristics for Manual Resuscitators will vary from user to user depending on a variety of factors: ambient temperature, patient lung compliance, ventilator frequency, size of operator's hands. The following data has been derived from Manual Resuscitators in accordance with ISO standards.

Temperature

-18°C/0°F 22°C/72°F 50°C/122°F

Cycle rate

| | | | |
|----------|----|-----|-----|
| Adult | 38 | 45 | 45 |
| Child A | 92 | 105 | 105 |
| Child B | 72 | 78 | 76 |
| Infant A | 95 | 95 | 95 |
| Infant B | 60 | 60 | 60 |

These results were obtained under the following conditions:

| | | | |
|----------|--------------------|---|---|
| Adult | $V_T=600\text{ml}$ | compliance -0,02 1/cmH ₂ O, | resistance -20 cmH ₂ O(1/S) |
| Child A | $V_T=70\text{ml}$ | compliance -0,01 1/cmH ₂ O, | resistance -20 cmH ₂ O(1/S) |
| Child B | $V_T=300\text{ml}$ | compliance -0,02 1/cmH ₂ O, | resistance -20 cmH ₂ O(1/S) |
| Infant A | $V_T=20\text{ml}$ | compliance -0,001 1/cmH ₂ O, | resistance -400 cmH ₂ O(1/S) |
| Infant B | $V_T=70\text{ml}$ | compliance -0,01 1/cmH ₂ O, | resistance -400 cmH ₂ O(1/S) |

PERFORMANCE CHARACTERISTICS**ADULT**

Values in parentheses are without an oxygen reservoir
O₂ flow in Tidal vol. (ml) x Ventilation rate

| | | | | | | |
|--------|--------|--------|--------|--------|---------|---------|
| 1/min. | 600x12 | 600x20 | 750x20 | 750x20 | 1000x12 | 1000x20 |
| 05 | 82(33) | 58(33) | 65(33) | 99(30) | 99(30) | 85(33) |
| 10 | 99(37) | 80(37) | 99(37) | 99(37) | 99(37) | 99(37) |
| 15 | 99(45) | 99(45) | 99(45) | 99(45) | 99(45) | 85(45) |

Conditions:

compliance -0,02 1/cmH₂O, Resistance -20 cmH₂O/(1/S) I:E ratio 1:2**CHILD**

Values in parentheses are without an oxygen reservoir
O₂ flow in Tidal vol. (ml) x Ventilation rate

| | | | | | | |
|--------|--------|--------|--------|--------|--------|--------|
| 1/min. | 70x20 | 70x30 | 200x20 | 200x30 | 300x20 | 300x30 |
| 02 | 99(65) | 99(65) | 58(37) | 50(37) | 45(34) | 40(32) |
| 05 | 99(86) | 99(83) | 99(44) | 96(43) | 84(37) | 68(36) |
| 10 | 99(93) | 99(90) | 99(46) | 99(46) | 99(45) | 99(43) |

Conditions:



 $V_T >150 <600\text{ml}$ compliance -0,02 1/cmH₂O,Resistance -20 cmH₂O/(1/S) I:E ratio: 1:2**INFANT**

Values in parentheses are without an oxygen reservoir.

O₂ flow in Tidal vol. (ml) x Ventilation rate

| | | | | | | |
|--------|--------|--------|--------|--------|--------|--------|
| 1/min. | 20x30 | 20x60 | 40x30 | 40x60 | 70x30 | 70x60 |
| 02 | 99(72) | 99(70) | 99(70) | 99(58) | 99(55) | 85(50) |
| 05 | 99(77) | 99(75) | 99(74) | 99(72) | 99(65) | 85(62) |
| 10 | 99(94) | 99(90) | 99(90) | 99(80) | 99(79) | 99(71) |

Conditions:

 $V_T =20\text{ml}$ Compliance -0,001 1/cmH₂O,Resistance -400 cmH₂O/(1/S) I:E ratio: 1:1 $V_T =40\text{ml} \& 70\text{ml}$ Compliance -0,01 1/cmH₂O,Resistance -20 cmH₂O/(1/S) I:E 1:2**Inspiratory flow resistance****Expiratory flow resistance**

WARRANTY

This device has 2-year warranty from the purchase date which is certified by the receipt or invoice. The distributor reserves the right either to repair or replace necessary parts to correct defects in the materials or workmanship without any charge.

The warranty does not apply to damage resulting from failure to follow the operating instructions, accidents, abuse, alteration or disassembly by unauthorized individuals.

In the case you need to get your device repaired and you fulfil the conditions of warranty, send the device to the nearest retailer together with the purchase receipt.

WARRANTY CERTIFICATE

DEVICE MODEL _____

PURCHASING DATE _____

AUTHORIZED DEALER _____

STREET _____ PLACE _____

PURCHASED BY _____

STREET _____ PLACE _____