**NEW ASKIR 230/12V BR**

**CE 0123**

**ASPIRATORE NEW ASKIR 230/12V BR  
ASPIRATOR NEW ASKIR 230/12V BR  
ASPIRATEUR NEW ASKIR 230/12V BR  
ABSAUGER NEW ASKIR 230/12V BR  
ASPIRADOR NEW ASKIR 230/12V BR  
АСПИРАТОР NEW ASKIR 230/12V BR  
ASPIRATÖR NEW ASKIR 230/12V BR**

APPROVED BY

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**NEW ASKIR 230/12V BR** Surgical Aspirator is a portable unit designed for the aspiration of bodily fluids in adults and children. It is particularly suitable for nasal, oral and tracheal aspiration of mucus, catarrh or blood after minor surgical procedures and can be used in post-operative therapy at home or conveniently transported from one hospital ward to another. The device operates on a 230V ~ AC 12V DC mixed power supply and is designed for virtually continuous use. This is made possible by an advanced electronic system that manages the power supply and enables the unit to operate even while the battery is being charged. The LED lights on the front control panel signal when the unit is operating, as well as the charge and re-charge status of the batteries. The New Askir 230/12V Br comes equipped an aspiration regulator on the front panel and a polycarbonate autocavable jar complete with overflow valve. It's made of a highly heat-resistant electrically insulated plastic material that conforms with the latest European safety standards.

## GENERAL WARNING

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE**

**ONLY HIGHLY QUALIFIED STAFF USE RESERVED**

**THE INSTRUMENT MUST NOT BE DISASSEMBLED  
FOR A TECHNICAL SERVICE ALWAYS CONTACT CA-MI**

## IMPORTANT SAFETY RULES

1. Check the condition of the unit before each use. The surface of the unit should be carefully inspected for visual damage. Check the mains cable and **don't connect to power** if damage is present.
2. Before connecting to power, always check the manufacture's label to ensure the electrical design of the unit and style of plug comply with the electrical network to be used.
3. If the plug supplied with the appliance is incompatible with the mains electricity socket, contact qualified staff for replacement of the plug with a suitable type. The use of simple or multiple and / or extension adapters is not generally recommended. Whenever their use is indispensable, use those in compliance with safety regulations, however paying attention not to exceed the maximum power supply limits, which are indicated on the adapters and extensions.
4. Pay special attention to the following:
  - Only use original accessories and components;
  - The device can be used only with the bacteriological filter;
  - Never immerge the appliance into water;
  - Position the appliance on flat stable surfaces;
  - Position the device in a way that the air inlets on the back aren't obstructed;
  - Never use the device in environments which have anaesthetic mixtures inflammable with air, oxygen or nitric oxide;
  - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
  - Keep off the reach of children or not capable people without supervision;
  - Don't leave the appliance connected to the power supply socket when not in use;
  - Hold plug to remove from socket. Don't pull on cable
  - Preserve and use the medical device in environments protected from atmospheric factors and at a distance from heat sources.
  - Don't use the device thoracic drainage.
5. For repairs, exclusively contact CA-MI technical service and request the use of original spare parts. Failure to comply with the above can jeopardise the safety of the device.
6. **This medical device must be destined exclusively for the use for witch it has been designed ad described in this manual.** Don't use for anything other than the use defined by the manufacturer. The manufacturer will not be responsible for damage due to improper use or connection to an electrical system not complying with current regulations.
7. The device is intended for use in the electromagnetic environment outlined in this manual. Install and use as directed.
8. Instrument and accessories discharging must be done following current law regulations in every country of use.
9. The lead battery integrated in the device is not to be considered as an ordinary domestic waste. Such a component must be disposed of in a specific collection centre in order to be recycled.
10. None of electric or mechanical parts have been designed to be repaired by customers or end-users. Don't open the device, do not mishandle the electric / mechanical parts. Always contact technical assistance

**The manufactured cannot be held liable for accidental or indirect damages should the device be modified, repaired without authorization or should any of its component be damaged due to accident or misuse.**

**Any minimal modification / repair on the device voids the warranty and does not guarantee the compliance with the technical requirements provided by the MDD 93/42/EEC (and subsequent changes) and its normative.**

### DISPOSAL OF WASTE BATTERIES - (Directive 2006/66/EC)

This symbol on the battery or on the packaging indicates that the battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery.


The recycling of the materials will help to conserve natural resources. At the end of their life hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

## TECHNICAL CHARACTERISTICS

TYPOLGY (MDD 93/42/EEC)	Class IIa Medical Device
MODEL	<b>NEW ASKIR 230/12V BR</b>
UNI EN ISO 10079-1	HIGH VACUUM / LOW FLOW
POWER FEEDING	230V~ / 50Hz or Internally powered equipment (PB Battery 12V 3Ah) or 12Vdc (1.9A) from car battery
POWER CONSUMPTION	50VA
FUSE	F 1 x 630mA 250V
MAXIMUM SUCTION PRESSURE (without jar)	-75kPa (-0.75 Bar)
MINIMUM SUCTION PRESSURE (without jar)	Less - 40kPa (-0.40 bar)
MAXIMUM SUCTION FLOW (without jar)	16 l/min
WEIGHT	4.9 Kg
SIZE	350 x 210 x 180 mm
BATTERY HOLDING TIME	45 ÷ 50 minutes
BATTERY TIME CHARGE	120 ÷ 150 minutes
WORKING CONDITION	Room temperature: 5 ÷ 35°C Room humidity percentage: 30 ÷ 75% RH Altitude: 0 ÷ 2000m s.l.m.
CONSERVATION CONDITION AND TRASPORT	Room temperature: -40÷ 70°C Room humidity percentage: 10 ÷ 100% RH

## SYMBOLS

	Class II isolation equipment
<b>CE 0123</b>	CE marking in conformity with EC directive 93/42/EEC and subsequent changes
	Warning, consult the instruction manual
	To Preserve in place coolness and dry land
	Conservation temperature: -40 ÷ 70°C
	Type B equipment
	Fuse
~	Alternate Current
	Battery
	Direct Current
Hz	Mains Frequency
I	ON
O	OFF
	DEHP Phthalates (Suction catheter)

## ACCESSORIES SUPPLIED

DESCRIPTION
COMPLETE ASPIRATION JAR 1000cc
CONICAL FITTING
TUBES SET 6 mm x 10 mm
ASPIRATION PROBE CH20
ANTIBACTERIAL FILTER

The filter is produced with (PTFE) hydrophobic material to prevent fluids entering the pneumatic circuit. It should be changed immediately if it becomes wet or if there is any sign of contamination or discolouration. It should also be changed if the unit is used with a patient whose risk of contamination is unknown. **Don't use the suction unit without the protection filter.** If the suction unit is used in an emergency or in a patient where the risk of contamination is not known the filter must be changed after each use. Available under request with different versions with complete jar 2000cc.

**Suction catheter:** Single-use device to be used on a single patient. Do not wash or re-sterilize after use.

Reuse may cause cross-infections. Don't use after lapse of the sell-by date

**WARNING:** Suction tubes for insertion in the human body purchased separately from the machine should comply with ISO 10993-1 standards on material biocompatibility.

Guidance and manufacturer's declaration – Electromagnetic Emissions		
The surgical aspirator <b>NEW ASKIR 230/12V BR</b> is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator <b>NEW ASKIR 230/12V BR</b> should assure that it's used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment - guidance
Irradiated / Conducted emissions CISPR11	Group 1	The surgical aspirator <b>NEW ASKIR 230/12V BR</b> only used RF energy only for its internal functioning. Therefore its RF emissions are very low and are not cause interference in proximity of any Electronic appliances.
Irradiated / Conducted emissions CISPR11	Class [B]	The surgical aspirator <b>NEW ASKIR 230/12V BR</b> can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes or environments feeds to you from batteries.
Harmonic emissions IEC/EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

Guidance and manufacturer's declaration – Electromagnetic Emissions		
The surgical aspirator <b>NEW ASKIR 230/12V BR</b> is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator <b>NEW ASKIR 230/12V BR</b> should assure that it's used in such an environment.		
Immunity Test	Compliance	Electromagnetic environments - guidance
Electrostatic discharge (ESD) IEC/EN 61000-4-2	± 6kV on contact ± 8kV in air	Floors should be wood, concret or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN 61000-4-4	± 2kV power supply	Mains power quality should be that of a typical commercial environment or hospital or environments feeds to you from batteries.
Surge IEC/EN 61000-4-5	± 1kV differential mode	Mains power quality should be that of a typical commercial environment or hospital or environments feeds to you from batteries.
Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11	5%U <sub>T</sub> for 0.5 cycle 40%U <sub>T</sub> for 05 cycle 70%U <sub>T</sub> for 25 cycle <5%U <sub>T</sub> for 5 sec	Mains power quality should be that of a typical commercial environment or hospital If the user of the surgical aspirator <b>NEW ASKIR 230/12V BR</b> request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field IEC/EN 61000-4-8	3A/m	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Conducted Immunity IEC/EN 61000-4-6	3Vrms 150kHz to 80MHz (for appliances that aren't life - supporting)	-
Irradiated Conducted IEC/EN 61000-4-3	3V/m 80MHz to 2.5 GHz (for appliances that aren't life - equipment)	-
Note U <sub>T</sub> is the value of the power supply voltage		

## CLEANING OF ACCESSORIES

To clean the plastic housing of the device wear disposable latex gloves and clean with denaturated alcohol or hypochlorite solutions.

Washing and / or cleaning the autoclavable jar as to be carried out as follows:

- Wear protection gloves and apron (glasses and face mask if necessary) to avoid contact with contaminating substances;
- Disconnect the jar from the device
- Disconnect all tubes from the jar and the protection filter
- Empty and dispose of the content and of the suction catheter according to the laws in force in your country;
- Separate all parts of the cover (overflow valve, o-ring);

After disposing of disposable parts and disassembling the jar wash in running cold water and rinse thoroughly.

Then soak in warm water (temperature shall not exceed 60°C). Wash thoroughly and if necessary use a non-abrasive brush to remove incrustations.

Rinse in running warm water and dry all parts with a soft cloth (non-abrasive).

The jar and the cover can be autoclaved by placing the parts into the autoclave and running one sterilization stem cycle at 121°C (1 bar relative pressure) making sure that the jar is positioned upsidedown.

Mechanical resistance of the jar is guaranteed up to 30 cycles of sterilization and cleaning at the indicated conditions (EN ISO 10079-1).

Beyond this limit the physical-mechanical characteristics of the plastic may decrease and replacement of the part is therefore recommended.

After sterilization and cooling at environment temperature of the parts make sure that these are not damaged.

Assemble the jar as follows:

- Place the overflow valve into its seat in the cover (under VACUUM connector)
- Insert floating valve keeping the o-ring towards the opening of the cage
- Place the o-ring into its seat around the cover
- After completing assembling operations always make sue that cover seals perfectly to avoid vacuum leakages or liquid exit

The aspiration tubes can be sterilized on autoclave using a sterilization cycle at 120°C. The conical connector can be sterilized on autoclave using a sterilization cycle at 121°C. The device is ready for a new employment now.

**DO NOT WASH, STERILIZE OR PUT IN AUTOCLAVE THE ANTIBACTERIAL FILTER**

## RULES FOR RETURNING AND REPAIRING

**COMPLYING WITH THE NEW EUROPEAN RULES, CA-MI INDICATES THE IMPORTANT POINTS TO PROTECT INSTRUMENT AND OPERATORS HYGIENE. THESE RULES MUST BE RESPECTED IN ORDER TO GUARANTEE HYGIENE AND SAFETY TO ALL THE PEOPLE OPERATING WITH THE INSTRUMENT TO OBTAIN QUALITY AND WELL BEING.**

Before returning an instrument for repair, the external surfaces and all accessories **MUST** be carefully disinfected with a cloth soaked in methylated spirits or hypochlorite-based solution. If CA-MI finds instrument not suitable for repairing due to clear signs of internal or external contamination, the same will be returned to customer with specification of **NOT REPAIRED INSTRUMENT**, accompanied by an explanation letter. CA-MI will decide if contamination is due to bad functioning or misuse. If contamination is due to bad functioning, CA-MI will substitute the instrument, only if a **SALE RECEIPT** and **STAMPED GUARANTEE** accompany the same.

CA-MI is not responsible for contaminated accessories, they will be substitute at customer's expenses.

For this reason it is **COMPULSORY** to carefully disinfect the external part of the instrument and accessories with a cloth soaked in methylated spirits or hypochlorite-based solutions. Put the instrument and accessories in a bag with indication of disinfecting.

We also request to specify the kind of fault, in order to speed up repairing procedures.

To this end, please read the instructions carefully in order to avoid damaging the equipment through improper use.

Always specify the fault encountered so that CA-MI can establish whether it falls into the category of the faults covered by the guarantee.

**CA-MI Srl cannot be held liable for accidental or indirect damages should the device be modified, repaired without authorization or should any of its component be damaged due to accident or misuse.**

**Any minimal modification / repair on the device voids the warranty and does not guarantee the compliance with the technical requirements provided by the MDD 93/42/EEC Directive and its normatives.**

## MAINTENANCE

The **NEW ASKIR 230/12V BR** suction equipment does not need maintenance or lubrication.

It is, however to inspect the unit before each use. Unpack the instrument and **always check** integrity of plastic parts and feeding cable, they might have been damaged during previous use. Connect cable to electrical network and turn switch on.

Close the aspirator outlet with your finger and with suction regulator in maximum vacuum position check that the vacuum indicators reaches – 75 kPa (-0.75 bar) minimum (internal battery). Rotate the knob from right to left and check the aspiration regulating control.

The vacuum indicator should go down – 40 kPa (-0.40 bar). Verify that loud noises are not present, these can indicate wrong functioning.

The equipment may also be used with the power cord hooked up to the electricity mains (the yellow light switches over to the fixed green light).

Ensure that the green light comes on and that, after a few seconds it switches over to the yellow light (meaning that the charged battery is operating).

A protection fuse (**F 1 x 630mA 250V**) reachable from exterior and it situated in the plug protects the instrument.

For fuse replacing, always the type and the range. The device is made up of a lead battery which cannot be accessed by outside. In order to replace it, consult the technical staff authorised by the manufacturer.

Fault type	Cause	Solution
1. Red light on	Battery run down	Hook up the power cord to the electricity mains, positioning the equipment power switch on 0. Charge it for 8 hours, until the yellow light switch over to the green light.
2. Red light on	Battery run down	The device may be also be supplied by an power cord hooked up to the electricity mains: the yellow light comes on and then switch over to the green light.
3. No aspiration	Jar Cap badly screwed down	Unscrewed the cap, then rescrew it correctly
4. No aspiration	Lid seal not in its seat	Unscrew the cap and insert the seal properly in its seat
5. The float doesn't close	If the cap has been washed, ensure that the float is not partially detached	Insert the float into it's place
6. The float doesn't close	The float it's covered by dirty material	Unscrewed the cap, leave the and put in on autoclave
7. Low suction	Foam inside the jar	Fill the jar to 1/3 full of ordinary water
8. No aspiration due to flow leakage of mucus	Filter blocked	Replace filter
9. The Vacuum power on the patient side is either very low or absent	<ul style="list-style-type: none"> <li>• Vacuum regulator set to minimum</li> <li>• Protection filter blocked or damaged</li> <li>• Connection tubes blocked, kinked or disconnected</li> <li>• Shut-off valve blocked or damaged</li> <li>• Pump motor damaged</li> </ul>	<ul style="list-style-type: none"> <li>• Turn the vacuum regulator clockwise and check the value of the vacuum on the gauge</li> <li>• Replace the filter</li> <li>• Replace or reconnect the tubes, check the jar connections</li> <li>• Empty the jar, or disconnect the tube from the jar and unblock the shut-off valve. The unit will only work in the upright position</li> <li>• Refer to authorised service personnel</li> </ul>
<b>Faults 1 - 2 - 3 - 4 - 5 - 6 - 7 - 8 - 9</b>	<b>None of the remedies has achieved the desired results</b>	<b>Contact the seller or CA-MI After-sales Assistance Service</b>

If the overflow security system it's activated, don't proceed with the liquid aspiration.

If the overflow security system doesn't work there are two cases:

1° case – If the overflow security system doesn't work the aspiration will be stopped by the bacteriological filter who avoid the liquid penetration inside the device.

2° case – If both the security system doesn't work, there is the possibility that liquid comes inside the device, in this case return the device to CA-MI technical service.

**CA-MI Srl will provide upon request electric diagrams, components list, descriptions, setting instructions and any other information that can help the technical assistance staff for product repair.**

**BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT CA-MI TECHNICAL SERVICE.**

**CA-MI DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.**

## INSTRUCTIONS

- Place the unit on a flat, horizontal surface.
- Connect the short silicon tube 3, with antibacterial filter 1, to the suction connector 4. The other tube, with one end connected to the antibacterial filter with the other end to jar's lid connector 2 where has been fixed the red float. When the 90% of the volume of the jar is reached there is the activation of the security float (the float close the aspiration connector on the jar) to the avoid liquid penetration inside the device.
- Connect the long silicon tube 5 to the other jar's lid connector 6
- Connect the other end of the long silicon tube 5 to the probe plastic connector 8 then connect the suction probe to it.
- Connect the power cord to the device then connect the plug to the electrical mains supply.
- Push switch 7 on position I to start suction
- Unscrew the jar's lid and fill the jar 1/3 full or ordinary water (this for an easy cleaning operations and an rapid reaching of the functionally vacuum) then rescrew the lid on the jar correctly.
- Once finished push switch on 0 position and unplug.
- To extract the accessories and start with cleaning.

**WARNING:** Ensure that the IN marker on the filter is on the side facing the collection jar lid and fitted into the "VACUUM".

A wrong connection causes immediate destruction in case of contact with sucked liquids.

### Battery functions

The internal battery provides 45 to 50 minutes of continuous use.

After use, the unit must be recharged for 120 to 150 minutes. A standard cycle would be to use the equipment for 45 minutes and then recharge it for 120 minutes. A complete recharge takes 8 hours.

To recharge, connect the power cable to a 230V outlet and ensure the power switch on the unit is turned on.

Battery function is indicated by the LED lights at the top right of the unit.

The green LED is on if either the battery is fully charged or the mains power has been connected. The yellow LED indicates the unit is charging.

The red LED indicates the battery is fully discharged.

**NEVER USE THE DEVICE WITHOUT JAR AND / OR PROTECTION FILTER**

**MAKE SURE THAT CHILDREN AND/OR MENTALLY ILL PEOPLE DO NOT USE THE DEVICE WITHOUT ADULT SURVEILLANCE**

**ALWAYS PLACE THE DEVICE IN POSITIONS FOR EASY DISCONNECTION**

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**IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC:**

At the end of its working life, the product must not be disposed of as urban waste. It must be taken a special local authority differentiated waste collection or to a dealer providing this service. Disposing of a household appliance separately avoids possible negative consequences for the environment and health deriving from inappropriate disposal and enables the constituent materials to be recovered to obtain significant saving in energy and resources. As a reminder of the need to dispose of separately the product is marked with a crossed-out wheeled dustbin.

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### WARRANTY

CA-MI warrants its products for **24 months** after purchasing date. In front of this warranty, CA-MI will be obliged only to repair or substitute free of charge the products or parts of them that, after verification effected on our factory, or our authorized Service Center, by the Technical Service, results defective. The product must be accompanied by a description of the defect. The warranty, with exclusion of responsibility for direct and indirect damages, is thought limited to the solos defects of material or workmanship and it stops having effect when the device results however gotten off, tampered or sheltered out of the Factory or from the Authorized Service center.

The commodity always travels to risk and danger of the buyer, without any responsibility of CA-MI for damages caused by the transport or dismay from the vector.

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