
GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 - fax 199 400 403
Export: tel. +39 02 953854209/221/225 fax +39 08 95380056
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

SCOLIOMETRO - SCOLIOMETER SCOLIOMÈTRE - DAS SKOLIOMETER ESCOLIÓMETRO - ESCOLIÔMETRO ΣΚΟΛΙΟΜΕΤΡΟ - ال ج ن ف م ق ياس

Manuale utente - User manual - Notice d'utilisation
Betriebsanweisungen - Manual de uso - Manual de uso
Εγχειρίδιο - دليل الإستعمال والرعاية

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SCOLIOSIS is a deviation (morphological or functional) of the spinal column to the side, to the left or to the right. Often, scoliosis is first diagnosed by finding:

- 1) irregular vertical alignment of spinal hypophyses of the vertebrae (misalignment);
- 2) asymmetry of shoulders and hips;
- 3) prominence of one shoulder blade;
- 4) a posterior rib hump in the chest wall.

The first three checks above can be easily done by observing the patient from behind and in erect position. It is more difficult to observe the entity of the posterior rib hump in the chest wall, which is formed as vertebrae rotate around the spinal column. In fact, scoliosis causes the patient to modify his posterior rib position, so that both sides are no longer on a horizontal plane.

The concave side of the scoliosis (**Fig. 1**, left side for observer) is curved inwards, therefore the posterior ribs are lowered, while on the convex side (**Fig. 1**, right side for observer) it is more protruding.

The hump form can be seen more clearly by examining the patient's back view and having him bend his chest forward over his hips.

The non-horizontal alignment is thus visible. On this subject, for an accurate overall evaluation, it is most important that the patient bends forward correctly, feet together, arms hanging beside his legs and in a relaxed position (**Fig. 2**).

Various studies have shown that about 5% of school-age children have sideways bends.

So careful screening by the family doctor, from the age of 6 till puberty, should be one of the routine checkups carried out.

The **Scoliometer**, a useful instrument for doctors in their everyday practice, diagnoses and measures rib hump. Based on the spirit-level principle,


fig. 1

fig. 2

it has been adapted to allow measurement of the different heights of both parts of the thorax, with respect to horizontal, of scoliosis sufferers.

In order to quantify the amount of scoliosis, expressed in degrees, simply place the **Scoliometer**, without pressure, upon the patient's back, with the concave part of the instrument on the spinal apophysis. Then read on the scale the torso's angle of rotation (the angle between horizontal and an axis touching the back of the torso on top of the rib hump).

Fig. 3 shows how the upper torso is measured, with the patient bending over slightly.

Fig. 4 shows the correct position of the thoracic area over the hips for measuring the lumbar part of the spinal column, which requires a full forward bend.

The instrument contains a small sphere that slides within a track, directly indicating on the scale divergence from horizontal in degrees.

This screening makes it possible to identify the defect and quantify it.

Results greater than or equal to 5 degrees require attention.

It is also possible to track results over time,

thus obtaining a precise idea

of the advancement, if any,

of the pathological process,

while also reducing X-ray exposure.


fig. 3

fig. 4

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA. During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included. All components subject to wear are not included in the warranty. The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic elds, radio interferences, etc. The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed. The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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