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

## Warranty

We guarantee this device for five years from date of purchase. We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. In the case of a complaint of a defect in the product, the buyer must prove that any such defect was present at the time the product was delivered. This guarantee excludes any defects which occur because of wear and tear, incorrect use, the use of non-original HEINE parts. Also excluded are repairs or modifications made by persons not authorised by HEINE or cases where the customer does not follow the instructions supplied with the product. Further claims, in particular claims for damage not directly related to the HEINE product, are excluded.

For **repairs**, please contact your supplier, who will return the goods to us or our authorized agent.

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## Instructions for use

### HEINE Disposable Sigmoidoscope/ Proctoscope/Anoscope

Read these instructions carefully before using the HEINE Disposable Sigmoidoscope/Proctoscope/Anoscope attachment and keep them in a safe place for future reference.					

#### Safety information

**Range of applications:** The HEINE disposable sigmoidoscopes, proctoscopes and anoscopes should be used once only. The should be fitted to the instrument head for UniSpec tubes with fibre optics illumination and used only for examination and treatment of the anal passage.

**Safety in use:** Use only light sources for the instrument head which have been approved for use in a medical environment. This tubes should only be inserted when an obturator is fitted and when attached to the instrument head.

#### 1. Preparing the Instrument

This is very simple, if the following order is observed:

**Loosen the sealing screw (2) and open the window (1), which should be removed by extracting it upwards out of the instrument head.**

**Attach the instrument to the handle of the light source.**

**Insert the extension-piece (4) as far as possible into the instrument head and secure this by twisting to the right.**

**If required connect insufflation bulb to the window (5) and attach window to head.**

**Insert the obturator.**

**Apply a suitable lubricant to the tip of the obturator and the tube.**

#### 2. Practical hints for use during the examination

The plane glass window (1) can be opened or closed at any time. The instrument can be made airtight by gently tightening the screw (2). If an insufflation bulb without a valve is used, the insufflated air can be regulated by gently loosening or tightening the screw.

If magnification of the object is required, insert the loupe into one of the two holes provided (3). The loupes are marked P for procto- and anoscope and R for recto- or sigmoidoscope. Both swivel to allow the window to be opened or closed. The telescope is inserted into both of the location holes (3) and clamped by means of a clamping screw. During instrumentation, the telescope can be folded downwards. (The telescope cannot be used with anoscopes).

The vacuum tube accessory can also be located in the holes (3) by means of the attachment. It can then be moved backwards or forwards along the tube and remains steady in any given position. This does not limit the use of a loupe in any way.

#### 3. Care and maintenance

##### 1. Instrument head

###### 1.1 Cleaning

We recommend washing with a soft brush and soapy or plain water. If machine washing or chemical cleaners are used, the manufacturer's instructions regarding duration and concentration must be followed.
**No ultrasonic cleaning.**
**After cleaning, rinse thoroughly with demineralized water and dry at 65 °C max.**

**1.2 Disinfection**
We recommend soaking in disinfectant or thermochemical disinfection at up to 65 °C in a washer sterilizer. Adhere strictly to manufacturer's instructions regarding concentration and duration.
After disinfection, rinse with sterile water and wipe dry with a sterile cloth.

**1.3 Sterilization**
After cleaning as in 1.1 the components can be gas-sterilized with ethylene oxide at up to 65 °C.
Autoclaving can be used occasionally, but repeated autoclaving will reduce light transmission and life expectancy of the optical fibers.
**Flash autoclaving and hot-air sterilization are not recommended.**

##### 2. Disposable tubes

The tubes are delivered clean and sealed in a transparent plastic sleeve. If desired, they can be sterilized with Ethylene Oxide. These tubes are for single use only. An attempt to clean them could result in damage to the tube or obturator and injury to the patient.

##### 3. Accessories

**Sealing window, loupes and vacuum tubes** can be sterilised by any normal method.
**The telescope** may only be gas sterilized.
**Attention:** The polished end-surfaces of the optical fibers in the instrument head must be handled carefully to avoid damage to them which might lead to a reduction in light intensity.

#### 4. Accessories

<b>Biopsy Forceps</b>	
straight jaws, working length 430 mm	<b>E-00.18.921</b>
30° offset jaws, 430 mm	<b>E-00.18.922</b>
straight jaws, working length 350 mm	<b>E-00.18.923</b>
30° offset jaws, 350 mm	<b>E-00.18.924</b>

<b>Prism Telescope</b> with adjustable eyepiece, 3–5, 5 x magnification	<b>E-00.18.908</b>
<b>Suction Tube</b> , for all tubus length	<b>E-00.18.907</b>
<b>Sponge holder</b> , triple metal jaws	<b>E-00.18.906</b>
<b>Rubber dilation tube (double)</b> for connection to the end-window	<b>E-00.18.105</b>

<b>Disposable UniSpec tubes, dia. 20 mm</b> complete with obturator in foil pack	
<b>Sigmoidoscope</b> (250 mm) sales pack with 25 tubes, boxed	<b>E-03.18.825</b>
Case with 4 sales packs = 100 tubes	<b>E-03.18.811</b>
<b>Proctoscope</b> (130 mm) sales pack with 25 tubes, boxed	<b>E-03.19.825</b>
Case with 4 sales packs = 100 tubes	<b>E-03.19.811</b>
<b>Anoscope</b> (85 mm) sales pack with 25 tubes, boxed	<b>E-03.19.025</b>
Case with 4 sales packs = 100 tubes	<b>E-03.19.911</b>

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