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We would like to introduce you to the Rydel-Seiffer Neurological Tuning-Fork with a C 64 damper and a c 128 base.

This tuning-fork is for qualified medical personnel use only.

The AB-125 model can be used to diagnose those nervous affections defined as polyneuropathy and manifesting themselves in the form of reduced sensitivity to vibrations. This symptom is found in the following diseases:

- Diabetes
- Toxic nerve damage
- Alcoholism
- Neuritis of bacterial origin

The Rydel-Seiffer Tuning-Fork produces two different frequencies: note C = 64 Hz with the dampers (as marked on the dampers) and note c = 128 Hz without the dampers (as marked on the tuning-fork itself). The dampers must be fitted with the lower edge flush with the gradation on the prongs. The imprint "C 64" on the dampers, the imprint "c 128" on the front of the tuning-fork and the gradation must all be facing the user when the damper is fitted. The knurled-head screws are tightened to fix the dampers and slackened to remove them. They should be tightened firmly when the dampers are fixed.

Practical to use:

There is a black and a white triangle on the dampers, each with the scale from 2 to 8.

When the tuning-fork is struck against the ball of the thumb, the prongs start to oscillate and two triangles are visible on each weight.

Never strike the tuning-fork against a hard object!

As the oscillation decreases, the two triangles gradually move closer together again and their point of intersection moves slowly upwards. The scale shows the intensity of the vibrations. This can be read off either the white or the black triangle, depending on the amount of light, available or which scale is more easily readable. For the assessment of insensitivity to vibrations, the tuning-fork is struck and its plastic base is placed against the olecranon or the *processus mastoideus*.

The patient must then indicate the point at which he ceases to sense vibration. The intensity of the vibration at the moment is read off the scale. The following values can be regarded as normal:

Patients under 60 years of age: between 6 and 8

Patients over 60 years of age: between 4 and 8

The sensitivity to vibration in the region of the feet can then be checked by placing the tuning-fork against the head of the fibula, the inner malleolus, the head of the os metatarsal and the tip of the big toe. The patient must once again indicate the point at which vibrations cease to be felt. When neuropathy is present, one generally finds a symmetrical decrease in sensitivity to vibrations at the periphery as compared with the values obtained proximally, especially on the forefeet. The extent of the neuropathy can be determined more exactly by placing the tuning-fork at successive points along the edge of the tibia.

Patients with diminished sensitivity caused by sensory neuropathy are especially prone to foot injuries because these are generally not noticed until too late.