Test Rapido Chlamydia
Chlamydia Rapid Test Card
Test rapide Chlamydia
Chlamydia Schnell Testkassette
Prueba Rápida de Chlamydia
Σύστημα Ταχείας Εξέτασης Χλαμυδίων

MANUALE D'USO
USER MANUAL
MANUEL D'UTILIZATION
BEDIENUNGSANLEITUNG
MANUAL DE USO
Εγχειρίδιο χρήσης

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BIBLIOGRAFIA

1. Sanders J.W. et al *Evaluation of an Enzyme Immunoassay for Detection of Chlamydia trachomatis in Urine of Asymptomatic Men*. J.Clinical Microbiology, 32,24-27, (1994).
2. Jaschek, G. et al *Direct Detection of Chlamydia trachomatis in Urine Specimens from*

Symptomatic and Asymptomatic Men by Using a Rapid Polymerase Chain Reaction Assay. J. Clinical Microbiology, 31, 1209-1212, (1993).

3. Schachter, J *Sexually transmitted Chlamydia trachomatis infection.* Postgraduate Medicine, 72, 60-69, (1982).

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0197

Tamponi:

 Puritan Medical Products
Company LLC
Fabbricante Guildford, ME 04443-0149, USA

0086

Methodology

Rapid test for the qualitative determination of the Chlamydia antigen on the cervical swab in women and urethral swab or urine sample in men.

Only for professional diagnostic in vitro use.

CORRECT USE

The Chlamydia Rapid Test Card (Swab/Urine) is a rapid immunochromatographic test for the qualitative determination of Chlamydia trachomatis on the cervical swab in women and the urethral swab or urine sample in men, to be used during the diagnosis of Chlamydia infection.

SUMMARY

Chlamydia trachomatis is the main cause of STD in the world. It is formed by elementary bodies (infective form) and reticulate bodies (non-infective form).

Chlamydia trachomatis has a very high transmission level and often shows no symptoms, with frequent complications both in women and newborns. The complications due to Chlamydia in women include: cervicitis, urethritis, endometritis, pelvic inflammation (PID) and an increase in the number of extra-uterine pregnancies and infertility.¹ Vertical transmission during birth, between the mother and newborn, can cause inclusion conjunctivitis or pneumoniae. In men, the complications due to Chlamydia infection include: urethritis and epididymitis. At least 40% of urethritis nongonococcal are associated to Chlamydia infection. About 70% of women with endocervical infections and over 50% of men with urethral infections are asymptomatic. Generally, Chlamydia infection was diagnosed through the determination of the Chlamydia inclusion in tissue cultures. The cultural method is the most sensitive and specific of all laboratory methods, but it is demanding for the laboratory, expensive, long (48-72 hours) and it is not a routine method in most laboratories. The Chlamydia Rapid Test Card (Swab/Urine) is a rapid test for the determination of the qualitative determination of the Chlamydia antigen on the cervical swab in women and urethral swab or urine sample in men, able to produce results in 10 minutes. The test uses a specific antibody for Chlamydia for the selective determination of the Chlamydia antigen on the cervical swab in women and urethral swab or urine sample in men.

PRINCIPLE

The Chlamydia Rapid Test Card (Swab/Urine) is a qualitative test with lateral flow for the determination of the Chlamydia antigen on the cervical swab in women and urethral swab or urine sample in men, to be used as aide in the diagnosis of Chlamydia infection.

In this test, a specific antibody directed towards the Chlamydia antigen is placed on the test area of the card. During the test, the antigen solution reacts with a Chlamydia antibody stuck to the particles. The mixture migrates until it reacts with the anti-Chlamydia antibody on the membrane and produces a colored strip in the test area. The presence of this colored strip in the test area indicates a positive result, whereas a lack of this colored strip in the test area indicates a negative result. As an internal control of the procedure, there will always be a colored strip in the control area (C) of the test which shows that a correct quantity of sample was used and that the migration towards the membrane took place correctly.

REACTANTS

The membrane of the card contains antibodies linked to the Chlamydia particles and antigen adhered to the membrane.

CAUTION

- Only for professional in vitro diagnostic use. Do not use after expiration date..
- Keep card in closed package until use.
- Do not eat, drink or smoke in the area where the kits and samples are handled.
- Do not use card if package is damaged.
- Handle all the samples as potentially infective. Take all necessary measures against microbiological risks, during all analysis phases, and follow the standard procedures for the correct disposal of the samples.
- During the analysis of the samples wear personal protection equipment: laboratory coats, disposable gloves and eye protection.
- After use, the test must be eliminated according to the local regulations in force.
- Humidity and temperature can negatively affect the results.
- For endocervical samples, use **sterile** swabs only.

PRESERVATION AND STABILITY

Keep in closed package and room temperature or refrigerated (2-30°C). The card is stable up to the date of expiration indicated on the label of the package. The card must be kept in the closed package up to the moment of use. **DO NOT FREEZE**. Do not use after the expiration date.

TAKING AND PREPARING THE SAMPLES

- The Chlamydia Rapid Test Card (Swab/Urine) can be performed on a cervical swab in women and urethral swab or urine sample in men.
- The quality of the sample taken is extremely important. The determination of Chlamydia requires a complete sample collection technique with firm movements, able to take cellular material and not only bodily fluids.
- How to take the **Sample from Female Cervical Swab**:
 - Only use the swab supplied in the kit. Alternatively it is possible to use a plastic swab.
 - Before collecting the sample remove the excess mucus in the endocervical area with a cotton swab and eliminate. The swab must be inserted in the endocervical channel, beyond the squamous-columnar junction until most of the swab is no longer visible. It will be possible to take columnar cells or cubic epithelium, which are the main reserve of Chlamydia of the body. Firmly rotate the swab by 360° in one direction (clockwise or counter-clockwise), wait 15 seconds and remove swab. Avoid contamination with extra-cervical or vaginal cells. Do not use sodium hydrochloric 0.9% to treat the swab, before taking the sample.
 - If the test must be carried out immediately insert the swab in the extraction test tube.
- How to take the **Sample from the Male Urethral Swab**:
 - Standard plastic swabs or sterile metal swabs can be used to collect the urethral sample. Tell patient not to urinate for at least one hour before the collection of the sample.
 - Insert the swab in the urethra for 2-4 cm, rotate the swab by 360° in one direction (clockwise or counter-clockwise), wait 10 seconds and remove swab. Do not use sodium hydrochloric 0.9% to treat the swab, before taking the sample.
 - If the test must be carried out immediately insert the swab in the extraction test tube.
- How to take the **Male Urine Sample**:
 - Collect 15-30 ml of the first morning urine in a sterile urine container. The first morning urine is preferable because it contains a high concentration of the Chlamydia antigen.
 - Mix the urine sample by turning the container upside-down. Transfer 10 ml of the urine sample in the test tube for centrifugation, add 10 ml of distilled water and spin up to 3,000 rpm for 15 minutes.
 - Carefully eliminate the supernatant, keep the test tube upside-down and remove all the

supernatant from the edges of the test tube, tapping on blotting paper.

- If the test must be performed immediately treat the spun urine as indicated in the instructions for use.
- We recommend performing the test as soon as possible after collecting the sample. If it is not possible to perform the test immediately, the samples must be kept in a dry test tube for preservation or for transportation. The swabs can be kept for 4-6 hours at room temperature (15-30°C) or 24-72 hours if refrigerated (2-8°C). The urine samples can be refrigerated (2-8°C) for 24 hours. Do not freeze. All the samples must be brought to room temperature (15-30°C) before carrying out the test.

COMPLETE KIT

Supplied Material

- Card
- Test Tubes
- Dropping cap
- Test Tube holder/reactant
- Sterile swabs for female cervical sample
- Methodology
- Measuring pipette
- Reactant A (0,2M NaOH)

Xi

R36/38 Irritating for eyes and skin.

(S2 Keep out of reach of children.)

(S46 If swallowed seek immediate medical attention and show container or label).

S60 This material and its container must be disposed of as dangerous waste.

- Reactant B (0,2N HCl)

C

R35 Causes serious burns.

S1/2 Keep locked and out of reach of children.

S45 In case of accident or sickness seek immediate medical attention (if possible, show label).

S60 This material and its container must be disposed of as dangerous waste.

Necessary Material but Not Supplied

- Timer
- Test Tube for spinning (only for male urine samples)
- Positive control
- Container for urine (only for male urine samples)
- Sterile swab for male urethral samples
- Negative control

PROCEDURE

Bring card, the samples, the reactants and/or the controls to room temperature (15-30°C) before the analysis.


1. Remove the card from the sealed envelope and use it as soon as possible. The best results are obtained if the analysis is performed immediately after opening the package.
2. Extract the Chlamydia antigen according to the type of sample taken.

For swabs with Female Cervical Samples or Male Urethral Samples:

- Keeping the bottle of Reactant A vertically, add **5 drops of reactant A** (about 300 µl) in a test tube for the extraction of the sample. Reactant A has no color. Immediately add the swab, press the bottom of the test tube and rotate the swab 15 times. Wait 2 minutes.
- Fill the measuring pipette for Reactant B up to indicated level (about 220 µl), then **add Reactant B** in the sample extraction test tube. The solution will become cloudy. Press the bottom of the test tube and **rotate the swab 15 times**, until the solution is clear again, with a green or blue color. If the swab shows traces of blood, the color of the solution will become yellow or brown. Wait 1 minute.
- Press the swab against a side of the test tube and extract it by pressing the test tube. Put the dropping cap onto the test tube to extract sample.

For Male Urine Samples:

- Fill the measuring pipette for **Reactant B** up to indicated level (about 220 μl), then **add Reactant B** in the spun urine in the test tube, then shake fluid well, taking and releasing fluid with a pipette until the suspension is homogeneous.
 - Transfer the content of the test tube used for centrifugation to a test tube for the extraction of the sample. **Wait 1 minute**, keep the bottle of **Reactant A** in vertical position, **add 5 drops of Reactant A** (about 300 μl) in the test tube for the extraction of the sample. Rotate or tap the bottom of the test tube to mix solution. **Wait 2 minutes**.
 - Apply the dropping cap to the test tube for the extraction of the sample.
3. Position the card on a flat and clean surface. **Transfer 3 drops of the extracted solution** (about 100 μl) in the pit of the sample (S) of the card and start timer. Avoid the formation of air bubbles in the sample pit (S).
4. Wait for the red strip/s to appear. **Read the results after 10 minutes**. Do not read results after 20 minutes.

INTERPRETATION OF RESULTS

(See previous illustration)

POSITIVE:* **Two separate colored strips appear.** One in the control area (C), the other in the reactive area (T).

***NOTE:** The color intensity of the strip in the reactive area (T) could change, but the test must be considered negative if just a colored shadow appears.

NEGATIVE: **A colored strip appears in the control area (C).** No colored strip appears in the reactive area (T).

NOT VALID: **No colored strip appears in the control area.** The most likely reasons for this are an insufficient volume of sample or an incorrect procedure. Review procedure and repeat the test with a new card. If the problem continues, immediately stop using the kit and contact local distributor.

QUALITY CONTROL

Every test includes an internal procedure control. The colored strip that appears in the control area (C) is a positive internal control of the procedure and it means a sufficient volume of sample was taken, the migration along the membrane is correct and the procedure is correct.

No control standards are supplied with the kit; we recommend analyzing a positive and a negative control as proper laboratory practice to confirm the analytical procedure and verify the quality of the test.

LIMITATIONS

1. The Chlamydia Rapid Test Card (Swab/Urine) is exclusively for in vitro diagnostic use. The test must be used for the determination of the Chlamydia antigen on the cervical swab in women and urethral swab or urine sample in men. This qualitative test is not able to determine the quantitative value nor the percentage of increase of the Chlamydia antigen.
2. This test solely indicates the presence of the Chlamydia antigen in the vital and non-vital Chlamydia samples. Do not perform test with other types of samples.
3. The determination of Chlamydia depends on the number of organisms in the sample. This can be influenced by the collection method of the sample and by other patient factors such as: age, past history of STDs, presence of symptoms, etc. The minimum level that can be determined with this test varies according to the serotype. Consequently, the results must be interpreted along with other clinical and laboratory data available to the doctor.
4. It is not possible to determine the success or failure of a therapy, because the antigen could persist after an adequate antimicrobial therapy.
5. An excessive quantity of blood on the swab could produce falsely positive results

EXPECTED VALUES

For women, patients of STD clinics, and for other high risk populations, the prevalence of Chlamydia infection could be between 20% and 30%. For low risk population, such as patients of ob/gyn clinics, the prevalence is about 5% or less.

The data shows that in men, patients of STD clinics, the prevalence of Chlamydia infection is approx. 8% in asymptomatic men and 11% in symptomatic men.^{1,2} The normal rate of Chlamydia in asymptomatic men is less than 5%.³

PERFORMANCES

Sensitivity

The Chlamydia Rapid Test Card (Swab/Urine) was evaluated with samples taken from patients of STD clinics. PCR was used as reference method for the Chlamydia Rapid Test Card (Swab/Urine). The samples were considered positive if the PCR analysis was positive. The samples were considered negative if the PCR analysis was negative. The results show that the Chlamydia Rapid Test Card (Swab/Urine) has a high level of sensitivity compared to PCR.

Specificity

The Chlamydia Rapid Test Card (Swab/Urine) uses a highly specific antibody for the Chlamydia antigen on cervical swabs in women and urethral swabs in men or male urine samples. The results show that the Chlamydia Rapid Test Card (Swab/Urine) has a high level of specificity compared to PCR.

For Samples from Female Cervical Swabs:

Method		PCR		Total Results
Chlamydia Rapid Test Card	Results	Positive	Negative	
	Positive	46	3	49
	Negative	6	87	93
Total Results		52	90	142

relative Sensitivity: 88.5%
relative Specificity: 96.7%
relative Accuracy: 93.7%

For Samples from Male Urethral Swabs:

Method		PCR		Total Results
Chlamydia Rapid Test Card	Results	Positive	Negative	
	Positive	40	8	48
	Negative	11	104	115
Total Results		51	112	142

relative Sensitivity: 78.4%
 relative Specificity: 92.9%
 relative Accuracy: 88.3%

For Samples of Male Urine:

Method		PCR		Total Results
Chlamydia Rapid Test Card	Results	Positive	Negative	
	Positive	20	0	20
	Negative	2	35	37
Total Results		22	35	57

relative Sensitivity: 90.9%
 relative Specificity: >99.0%
 relative Accuracy: 96.5%

Cross-Reactivity

The antibody used for the Chlamydia Rapid Test Card (Swab/Urine) proved to be able to determine all the serotypes of Chlamydia known. The group Chlamydia psittaci and Chlamydia pneumoniae were tested with the Card for the Chlamydia Rapid Test Card (Swab/Urine). Their cross-reactivity with other organisms was evaluated using suspensions of 109 CFU/ml.


The following organisms resulted negative when tested with the Chlamydia Rapid Test Card (Swab / Urine):

Acinetobacter calcoaceticus	Pseudomona aeruginosa	Proteus mirabilis
Acinetobacter spp	Neisseria meningitides	Neisseria gonorhea
Enterococcus faecalis	Salmonella choleraesius Group B/C	Streptococcus
Enterococcus faecium	Candida albicans	Hemophilus influenzae
Staphylococcus aureus	Proteus vulgaris	Branhamella catarrhalis
Klebsiella pneumoniae	Gardnerella vaginalis	

BIBLIOGRAFY

1. Sanders J.W. et al *Evaluation of an Enzyme Immunoassay for Detection of Chlamydia trachomatis in Urine of Asymptomatic Men*. J.Clinical Microbiology, 32,24-27, (1994).
2. Jaschek, G. et al *Direct Detection of Chlamydia trachomatis in Urine Specimens from Symptomatic and Asymptomatic Men by Using a Rapid Polymerase Chain Reaction Assay*. J. Clinical Microbiology, 31,1209-1212, (1993).
3. Schachter, J *Sexually transmitted Chlamydia trachomatis infection*. Postgraduate Medicine, 72, 60-69, (1982).

Index of Symbols

	Attention, ^{see} instructions for use		Tests per kit		Manufacturer
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2-30°C		Lot Number		Catalog 29117
	Irritating		Corrosive		

 **Manufacturer** Innovacon, Inc.
4106 Sorrento Valley Boulevard
San Diego, CA 92121, USA

0197

Swabs:

 **Manufacturer** Puritan Medical Products
Company LLC
Guildford, ME 04443-0149, USA

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BIBLIOGRAPHIE

1. Sanders J.W. et al Evaluation of an Enzyme Immunoassay for Detection of Chlamydia trachomatis in Urine of Asymptomatic Men. J.Clinical Microbiology, 32,24-27, (1994).
2. Jaschek, G. et al Direct Detection of Chlamydia trachomatis in Urine Specimens from Symptomatic and Asymptomatic Men by Using a Rapid Polymerase Chain Reaction Assay. J. Clinical Microbiology, 31,1209-1212, (1993).
3. Schachter, J Sexually transmitted Chlamydia trachomatis infection. Postgraduate Medicine, 72, 60-69, (1982).

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Écouvillons:

Fabricant

Puritan Medical Products
Company LLC
Guildford, ME 04443-0149, USA

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Puritan Medical Products Company LLC