SterilAir Plus - IgenHand

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Why sterilization devices should be used

Aerogenic contamination has always been a serious problem for our society, in fact, many diseases are transmitted through the air.

Meningitis, diphtheria, tuberculosis, measles, German measles, mumps, right up to simple influenza, the common cold and the new varieties of influenza which have appeared over the last few years, such as avian and swine 'flu.

In some environments the problem of aerogenic spread is particularly critical; as a rule, we can take into consideration all those environments where there is **an intensive gathering or flow of people** (public places, meeting points, schools, nurseries etc.), **or medical or surgical environments**. Some medical environments are especially critical, for example dentists' surgeries, where the **sprays and ultrasound instruments used create a greater diffusion** of micro-particles which are then carried into the surroundings.

All of these infections are more and more diffuse, so much so that in some cases they are even defined as pandemic, and this generates heavy costs and enormous social problems. **The use of environmental sterilizers finds immediate application in medical environments to reduce the risk of transmitting air-borne viruses**, on the basis of the recommendations set down by Disease Control and Prevention Centres, however **these Devices can be used in all environments**.

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Air sterilization: SterilAir Plus

Tecno-Gaz S.p.A. has available a system called **SterilAir Plus, certified for the sterilization of the air**, via irradiation and **elimination of all the microbial agents including pathogenic ones**.

Using a disinfection system means:

- **Reducing the risk of contaminating workers**, as also requested by the consolidated act on safety.
- **Reducing or eliminating costs and lost profits**, due to long absences caused by sickness.
- **Reducing or eliminating risks** for patients or the population.
- **Having an environment which is** bacteriologically pure.

SterilAir Plus is an indispensable product, which is low cost and generates a service of enormous social usefulness.

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
The controlled ultraviolet radiations

There are a lot of methods for decontaminating the air, for example by means of free ultraviolet radiation, by means of chemical methods which use iodoforms, chlorhexidine, or quaternary ammonium derivatives sprayed with special appliances. These techniques although considered efficient, have only a temporary effect; they are used just when operators are not present and there is no activity at all. **The most suitable and efficient method is the one that uses ultraviolet radiation having controlled wavelength.** This radiation has an intensive microbicide activity against numerous microorganism. This radiation blocks the microorganism's reproductive capacity by altering their chromosomes. There are no immune microorganisms to germicide radiation: **spores, bacteria, mould, yeast and even virus DNA are destroyed by this radiation.**

The great advantage of our system consists in **the absolute lack of danger (certified) for man** since there is not even the slightest leakage of UV-C radiation from the module.

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Microbiological Analysis

Figure 1 shows the comparison between the data concerning the air outside and the air inside a refuse dump chamber (the latter considered to be a model of acute contamination). The external values were measured to consider their influence on the internal ones. The initial contamination values are those relating to 9.30 a.m.

The first sample at 11.45 a.m. after the installation of the lamp already shows a notable reduction in the load of micro-organisms in the air with respect to the initial load. This reduction, calculated net of the variation in the external microbial load, represents 50%. The final results underline an environmental microbiological situation different from the initial one. If initially the levels of internal microbial load are greater than the external ones, after the treatment with the Spectral UV Lamp the situation has been completely modified, obtaining negligible internal pollution values.

Microbial Reduction

Another significant trial was carried out, taking air corresponding to the input aperture of the SPECTRAL UV lamp and the exit aperture of the same lamp (fig. no.2). From the values shown it can be noticed that the levels of microbiological pollution in the air before undergoing the photolytic treatment are very when compared with those of the air after treatment. In order to assess the sanitizing capacity of the lamp as regards the air, a room of around 70m³ was contaminated in advance with a suspension of Bacillus Subtilis spores, by nebulization. Later on, checks were made on the degree of air pollution at different intervals of time, via the S.A.S. system (air sampler), before and after turning on the Spectral UV Germicide Lamp. The results obtained, (fig. 3) show a drastic reduction in the microbial load in the air over time, as testimony of the effectiveness of the germicidal action of the device.

Fig. 1

Abbattimento microbico UV lampada germicida Spectral UV
Microbial depression UV germicidal lamp Spectral

■ Riferimento esterno UFC/m³
External reference UFC/m³
■ Area interna trattata con UV UFC/m³
Internal treated area with UV UFC/m³


Fig. 2

Abbattimento microbico con trattamento lampada germicida Spectral UV | Microbial depression with germicidal lamp UV Spectral treatment

■ Aria ambiente prima del trattamento UFC/m³
room air before the treatment UFC/m³
■ Aria trattata UFC/m³, abbattimento microbico:72,9%
treated air UFC/m³, microbial depression:72,9%



Fig. 3

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UV lamps against aerogen infections

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Produced with an innovative design and studied to avoid problems of direct and indirect exposure to the shortwave ultra-violet rays (UVC 254 nm) offering the possibility to use it continuously also in the presence of people without any risk for the latter. **The most effective device for the destruction of micro-organisms present in the air.**

They work via a system of closed-cycle forced ventilation. The air sucked into the module first passes through a mechanical filter located in the entrance aperture. Here the larger pollutants are blocked to avoid staining the germicidal lamps. **After this, the air is forced to pass into direct contact with tubes of mercury vapour which, thanks to the emission of UV-C radiation** accomplish the maximum germicidal action. A high-powered reflecting specular screen concentrates the reflections of the UV-C radiation. The air is expelled from the exit aperture to determine the microbiological reduction.

The great advantage of this system consists in the absolute lack of danger for man since there is not even the slightest leakage of UV-C radiation from the module.

It is therefore possible to carry out continuous and constant disinfection of the air in any environment, during every working phase.

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Keep safe your surgery and your business

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Safeguard your professional image and your prestige, while looking after your health, and that of your associates and patients. **STERILAIR PLUS** is the first system to sterilize air.

STERILAIR PLUS is a device to sterilize air which **ensures the elimination of all micro-organisms including spores**. The air is transported into the sterilization chamber inside the device, which, by irradiation eliminates the entire microbial load present in the air. This device can run continuously also in the presence of people, ensuring maximum safeguarding of the operators. **STERILAIR PLUS is simple, silent and ergonomic.** Tecno-Gaz is the only company able to offer a complete sterilization programme which is rational, ergonomic, and safe.

Some of the possible applications of SterilAir Plus systems

Germicidal lamps to disinfect air can be used:

- in hospitals
- in operating theatres
- In Accident & Emergency
- in dentists' surgeries
- in outpatient clinics
- in veterinary surgeries
- in analytical laboratories
- in pharmaceutical industries
- in animal farms
- in premises for the production of food and beverages
- in cold stores
- in dwellings
- in offices
- in public places
- in premises with air conditioning
- etc.

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
The hands are the primary means of spreading cross infection, thus the washing and disinfection of the hands is normally the first and most important procedure to be made official in all preventative medical and social protocols.

“According to a multi-centre study published in the British Medical Journal, **to keep the spread of influenza viruses under control it would be sufficient to wash the hands.** By inspecting the results of 51 previous investigations which analyzed the efficacy of a series of measures to limit the transmission of respiratory viruses both from animals to man and from man to man, the researchers underlined how, within specific national prevention plans, **priority should be given to the realization of simple and inexpensive measures of physical protection** (such as frequent hand-washing).

In fact, from the remarks of the experts it emerges that the use of vaccines and antiviral medicines could be insufficient to interrupt the spread of influenza and that hand-washing, instead, generates an efficient “barrier” effect.

Also the Ministry of Health, as well as advising the use of paper handkerchiefs, to be thrown away after use, and good respiratory hygiene (covering the mouth and nose when sneezing or coughing), has recently reaffirmed the importance of frequent hand-washing as a preventative measure against influenza contagion. Some studies had already demonstrated that the “promotion” of hand-washing amongst children at nursery school was able, on its own, to reduce the incidence of respiratory infections by 14 per cent in Canada, 12 per cent in Australia and a good 32 per cent in the United States. And what is more: **it has been demonstrated that washing the hands at least 5 times a day reduces paediatric visits for respiratory infections by 45%.** Moreover, in developing countries, where the highly contagious infections are a cause of a high mortality rate, programmes to “encourage” the use of these hygienic regulations have been launched. Good. In one of these countries, Pakistan, it has been demonstrated that **a project to encourage hand-washing with soap has reduced by 50 per cent incidences of pulmonitis, gastroenteritis and also pyodermitis in children under 5 years of age”**

Taken from: http://it.health.yahoo.net/c_special.asp?id=23914&c=13&s=1

Obviously the washing and disinfection of the hands play an even more important role for health workers, dentists and all those who practice therapy and other operations on people.

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Tecno-Gaz S.p.A. offers IgenHand, a tank loaded with liquid disinfectant a device of the utmost usefulness, which helps to rationalize the disinfection of the hands.

This device permits the **best possible disinfection, without wasting time, without activating levers, buttons or using gel or liquids to be spread. IgenHand is a device which must be pre-loaded with a suitable liquid disinfectant**, which is then finely vaporized and scattered over the hands in a far-reaching way, when they are inserted under the control photocell. **IgenHand is a comprehensive device**, which can count on **solutions designed to be used in all environments**. The vaporization of the liquid is micro-metric and this allows **far-reaching disinfection and rapid drying**.

1 Serbatoio di carico liquido disinfettante
Tank loaded with liquid disinfectant

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All rights reserved. Variations can be done without notice. Tecno-Gaz S.p.A. is to be considered not responsible for damages caused by the lack or the wrongness of the information here mentioned.

www.tecnogaz.com