Istruzioni per l'uso
Instructions for use
Mode d'emploi
Bedienungsanleitung
Instrucciones para el uso

TECNO40

TECNO70

TECNO100

VASO/POT 2L

VASO/POT 4L

I

GB

F

D

E

	PAGINA - PAGE PAGE - SEITE - PÁGINA
I GB F D E	1 - 7
I	8 - 16
GB	17 - 25
F	26 - 34
D	35 - 43
E	44 - 52

**AVVERTENZE GENERALI E INFORMAZIONI AL DESTINATARIO
CONTENT ORGANIZATION PRESENTATION AND CONSULTATION MODES
GÉNÉRALITÉS ET INFORMATIONS DESTINÉES AU CLIENT
ALLGEMEINE HINWEISE UND INFORMATIONEN FÜR DEN ANWENDER
ADVERTENCIAS GENERALES E INFORMACIONES AL DESTINATARIO**

1.1

**DICHIARAZIONE DI CONFORMITÀ / CONFORMITY STATEMENT
DÉCLARATION DE CONFORMITÉ / KONFORMITÄTSERKLÄRUNG DECLARACIÓN DE
CONFORMIDAD**

*La sottoscritta - The company- La soussignée - Die unterzeichnende - La que suscribey:
TECNO-GAZ S.p.A. - Strada Cavalli n.4 - 43038 Sala Baganza - Parma - Italia.*

*Dichiara sotto la propria responsabilità, che il Dispositivo Medico :
Declares under its own responsibility that the Medical Device :
Déclare sous sa propre responsabilité que l'appareil à usage médical :
Erklärt auf eigene Verantwortung, dass das Medizinprodukt :
Declara bajo la propia responsabilidad, que el Dispositivo Médico:*

*ASPIRATORE CHIRURGICO CARRELLATO,costruito nella versione:
SURGICAL SUCTION UNIT WITH TROLLEY ,available in the following version :
ASPIRATEUR CHIRURGICAL SUR ROUES,dans la version :
CHIRURGISCHER ASPIRATOR AUF ROLLEN, in der Ausführung :
ASPIRADOR QUIRÚRGICO SOBRE RUEDAS,fabricado en la versión:*

INCOLLARE ETICHETTA LOTTO
LOT PLATE
COLLER L'ÉTIQUETTE DU LOT
CHARGENETIKETT EINKLEBEN
APLICAR LA TARJETA LOTE

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- Manufactured by TECNO-GAZ S.p.A. complies with Medical Equipment Directive 93/42/EEC of 14 June 1993, which corresponds to the Italian Decreto Legislativo N.46 of 24 February 1997.

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WARNING

- BEFORE USING THE DEVICE, READ THE INSTRUCTIONS
- BEFORE USING THE DEVICE, CHECK THE CORRECT WORKING OF ITS FUNCTIONS
- In case of possible contamination or discolouration, change the filter immediately. Don't use the suction unit without the protection filter fitted. If the suction unit is used in an emergency or in a patient where the risk of contamination is not know the filter must be changed after each use.
- DURING ITS WORKING, CHECK THE ASPIRATION CIRCUIT
- THE SUCTION UNIT MUST BE USED BY SKILLED STAFF
- DO NOT OPEN THE DEVICE, ADDRESS TO A CUSTOMER SERVICE FOR ASSISTANCE.
- DO NOT USE IT IF THERE ARE ELECTROMAGNETIC INTERFERENCES WHICH COULD INFLUENCE THE CORRECT WORKING OF THE ASPIRATOR
- The machine contains a sealed leaden – calcium battery, it doesn't need maintenance. For the replacement address to the customer services.
- Check the integrity of the device and of its parts, possible failures, cracks or breaking.
- Make sure that the feeling cables (12V, 230V) are whole.
- Don't use the plugs when your hands are wet.
- Don't plunge the device in water.
- Keep away from heat sources and don't expose to atmospheric agents.
- Use the device for the aims expected by the manufacturer.
- DON'T USE THE DEVICE IN POTENTIALLY EXPLOSIVE ATMOSPHERE.

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WARRANTY

THE SUCTION UNIT WITH TROLLEY is under warranty for twelve months except for the parts subject to wear.

This warranty starts from the delivery to the customer and is proved by the return of the coupon properly filled in, stamped and signed by the seller.

If there are any objections, the date on the sale bill is considered as valid.

The reparation or the replacement of parts under warranty is made after unquestionable permit of the manufacturer and doesn't include the staff travelling expenses, the packing expenses and the delivery expenses.

The warranty doesn't include the failures or damages caused by a wrong maintenance, lack of carefulness, negligence or other reasons which cannot be attributed to the manufacturer.

Any kind of compensation for people or things is therefore excluded.

The whole suction unit cannot be replaced.

No compensation for the not working of the device is possible.

The warranty automatically stops if the customer or not authorized people repair or modify the device.

The customer must address only to the dealer or the customer services indicated by the manufacturer.

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INDEX

PAR	DESCRIPTION	PAGE
01	INFORMATION FOR USERS	1-7
02	INTRODUCTION OF THE PRODUCT	18
03	PRESENTATION OF THE CONTENT AND WAY OF CONSULTATION	19
04	DATA AND TECHNICAL FEATURES	20
05	INSTALLATION	21
06	WORKING AND USE	22
07	INSTRUCTIONS FOR THE USER	23
08	ORDINARY AND PLANNED MAINTENANCE	24
09	ANALYSIS AND RESEARCH OF FAILURES	25
10	COMMERCIAL DATA, SPARE PARTS AND DOCUMENTATION	25

2.1 COMPOSITION

The SUCTION UNIT WITH TROLLEY has a structure made of thermoplastic material (PTS) equipped with control board and place for 2 pots .

The surface in thermoplastic material (PTS) is equipped with 3 pivoting wheels made of grey-nylon (anti-marks), the two rear wheels are equipped with brake.

The pump doesn't unload anything in the environment.

The control board is easy to use.

The switches have a IP65 protection and the device is IPX1 (protected against dripping).

The device is equipped with rubber feet ensuring stability.

The pot is made of transparent polycarbonate(in polysulfone upon request) and its capacity is 2L or 4L on request, marked with their own value, expressed in millimetres.

Thanks to a gradation the user has an indication of the volume, with intervals of 100 ml.

DIFFERENT VERSIONS:

SURGICAL SUCTION UNIT WITH TROLLEY

ITEM	REFERENCE	DESCRIPTION
TECNO40	121-A	SUCTION UNIT WITH TROLLEY 40 L/MIN WITH 2 L POT
TECNO40-4L	144-A	SUCTION UNIT WITH TROLLEY 40 L/MIN WITH 4 L POT
TECNO70	122-A	SUCTION UNIT WITH TROLLEY 70 L/MIN WITH 2 L POT
TECNO70-4L	145-A	SUCTION UNIT WITH TROLLEY 70 L/MIN WITH 4 L POT
TECNO100-4L	146-A	SUCTION UNIT WITH TROLLEY 100 L/MIN WITH 4 L POT

PORTABLE SURGICAL SUCTION UNIT

ITEM	REFERENCE	DESCRIPTION
TECNO16-B	111-A	PORTABLE SUCTION UNIT WITH BATTERY 16L/MIN WITH 1 L POT
TECNO16-B-2L	114-A	PORTABLE SUCTION UNIT WITH BATTERY 16L/MIN WITH 2 L POT
TECNO15	112-A	PORTABLE SUCTION UNIT 15 L/MIN WITH 1 L POT
TECNO15-2L	115-A	PORTABLE SUCTION UNIT 15 L/MIN WITH 2 L POT
TECNO25	113-A	PORTABLE SUCTION UNIT 25 L/MIN WITH 1 L POT
TECNO25-2L	116-A	PORTABLE SUCTION UNIT 25 L/MIN WITH 2 L POT

* Some versions have been omitted, for further information please contact TECNO-GAZ S.p.A

**3.1
GLOSSARY**

VACUUM	Pressure lower than the atmospheric pressure
FREE AIR FLUX	Air flux without limits passing through an appointed entry
LOW VACUUM	No more than 20kPa under the atmospheric pressure
MEDIUM VACUUM	Between 20kPa and 60kPa under the atmospheric pressure
HIGH VACUUM	More than 60kPa under the atmospheric pressure
LOW FLUX	Less than 20 l/min
HIGH FLUX	More than 20 l/min
VACUUM PUMP	Device for generating vacuum
VACUUM INDICATOR	Device for checking the vacuum level
VACUUM REGULATOR	Device for checking the highest vacuum applied to the patient
FILTER	Device for holding particles
ENDING PART	Part of the aspirator applied to the patient from the ending part to the first removable joint
INTERMEDIATE TUBES	Tube between the pot and the vacuum source
ASPIRATION TUBE	Tube for the liquids and/or solid particles transport from the ending part to the pot
ASPIRATION	Application of vacuum for the aspiration of fluids and/or solid particles
POT	Container for liquids and solid particles

**3.2
REMARKS ON THE GRAPHIC CAPTION**

MEASUREMENT OF THE VACUUM LEVEL (VACUUM-METER)

REGULATOR OF THE VACUUM LEVEL

VERSION	TECNO40	TECNO70	TECNO100-4L
EN ISO 10079 - 1	HIGH VACUUM - HIGH FLOW		
Surgical Suction Unit	HIGH VACUUM - HIGH FLOW		
EN 60601 - 1			
Protection against electric dangers	Class 1		
Protection against direct/indirect contacts	TYPE B		
Level of protection against the penetration of substances	IPX1		
Level of protection in zone AD-FT (CEI64-4)	NOT SUITABLE		
Level of safety with anesthetic and /or explosive mixtures	NOT SUITABLE		
Use in field structures	NOT SUITABLE		
Use on means of transport	NOT SUITABLE		
Use for thoracic drainage	NOT SUITABLE		
TECHNICAL FEATURES			
Power supply	230V ~ 50 Hz		
Absorbed electric power	240 W	420 W	405 W
Intensità di corrente	1.02 A	1.06 A	1.09 A
Intensità di corrente di spunto	3.09 A	5.07 A	5.08 A
Electric protection	2 fuses Ø 5x20 F 3.15 AL 250V	2 fuses Ø 5x20 F 6 AL 250V	
ASPIRATION			
Pump	WITH PISTON-NOT OILED		
Pump nominal flow	50 Lt/min	90 Lt/min	120 Lt/min
Flux at use- max. established rate	40 Lt/min	70 Lt/min	100 Lt/min
Max. depression	-0.9 bar		
POTS			
N° pots	2		
Capacity	Version 2 Lt and version 4Lt		
Material	POLYCARBONATE OR POLYSULFONE UPON REQUEST		
Single use collection bag	UPON REQUEST		
Non-return valve system	VALVE WITH FLOAT		
HYDROFOBIC FILTER			
Filtering Capacity	0.3 µm (99.999%)		
Material	MICRO-GLASS FIBREPAPER		
TUBE			
Material,diameter	SILICON FOR MEDICAL USE ØINTERNAL= 8 ØEXTERNAL= 14		
VACUUM REGULATION			
Valve for vacuum regulation	PRESENT WITH SCALE		
Vacuum-meter of graduated control	Ø 63 0 / -1 bar		
DIMENSIONS			
Length (mm)	540 mm		
Depth (mm)	480 mm		
Height with 1 l pot	850 mm		
Weight	16 kg	19 kg	20 kg
ACOUSTIC LEVEL	/		
LIMITS OF USE	SEE INSTRUCTION MANUAL		
CONFORMITY	MDD 93/42 EEC - EN60601-1-2 (62-5) - EN ISO 10079-1		
EXPECTED LIFE OF THE DEVICE	5 YEARS		

05**TRANSPORT, PACKAGING, INSTALLATION**

The SURGICAL SUCTION UNIT WITH TROLLEY weighs 16/19/20 kg, the special handle is the only means to raise and move the machine.

Protected by a cellophane bag, placed in a foam polystyrene cover and inserted in a corrugated paper sealed with adhesive tape.

During the transport shakes and crashes must be avoided.

INSTALLATION:

1. Remove the packaging, make sure that the parts are whole.
2. Insert the handle in the bushes and close trough the two screw with screwdriver.
3. Set the filter and silicon pipes (to see paragraph 07).
4. Connect supply cables to the network and to the device, set to ON the general switch, verify if the pump works.

5. We don't supply suction catheters, if necessary, use only catheters according to norm ISO 8836: 1997

WARNING:

The packaging must be kept during the warranty time.

Do not disperse the packaging in the environment, recover the recyclable products and hand the others over to a company of rubbish disposal.

The packed machines must be kept in dry places and at a temperature between -10°C and +40°C.

**5.1
DEMOLITION**

When the device is ready to be thrown away, you have to respect the rules for the disposal of special rubbish in accordance with the current laws of destination countries for special waste disposal and hospital.

The Surgical Suction Unit with trolley in its different versions is suitable for the aspiration of body and extr-body liquids.

They've been conceived for displacements in hospital wards.

The small dimensions make this device versatile, handy and easy to move through three wheels with pads and a chromium-plated handle.

With the cable supplied with all the versions it is possible to feed the suction unit with electricity network 230V ~ 50 Hz.

The light on the switch indicates the switching on of the device.

The good working of the device depends on the correct use, maintenance and cleaning.

The improper use causes bad working, damages to the device and to the environment.

The position of work must allow to reach the control panel and to have a good view of the vacuum indicator, of the pot and filter.

We suggest to place the suction unit on a flat and solid surface in order to have easy movement and the complete volume of use of the pot and the major efficacy of the non return valve system.

The suction stops when the non return valve system starts to work. Switch off the system and clean the pot (see page 8).

ADJUSTMENT OF THE SUCTION UNIT WITH 2 L POTS TO 4 L POTS

The versions with 2 pots of 2 l. capacity have two bushes that reduce the dimension of the basic seat (for 4 l. pot) to receive 2 l. pot that has a smaller diameter.

These bushes are fixed on the seat basis (for 4 l. pot) with a screw placed in the middle of the bush.

In case of modifications please contact TECNO-GAZ.

Room temperature 5 – 35 °C

Remove the equipment and make ordinary/extraordinary maintenance in case liquids are sucked by the pump (Suction Rubber holder pad.12 pos. B) inside the unit.

We recommend a check up in TECNO-GAZ.

The correct use of the device is granted by the correct working of the vacuum meter and by an effective aspiration.

In case of bad working or for a check on parameters please contact Tecno-gaz.

INSTRUCTION FOR USE

Connect the alimentation cable supplied with the suction unit in the rear side.
 Connect the same cable to the electricity network 230 V AC.
 Operate on the switch on control panel to start the aspiration.


Follow carefully the instructions for tubes connections and insertion of filter described in the following page.

DESCRIPTION OF CONNECTIONS:

Important for a correct use of the device and to preserve the device to wire the external tubes of the suction unit following the indications.

Connect the tubes as described:

1. Fix the tube (A) to the tube connector of aspiration (B)
2. Insert the hydrophobic filter on the free side of the tube (A)
3. Fix the tube(C) to the free extremity of the filter (F)
4. Fix the tube (C) to the entering connection of the pot (E) (vacuum), characterized by the presence of a float to close the aspiration with full pot. This safety system allows to protect the filter from liquids (F).
5. Connect the tube (D) to the free connection (G) (patient) of the collection pot.
6. The final part (i.e. catheter) is connected to the tube (D) trough a connection (NOT INCLUDED)

IMP: the filter has to be mounted in vertical position above the connection of aspiration (vacuum)
FILTER

The filter has a great importance, it is protected from the liquids coming from the pot by a safety valve, but not from foam and steam that has to keep.

For this reason we suggest to replace it frequently and check its functionality.

- Switch on the suction unit with free flow and check that vacuum-meter doesn't move from ~0 if the degree of vacuum increases, the filter should be replaced
- A fall of the aspiration could indicate the necessity to replace the filter.
- **TO REVERSE THE WORKING OF THE FILTER CAUSES POLLUTION AND ENDANGERS THE FUNCTIONING OF THE SUCTION UNIT**
- **RESPECT THE CONNECTION OF TUBES, FILTERS AND TUBECONNECTOR**

VACUUM REGULATOR

A knob on the control board allows to adjust the max. vacuum level applied to the patient.

To check the correct working of the pump and of the vacuum regulator:

- Close the ending part and take the regulator to the highest degree of vacuum.
- Read the value on the vacuum indicator (~ -0.79 bar)
- Turn the regulator and check the correspondence of the values indicated by the index of the hand-grip and by the vacuum-meter.
- What is above mentioned refers to an efficient circuit.

THE ORDINARY MAINTENANCE CONSISTS FIRST OF ALL IN CLEANING:

Cleaning must be careful and done after every use

Important:

Cleaning the accessories and the device allows a correct working of the machine and the protection of both the patient's and the operator's health

INSTRUCTION FOR CLEANING AND DISINFECTION OF THE UNIT:

- Wear latex gloves

External:

- Remove all the accessories and clean the aspirator with a wet cloth with a not aggressive detergent

Pots:

- Empty the pot out (the pot's cleaning can be facilitated by putting some water inside: the aspired diluted liquid will be more easily eliminated)
- Disassemble every part,(the valve "too full" must be correctly reassembled, the float must have the gasket turned upward and the wrapping must be well inserted)

Tubes:

- Wash by running water
- Check that tubes don't present dirty parts with presence of corpuscles
- The tubes can be steam sterilized at 134°C, cold disinfected, or following the sterilization protocol of the site where the unit is used

Filter replacement:

- For a correct working, check the filter and replace it when necessary, in order to avoid every inner pollution which would make it unusable.
- The filter cannot be cleaned, it must be replaced when it is choked
- The accessories (pot, filter...) must be supplied by Tecno-gaz
- Temperature of stockage of the filter +5°C/ +40°C

THE EXTRAORDINARY MAINTENANCE IS NOT EXPECTED. IN CASE OF ANOMALIES CONTACT TECNO-GAZ FOR TECHNICAL ASSISTANCE AND REPARATION

WARNINGS FOR SHIPMENT:

The unit must be **SANITIZED** with particular attention to parts returned for assistance.

The package must be suitable, the best think is preserve the original package.

09

ANALYSIS AND RESEARCH OF FAILURES

N° FAULT	CAUSE	ASSISTANCE
1	IT DOESN'T SUCK	THE POT IS BADLY SCREWED DOWN UNSCREW AND SCREW THE TOP, CHECK GASKET'S INTEGRITY AND RIGHT POSITION
		EXTERNAL PIPES WITH FOLDING REMOVE NARROWING ON THE PIPES
		FILTER IS OBSTRUCTED FILTER MUST BE REPLACED
2	THE FLOAT DOESN'T	WRONG ASSEMBLING CHECK THE CORRECT POSITION OF NON RETURN VALVE
3	NOT WORKING	EXTERNAL FUSE INTERRUPTED  IN CASE OF NOT WORKING, ADDRESS TO TECNO-GAZ

10

COMMERCIAL DATA, SPARE PARTS AND DOCUMENTATION

The positions refer to the electro-pneumatic outline (following chapter).

DESCRIPTION	REFERENCE
POWER SUPPLY CABLE	CE06006
KIT OF PIPES FOR EXTERNAL WIRING	SB1A138
HYDROFOBIC FILTER	CPF0029
POT WITH NON RETURN VALVE 2 L	CM84151
POT WITH NON RETURN VALVE 4 L	CM84152

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