GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Monza, 102 - 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 - fax 199 400 403
Export: tel. +39 02 953854209/221/225 fax +39 02 95380056
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Lampada frontale WIDA
WIDA head light
Lampe frontale WIDA
Lámpa frontal WIDA

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FEATURES

The Wida head light (cod. 30850-30854) is an optical-fibre lighting device suitable for diagnostic examinations and minor surgery. The product shall be used by qualified and trained medical staff only.

The examination lamp Wida is assembled on a practical helmet, adjustable and adaptable to any kind of head and it is equipped with a high-flexible optical-fibre cable, assuring an optimal transmission of light.

The lamp is supplied with an adapter (cod. 30830) for the connection to GIMA sources of light (cod. 30801-30802), however adapters are available for spotlights of the most common brands. By means of a revolving ring located on the headlight it is possible to adjust focusing and therefore the diameter of the light beam.

Focus spot diameter WIDA 30850:

7 - 15 cm at a distance of 20 cm

10 - 22 cm at a distance of 30 cm

16 - 35 cm at a distance of 50 cm

Focus spot diameter WIDA PLUS 30854:

6 - 10 cm at a distance of 20 cm

9 - 15 cm at a distance of 30 cm

14 - 20 cm at a distance of 50 cm

The product or its components cannot be used for purposes different from the one specified in the present manual.

PRESCRIPTIONS

The optical-fibre cable is very delicate. It is advisable to carefully handle it. Avoid bending the cable or putting it in contact with acids or abrasive materials.

The product shall be used by qualified personnel only.

Do not use the equipment in case it is damaged. Apply to your retailer. Avoid precarious repairs. Repairs shall be carried out with original spare parts only, which shall be installed according to the intended use.

Since the product is made of corrosion-proof materials suitable for the environmental conditions foreseen for its normal use, does not require special care, however it is necessary to store it in a closed place making sure that is protected from dust and dirt to assure its hygienic conditions. Moreover, it is recommended to store the product in a place which can be reached easily by the personnel in case of necessity.

UNPACKING

Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled.

They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

INSTALLATION

Before using the product, it is possible to adjust the helmet housing the lamp Wida to the head size. To do this, make use of the apposite adjusting screws, located on the back and/or on the top (see picture) allowing tightening or loosening the headbands.

In case several spotlights are available, it is necessary to choose the most suitable optical-fibre cable terminal (adapter). It is sufficient to screw it in the apposite seat at the cable end.

FUNCTIONING

Once the lamp Wida has been firmly positioned on your head, to turn it on it is sufficient to insert the cable terminal into the socket of the source of light and turn it on. Then, it is possible to adjust the intensity of light according to your requirements.

It is also possible to direct the light beam by lifting and lowering the lighting unit. Once the most suitable position has been found, fix the unit by means of the apposite lock nut (see picture).

As already mentioned, it is possible to adjust the diameter of the light beam of the examination lamp Wida.

Clockwise rotate the headlight ring to reduce the diameter of the light beam and counter-clockwise rotate it to widen the diameter.

MAINTENANCE

Before the use, it is advisable to check that the optical-fibre is not damaged. In case it is damaged, remove and replace it. To this purpose, proceed as follows:

- Open all the automatic snap fasteners on the helmet and release the cable
- Unscrew the lock nut on the front and withdraw the lighting unit
- Clockwise rotate the fastening joint (cod. 30850)
- Counter-clockwise unscrew the lighting unit from the cable, until it can be drawn out
- Insert the new cable and screw it in the lighting unit
- Counter-clockwise turn the fastening joint to replace it in the original position (cod. 30850)
- Reassemble the fastening joint and tighten the lock nut
- Lock the optical-fibre cable by means of the apposite snap fasteners

The examination lamp Wida, cable included, can be cold cleaned with disinfectants and a sponge. Periodically clean the cable end to avoid dirt deposits.

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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GIMA S.p.A.