ELECTRODOS DE pH NAHITA pH ELECTRODE

- ◆ The electrode must be used with a high input impedance (= 10¹²) pH-meter.
- ◆ Do not submerge the electrode for a long period of time in distilled water, protein solutions or hydrofluoric acid. The electrode must not come into contact with silicone organic grease.
- ◆ The electrode slope value decreases with use and that increases time response and decreases sensibility of the electrode. Dip the electrode in 4% hydrofluoric acid solution for 3-5 seconds; then, clean it with distilled water and submerge it in KCl for approximately 1 hour before using it again.
- ◆ If the electrode bulb is contaminated or the diaphragm blocked, sensibility of the electrode will be affected. Clean the electrode with and adequate detergent. Detergents to be used for electrode cleaning:

CONTAMINANT	DETERGENT
Inorganic substances	HCl < 1 M
Organic substances	Weak base
Grease and oil	Mild detergent or methanol
Resin	Alcohol, acetone, ether
Protein deposits	Acidified enzymatic solution (1% pepsin in HCl 0.1 M)
Dye	H ₂ O ₂ diluted bleaching agent

NOTE: WITH THE AIM OF AVOIDING A WRONG OPERATION OF THE EQUIPMENT, ACCIDENTS, ETC., READ THE INSTRUCTIONS MANUAL OF THE PH-METER.

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1. APPLICATIONS

Nahita pH electrodes are combined electrodes characterized by including in the same body a reference and an indicator electrode, and are used for pH determination in several applications in chemistry pharmaceuticals, education, etc.

2. TECHNICAL SPECIFICATIONS

Features in common

BNC connector
 Reference system: Ag/AgCl
 pH measuring range: 0-14

Specific features

PLASTIC BODY

Model	Dimensions (ØxL)	Working temperature	Diaphragm	Membrane	Electrolyte
218	12x165 mm	5-60 °C	Fibre	Bulb	Gel, sealed
219	12x165 mm	5-60 °C ^o	Fibre	Bulb	Liquid, refillable
900	12x165 mm	5-60 °C	PTFE ring	Bulb	Gel, sealed
901	12x165 mm	5-60 °C	PTFE ring	Flat	Gel, sealed

GLASS BODY

Model	Dimensions (ØxL)	Working temperature	Diaphragm	Membrane	Electrolyte
651	12x148 mm	5-60 °C	Ceramic	Bulb	Liquid, refillable
652	12x145 mm	5-90 °C	Ceramic	Bulb	Liquid, refillable
953	12x160 mm	5-60 °C	Ground glass	Bulb	Gel, sealed
955	12x160 mm	5-90 °C	Ground glass	Bulb	Liquid, refillable

3. OPERATION AND MAINTENANCE

- ◆ Before measurement, both the electrode and the pH-meter must be calibrated with a standard buffer solution.
- ◆ The electrode bulb must not come into contact with hard objects or substances. Any scratches or cracks on the glass membrane will make the electrode useless.
- ◆ Once measurement has finished, put the bulb into the cap that should contain a KCl solution.
- ◆ Refillable electrodes present a 3 M KCl solution as a reference electrolyte. If the electrolyte level is 2.5 cm under the filling hole, the electrode must be refilled. Remove the small cap located at the upper part of the electrode and add new 3 M KCl solution.
- ◆ The connector must be kept always clean and dry.