LEVITAS®

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CE I Class Medical Device
according to Directive 93/42/EEC and further modifications

1. CODES

- LAD462** Kit mattress and pump with pressure adjustment Piuma UP (bubble pad)
LAD463 Kit mattress and pump with pressure adjustment Piuma UP (single cells pad)

2. INTRODUCTION

Thank you for purchasing a LEVITAS mattress and pump Piuma UP Kit (single cell pad) by Moretti. The pump and mattress kits by Moretti are designed and manufactured to meet all your needs for a practical, correct and safe use. This manual contains small suggestions for proper use of the device you have chosen and valuable advice for your security. You may read the complete manual before using the kit pump and mattress, in case of doubts please contact your dealer, who will help and advise you properly.

NOTE, CAUTION AND WARNING STATEMENTS:

NOTE: Indicate some tips.

CAUTION: Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

WARNING: Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

3. INTENDED USE

The kit pump and mattress have been designed for the bedsore treatment and prevention, which may occur during long hospital and/or house confinement.

CAUTION!

- Do not use the product for a purpose not indicated in this manual
- Only qualified personnel or trained for the treatment and prevention of bedsore, can operate with the device
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.

4. EC DECLARATION OF CONFORMITY

Moretti S.p.A. hereby declares, under its own exclusive responsibility, that the products manufactured and sold by Moretti S.p.A. and belonging to the kit pump and mattress category are in compliance with the applicable provisions outlined in the Medical Devices Directive 93/42/EEC and further modifications.

To this end Moretti S.p.A. guarantees and declares as follows, under its own responsibility:

1. The devices in question satisfy the essential requirements of Annex I of directive 93/42/EEC as prescribed by Annex VII of the above directive
2. The devices in question ARE NOT MEASUREMENT DEVICES
3. The devices in question ARE NOT INTENDED FOR CLINICAL INVESTIGATION
4. The devices in question are sold in NON STERILE packing
5. The devices in question should be considered as belonging to class I
6. DO NOT use/install the devices in question for any other purpose than that expressly indicated by Moretti S.p.A.
7. Moretti S.p.A. maintains and submits the documentation certifying conformity with Directive 93/42/EEC to the proper authorities and further modifications
8. The devices mentioned are manufactured using materials in accordance with Directive 2011/65/UE - ROHS2

5. STANDARDS

The system has been tested and successfully approved to the following standards: EN 60601-1:2005+A1:2012

EN 60601-1-2:2007

IEC 60601-1-11:2010

IEC 60601-1-6:2010+A1:2013

IEC 62366:2007+A1:2014

EN ISO 13485:2012+AC:2012

EN ISO 14971:2012

EN ISO 15223-1:2012

EN 1041:2008

EN 71-3:2013+A1:2014

6. GENERAL WARNINGS

- For proper use of the device, refer to the following manual carefully
- Advise with physician or therapist to determinate the correct adjustment and the correct use of the device
- Keep the packed kit away from heat sources
- SERVICE LIFE- the device use limit is defined by the wear of the parts
- DO NOT allow children to play on or operate the pump

6.1 Electromagnetic compatibility warnings

This equipment has been tested and it complies with medical devices limits belonging to EN 60601-1-2:2007. These limits have been designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guaranteed that the interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference following one or more of the suggestions below:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other devices are connected
- Consult the manufacturer or field service technician for help

7. SYMBOL DEFINITION

 Code	 Attention, should read the instructions	 Temperature Limitation	 BF type	 Follow instructions for use	 Class II Equipment	
 Production dates	 Waste Electrical & Electronic Equipment (WEEE)	 ON	 OFF	 Manufacturer	 Lot	 Protected against solid foreign objects of 12.5 mm and greater. Protection against vertically falling water drops

8. KIT CONTENTS


The Piuma UP series is composed by 2 kind of kits which can satisfy all the requirement about treatment and prevention of complications caused by long hospital and/or house confinement, in particular bedsores.

LAD462 Kit mattress and pump with pressure adjustment Piuma UP (bubble pad)

LAD463 Kit mattress and pump with pressure adjustment Piuma UP (single cells pad)

9 GENERAL DESCRIPTION

9.1 Pump unit

9.2 Front panel

9.2.1 Pressure Adjust Knob

Pressure adjust knob controls the air pressure output. Please consult the physician for a suitable setting

9.2.2 Main Power Switch

To turn the pump unit on/off

10. INSTALLATION

Unpack the box and check the package contents for completeness. If there are any damages, please immediately contact your retailer.

1. Place the mattress on a common mattress, with the air hose connectors oriented to the footboard

ATTENTION!

We recommend to place the mattress with interchangeable elements over a mattress with a height of at least 5 cm

2. Hang the pump onto the bed rail (footboard side), or place it on a flat surface
3. Connect air hose connectors from air mattress to the pump unit

ATTENTION!

Ensure the air hoses are not kinked or tucked under mattress

4. Plug the power cord into the electrical outlet

-  **NOTE:** 1. Make sure the pump unit is suitable for the local power voltage
2. The plug is also served to disconnect the device. Do not position the equipment so that it is difficult to disconnect the device

ATTENTION!

The pump should only be used with mattress recommended by the manufacturer. Do not use it for any other purpose

5. Turn the switch ON

Several installation tips are listed below:

After installation, the extra length of the power cord, if any, should be neatly arranged to avoid any tripping accidents. The EQUIPMENT should be firmly placed at position where users/doctors can access easily.

11. BEFORE USE

Check the correct assembly, paying particular attention to air hoses connector.

12. USING WARNINGS

DANGER - To reduce risk of electrocutione

1. Always unplug this product immediately after using
2. Do not use while bathing
3. Do not place or store this product where it can fall or be pulled into a tub or sink

4. Do not place in or drop into water or other liquid
5. Do not reach for a product that has fallen into water. Unplug immediately

WARNING - To reduce risk of burns, electrocution, fire or injury to people


1. Evaluate patients for entrapment risk according to facility protocol and monitor patient appropriately
2. This system is not intended to be used with patients who have spinal cord injury.
3. Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may occur if a child swallow a small part detached from the device
4. Use this product only for its intended use as described in this manual. Do not use attachment not recommended by the manufacturer
5. Never operate this product if it has a damaged cord or plug, if it is not working properly, or if it has been dropped, damaged, or dropped into water. Return the product to a service centre for examination and repair
6. Keep the cord away from heated surfaces
7. Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where their openings may be blocked. Keep the air openings free of lint, hairs, and other similar particles
8. Never drop or insert any object into any opening or hose
9. Do not modify this equipment without authorization of the manufacturer
10. Mattress cover has passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately
11. If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance to 3,3m between devices or turn off the mobile phone

13. OPERATION


 **NOTE:** Always read the operating instruction before use

13.1 General operation


1. Open the package: it is recommended to use a cutter with short blade or a tool not to damage the mattress.
2. The overlay mattress is rolled up: take it out of the box paying attention to the compressor. Unroll it on the bed mattress.
3. Fold and put the flaps of the overlay mattress between the mattress and the metal net, to secure it and prevent unintentional movements.
4. Take out and unpack the compressor: then hang it to the footboard of the bed by the hooks on the back of the compressor, or place it on a flat surface.

5. Make sure that the switch "ON/OFF" of the compressor is OFF.
6. Connect the tubes placed at the bottom of the overlay mattress to the compressor, checking they are connected correctly.

8. Turn the knob clockwise, putting the narrow's knob in proximity of "+" symbol. In this way the compressor will inflate the mattress more quickly.
9. Connect the compressor plug to the 220V outlet.
10. Switch on the compressor with the "ON/OFF" switch and wait about 25 minutes, so that the compressor inflate the mattress completely.

⚠ NOTE: every time the mattress is set up for use, it is recommended that the pressure first to be set to the max. The user/career can then adjust air mattress weight levels to the desired softness after set up has been completed

14. CLEANING

It is important to follow the cleaning procedures before using the equipment on human bodies; otherwise, patients and/or doctors may have the possibility of getting infection.

⚠ CAUTION: Do not immerse or soak pump unit.

Wipe the pump unit with a damp cloth and a mild detergent. If other detergent is used, choose one that will have no chemical effects on the surface of the plastic case of the pump unit. Wipe down the mattress with warm water containing a mild detergent. The cover may also be cleaned by using sodium hypochlorite diluted in water. All parts should be air dried thoroughly before use.

⚠ CAUTION: Do not use phenolic based product for cleaning.

⚠ CAUTION: After cleaning, dry the mattress without direct exposure sunlight.

WARNING!

Unplug the pump before starting the cleaning procedure

15. STORAGE

- Fold the mattress accurately to prevent creases
- Roll from the head end towards the foot end with the valve open

For single cells mattress:

1. Lay the bubble pad or mattress on a flat surface and upside down
2. Roll-up the mattress from the head end towards the foot end
3. Foot-end strap can then be stretched around the rolled pad/ mattress to prevent unrolling

 CAUTION: Do not fold, crease or stack the mattresses.

16. MAINTENANCE

The LEVITAS devices by MORETTI are checked carefully and supplied with EC mark, once launched on the market.

The products are intended to offer safe and reliable operation when used or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function. Otherwise, service and inspection of the devices generally should not be required.

1. Check main power cord and do not plug it if there is an abrasion or excessive wear
2. Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are stubbed together correctly
3. Check the airflow from the air hose. The airflow should alternate between each connector every half-cycle time
4. Check the air hoses if there is kink or breaks. For replacement, please contact Apex Medical Corp. or your suppliers

17. CONDITIONS OF DISPOSAL

General conditions of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the thermometer in the appropriate disposal areas for recycling.

 **Information for the correct disposal of the product in accordance with European Directive 2002/96/CE:**

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin the product bears the barred bin symbol.

18. SPARE PARTS/ACCESSORIES

For spare parts and accessories refer to main catalogue

19. TROUBLESHOOTING

If your questions can't be answered with above information, please contact your local agent directly. They might require a technician to take care of the problem

PROBLEM	SOLUTION
Power is not ON	<ul style="list-style-type: none"> • Check if the plug is connected to mains • Turn "ON" the compressor switch
The mattress doesn't inflate	<ul style="list-style-type: none"> • Check if the connection between mattress and pump unit is tightly secured • Check if there is any leakage on air cells
Mattress form is loose	<ul style="list-style-type: none"> • Check if all the snap buttons or straps of mattress are all securely fastened • Check if the mattress is fixed to the bed frame by straps
No air produced from some air outlets of the air tube connector	<ul style="list-style-type: none"> • This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time

⚠ NOTE: If the pressure level is consistently low, check for any leakage (tubes or air hoses). If necessary, replace any damaged tubes or hoses or contact your local qualified dealer for repair.

20. GUIDANCE AND MANUFACTURER'S DECLARATION

20.1 Guidance and Manufacturer's Declaration - Electromagnetic Emissions

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
RF emissions CISPR 11	Class B	
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	

20.2 Guidance and Manufacturer's Declaration - Electromagnetic Immunity

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	± 2kV for power supply line ± 1kV for input/out line	± 2kV for power supply line ± 1kV for input/out line	Mains power quality should be that of atypical commercial or hospital environment

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment - Guidance
Surge IEC61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of atypical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment

 NOTE: UT is the A.C. mains voltage prior to the application of the test level

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside ISM bands	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \frac{1.5 \sqrt{P}}{f}$ 150kHz to 80MHz $d = \frac{1.5 \sqrt{P}}{f}$ 150kHz to 80MHz 80 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

 NOTE: At 80 MHz and 800 MHz, the higher frequency range applies

 NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

 NOTE: The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land / mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

Recommended separation distances between portable and mobile RF communications equipment and this device

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2,5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

 NOTE: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

 NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people

21. TECHNICAL FEATURES

21.1 Pump

MODEL	LAD462	LAD463
Power supply	~ 220V 50/60Hz	~ 220V 50/60Hz
Fuse T1AL	250V	250V
Power cord	2 mt	2 mt
Pressure output	40~100mmHg	40~100mmHg
Air release	6/7LPM	6/7LPM
Cycle time	6 min	12 min
Noise	32-34 dD	32-34 dD
Dimensions	27x13x10,5 cm	27x13x10,5 cm
Weight	~ 1,1 Kg	~ 1,1 Kg
Environment temperature	Operation: +10°C/+40°C Storage: -18°C/+43°C Shipping: -18°C/+43°C	Operation: +10°C/+40°C Storage: -18°C/+43°C Shipping: -18°C/+43°C
Environment humidity	Operation: 10%/75% non condensing Storage: 10%/95% non condensing Shipping: 10%/95% non condensing	Operation: 10%/75% non condensing Storage: 10%/95% non condensing Shipping: 10%/95% non condensing
Atmospheric pressure	70/101,325 KPa	70/101,325 KPa
Classification	Class II, Type BF, IP21 Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)	Class II, Type BF, IP21 Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)

21.2 Mattress

MATTRESS	MATTRESS BUBBLE PAD (LAD662)	ELEMENT MATTRESS (LAD663)
Model	Bubble pad	Overlay mattress
Components	130 bubbles	18 cells
Material	Medical-grade PVC	Medical-grade PVC Cell: Nylon/PVC
Dimension	~ 200x90x7,6 cm	~ 190x80x11 cm
	-	~ 10,2 cm
Weight	2,3 Kg	3,8 Kg
Max user's weight	135 Kg	135 Kg
Thickness	0,3	-

22. WARRANTY

Moretti products are guaranteed from material or manufacturing faults for a period of (two) 2 years from the purchasing date, subject to any limitations and exclusions specified below. The warranty does not apply in case of damages caused by improper use, abuse or modification of the product, and the warranty will not be valid if the instructions for use are not strictly followed. The correct intended use of the device is specified in this manual.

Moretti is not responsible for damages resulting from personal injury or otherwise caused by or in relation to wrong installation or improper use. Moretti does not guarantee the products against damage or defects in the following conditions: natural disaster, not authorized maintenance or repair, damage due to power supply problems (where provided) use of spare parts not supplied by Moretti, improper use, not authorized alteration, shipping damage (other than original shipment by Moretti), or in case of insufficient maintenance as indicated in the manual. The warranty doesn't cover components subject to wear and tear during the correct use of the device.

23. REPAIR

23.1 Warranty repairs

In the event that a product Moretti to be defective in material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, in its sole discretion, can repair or replace the item under warranty at a Moretti dealer or at Moretti location . The labor cost for the repair may be charged to Moretti if the repair is covered by warranty. Repair or replacement does not renew or extend the warranty

23.2 Repair not covered by warranty

A product not covered by warranty can be returned for repairing only after receiving prior authorization from Moretti customer service. The labor and shipping costs for a repair not covered by warranty will be fully borne by the customer or dealer. Repairs of products not covered by warranty are warranted for six (6) months from the date you receive the product repaired.

23.3 Non-defective device

The customer will be informed if, after having examined and tested returned product, Moretti claims that the product is not defective. In this case the good will be returned to the customer, all shipping costs will be paid by custome.

24. SPARE PARTS

For spare parts only refer to MORETTI master catalogue.

25. EXEMPT CLAUSES

Except as expressly specified in this warranty and within the limits of the law, Moretti does not offer any representation, warranty or condition, express or implied, including any representations, warranties or conditions of merchantability, fitness for a particular purpose, non-infringement and interference. Moretti does not warrant that the use of the product Moretti will be uninterrupted or error free. The duration of any implied warranties that may be imposed by law are limited to the warranty period, within the limits of the law. Some states do not allow limitations on how long an implied warranty lasts or the exclusion or limitation of incidental or consequential damages in connection with products for consumers. In such states or countries, certain disclaimers and limitations of this warranty may not apply to you. This warranty is subject to change without notice.

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Acquistato in data / Purchasing date _____

Rivenditore / Authorized dealer _____

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Venduto a / Purchased by _____

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 **MORETTI S.P.A.**

Via Bruxelles, 3 - Melegnano 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com email: info@morettispa.com

MADE IN P.R.C.

Via Bruxelles, 3 - Melegnano
52022 Cavriglia (Arezzo)
Tel. +39 055 96 21 11

Fax. +39 055 96 21 200
www.morettispa.com
info@morettispa.com