
SINCE 1966
MADE IN ITALY

Aspira

by **FLAEM**®

IT> ISTRUZIONI PER L'USO - Aspiratore professionale medico chirurgico
EN> INSTRUCTIONS FOR USE - Professional medical surgical suction machine
FR> INSTRUCTIONS POUR L'UTILISATION - Aspirateur professionnel médico - chirurgical
NL> GEBRUIKSIINSTRUCTIES - Professioneel medisch chirurgisch afzuigapparaat
DE> GEBRAUCHSANWEISUNG - Professionelles medizinisches chirurgisches Absauggerät
ES> INSTRUCCIONES PARA EL USO - Aspirador profesional médico quirúrgico
RU> ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ - Профессиональный медицинский хирургический аспиратор
PL> INSTRUKCJA OBSŁUGI - Profesjonalny, chirurgiczny aspirator przenośny
EL> ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - Επαγγελματικός ιατρικός αναρροφητήρας χειρουργικής χρήσης
<AR> تعليمات الاستخدام - شفط مهني طبي جراحي

SCHEMA DI COLLEGAMENTO PER NOMENCLATURA - CONNECTION DIAGRAM FOR NOMENCLATURE
SCHÉMA DE CONNEXION POUR NOMENCLATURE - NOMENCLATUUR VERBINDINGSSHEMA
SCHALTPLAN FÜR DIE BEZEICHNUNGEN - ESQUEMA DE CONEXIÓN PARA NOMENCLATURA
СХЕМА СОЕДИНЕНИЯ СОСТАВНЫХ ЧАСТЕЙ - SCHEMAT POŁĄCZEŃ DLA NOMENKLATURY
ΣΧΕΔΙΟ ΣΥΝΔΕΣΗΣ ΚΑΙ ΟΝΟΜΑΤΟΛΟΓΙΑ - مخطط توصيل خاص بالتسميات - 82

ITALIANO
Pag. 1

ENGLISH
Pag. 10

FRANÇAIS
Pag. 19

NEDERLANDS
Pag. 28

DEUTSCH
Pag. 37

ESPAÑOL
Pag. 46

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min. 690 hPa; max. 1060 hPa

Mod. P1211EM/11 I - Mod. P1211EM/13.5 I - Mod. P1211EM/20 I

We are pleased you have purchased our product and we thank you for your trust in us.
It is very important that the health care worker, the person trained in use or the patient read and understand the information for use and maintenance.

INTENDED USE

ASPIRA is an aspirator designed for medical and surgical use in healthcare centres such as hospitals, for home care and for pharyngeal suction for permanent tracheotomy. The device generates a negative pressure (suction) which allows fluids to be extracted through a single use tube connected to a collection canister which retains the fluids until they can be properly disposed of. Use of the device must be prescribed by a doctor.

PRODUCT DESCRIPTION

ASPIRA (1) is equipped with vacuum level adjustment knob (4), vacuum gauge (3) and 1000 ml or 800 ml canister (9) with overflow safety device (10) to prevent fluids from entering the suction pump. This device interrupts the suction flow via a float-controlled valve housed in the canister lid. It does not require lubrication, is manageable, simple to use, reliable, resistant and noiseless. ASPIRA is equipped with the following FLAEM accessories: 1000 ml or 800 ml canister (9) with overflow safety device (10), connection tubing (6/8/14), single use antibacterial/antiviral hydrophobic filter (7), single use manual vacuum flow control (15).

NOTE: Only use original FLAEM accessories and spare parts, we disclaim any liability in the event of using non original spare parts or accessories.

DEVICE NOMENCLATURE

1) SUCTION MACHINE	10) SAFETY DEVICE
2) ON/OFF SWITCH	11) CANISTER LID
3) VACUUM GAUGE	12) CANISTER "VACUUM" INLET
4) VACUUM FLOW ADJUSTMENT KNOB	13) CANISTER "PATIENT" INLET
5) AIR INLET	14) CONNECTION TUBING
6) CONNECTION TUBING	15) MANUAL VACUUM FLOW CONTROL
7) HYDROPHOBIC ANTIBACTERIAL/ANTIVIRAL FILTER (SINGLE USE)	16) ASPIRATOR TUBE (NOT SUPPLIED WITH THE DEVICE). This tube is not supplied and is not available as a spare part for the medical device. Any tubes to be used with our aspirators must feature a shape that can adapt to the fitting of the "MANUAL VACUUM FLOW CONTROL" supplied with the aspirator.
8) CONNECTION TUBING	17) POWER CABLE
9) COLLECTION CANISTER - ALTERNATIVE STERILISABLE MODEL - SINGLE PATIENT MODEL	

IMPORTANT NOTICES

The manufacturer makes every effort to ensure that every product is of the highest quality and safety; however, as for any electrical device, basic safety regulations must always be followed in order to avoid harming persons and things.

- **The medical device must NOT be used in the operating theatre, for drainage or for thoracic drainage or on vehicles such as ambulances or emergency vehicles.**
- The Medical Device is not suitable for use in MRI premises (Magnetic Resonance Imaging)
- Carefully follow the user and maintenance instructions to assure proper operation and to extend the service life of the device.
- The supplied cable is equipped with a protective gasket to prevent the penetration of liquids into the device. Do not separate this gasket from the cable.
- Before using the product for the first time, and periodically during its lifetime, check the power cable to make sure it is not damaged. In the event it should be damaged, do not plug in the cable and immediately bring the product to an authorised FLAEM service centre or to your trusted dealer.
- In the presence of children and non self-sufficient individuals, the device must be used under the close supervision of an adult who has read this manual.
- Some parts of the device are small enough to be swallowed by children; keep the device out of the reach of children.
- Do not use the provided tubing and cables for any other purpose than those specified, as they can cause risk of strangulation. Be particularly attentive with children and individuals with disabilities because they are often unable to correctly assess risk.
- The suction machine is intended exclusively for the collection of NON-flammable fluids. It is not suitable for use in the presence of a flammable anaesthetic mixture with air, or with oxygen or nitrous oxide.
- Always keep the power cable away from hot surfaces.
- Keep the power cable away from animals (for example, rodents) which could damage the insulation.
- Do not handle the plug with wet hands. Do not use the device in damp conditions (for example when taking a bath or shower). Do not immerse the device in water; in the event of immersion immediately disconnect the plug. Do not remove or touch the immersed device; unplug the power cable first. Immediately bring the device to an authorised

FLAEM service centre or to your trusted dealer.

- Use the device only in dust-free conditions, otherwise treatment could be compromised.
- The device casing is not waterproof. Do not wash the device under running water or by immersion and keep it safe from being sprayed by water or other liquids.
- Do not expose the device to particularly extreme temperatures.
- Keep the device away from sources of heat, direct sunlight or excessively warm locations.
- Never obstruct the air inlets positioned on both sides of the device.
- Always operate the device on an unencumbered hard surface.
- Always unplug the power cable immediately after use.
- The device is equipped with an easy to inspect safety fuse in case of malfunction. Before doing so, unplug the power cable. Repairs may only be performed by authorised FLAEM personnel following the information in the restricted area of the website www.flaem.it. Unauthorised repairs void the warranty and may pose a danger for the user.

WARNING: Do not modify this device without authorisation from the manufacturer.

- The Manufacturer, the Dealer and the Importer shall be held responsible for safety, reliability and performance only if: a) the device is used in compliance with the instructions for use b) the wiring where the device is being used is in compliance with safety regulations and current laws.
- The Manufacturer must be contacted about any problems and/or unexpected events concerning operation and if any clarifications on use, maintenance or cleaning are needed.
- Should your device fail to provide the expected performance, contact the authorised Flaem service centre for clarifications.
- Ensure the connections and the canister lid have been carefully sealed in order to prevent loss of suction.
- Do not overturn the canister while connected to the operating device as the liquid may be aspirated into the device, thus damaging the pump. Should this occur, immediately switch off the suction machine, then empty and clean the canister (send it to an authorised FLAEM service centre).
- When the safety device (10) is triggered, suction is halted. Switch off the device, empty the canister (9) and carry out the cleaning procedure.
- We recommend personal use of the accessories, the collection canisters and the connecting tubing to prevent risk of contagious infection.
- The manual vacuum flow control (15) is a sterile disposable product and must be replaced after each use.
- Check the expiry date on the original package of the manual vacuum flow control (15) and ensure the sterile package is undamaged. Should the package not be intact contact the authorised Flaem service centre or your trusted dealer. FLAEM NUOVA disclaims any liability for harm to the patient related to deterioration of the aforementioned sterile package due to third parties handling the original packaging of the entire device.
- **The hydrophobic/antiviral and antibacterial filter supplied with Flaem aspirators must be strictly replaced for each new patient, or if the said filter is saturated. Should the aspirator be used directly on the patient, the filter for Single-patient use must be replaced within a maximum period of TWO MONTHS, or if the said filter is saturated.**
- The amount of time required between preservation conditions and use is approximately 2 hours.

INSTRUCTIONS FOR USE


Before each use, the accessories must be carefully inspected to ensure the absence of dust, incrustation, clots or liquid substances both inside the connecting tubing and in the canister and its respective lid. Furthermore, they must be cleaned following the instructions rigorously as stated in the "CLEANING, SANITISATION, DISINFECTION, STERILISATION" paragraph. We recommend personal use of the accessories, the collection canisters and the connecting tubing to prevent risk of contagious infection.

Instructions for suction on the patient:

1. Connect the accessories, referring to the "CONNECTION DIAGRAM" on the cover.

ATTENTION: in addition to being antibacterial/antiviral, the hydrophobic filter also acts as secondary line of safety and stops any fluids that the primary safety device in the canister lid accidentally does not stop.

Follow the directions below for proper assembly:

2. Plug the supply cable (17) into a power socket matching the device's voltage specifications. After using the device unplug it, the socket must be placed so as to allow for convenient disconnection from the power mains.
3. Switch on the device using the on/off switch (2), keeping in mind that the device is intended to be used on a horizontal surface in order for the safety device (10) to function correctly and prevent liquids from entering the suction pump. The device must be kept away from walls to allow for adequate engine cooling.
4. For more comfortable aspiration, set the desired vacuum value (bar) using the vacuum adjustment knob(4). Vacuum value is increased by turning the knob clockwise and decreased by turning the knob counter-clockwise; these values can be read on the "vacuum gauge" (3). To accelerate aspiration and simplify cleaning, we recommend putting approximately 400 ml of water in the canister (16).
5. Place your finger on the manual vacuum flow control (15) and, activating it in pulses, begin suction on the patient using the aspirator tip (NOT SUPPLIED WITH THE DEVICE).
6. Once the process is complete, switch off the device.
7. Remove the plug, then empty and clean the canister and connecting tubing.

CLEANING

Switch off the device before any cleaning procedure and unplug the power cable from the socket.


DEVICE CLEANING

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).


COLLECTION CANISTER AND CONNECTING TUBING

- Detach the tube (16) (NOT SUPPLIED WITH THE DEVICE), the manual vacuum flow control (15) and the tubing (14) from the canister (9). Disconnect the tubing (6/8) from both the canister and the filter (7). Remove the canister from its holder, keeping it upright, and empty it (in the WC at home, in the biological waste container in the hospital) and clean it, disassembling it as shown below:


STERILISABLE CANISTER


1. Remove the safety device (.10) from the canister lid (.11)


2. Remove the float (.10a) from the safety device (.10)


3. Remove the silicone gasket (11a) from the canister lid (.11)


SINGLE PATIENT CANISTER


1. Remove the safety device (.10) from the canister lid (.11)

2. Remove the float (.10a) from the safety device (.10)

3. Remove the silicone gasket (11a) from the canister lid (.11)

STERILISABLE CANISTER

SANITISATION

Before and after each use, sanitise the collection canister and the connecting tubing choosing one of the methods described below.

(method A): rinse each part (**6, 8, 9, 10, 10a, 11, 11a, 14**) with hot potable water (approximately 40°C) with gentle dish washing soap (non-abrasive) or in the dishwasher set to a hot cycle

(method B): immerse each part (**6, 8, 9, 10, 10a, 11, 11a, 14**) in a solution of 50% water and 50% white vinegar. Lastly, rinse abundantly with hot potable water (approximately 40°C).

(method C): boil each of the parts (**6, 8, 9, 10, 10a, 11, 11a, 14**) in water for 10 minutes; it is preferable to use de-mineralised or distilled water to prevent calcium deposits.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

DISINFECTION

The accessories which can be disinfected are (**6, 8, 9, 10, 10a, 11, 11a, 14**)

The disinfection procedure described in this paragraph is to be followed before using the accessories. It is effective on the parts that undergo this treatment only if each step is followed and only if the parts to be treated have previously been sanitised.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

STERILISATION

The accessories that can be sterilised are (**6, 8, 9, 10, 10a, 11, 11a, 14**)

The sterilisation procedure described in this paragraph is effective on the parts that undergo this treatment only if every step is followed and only if the parts being treated have previously been sanitised. The procedure is validated in compliance with ISO 17665-1.

Device: Fractionated vacuum overpressure steam steriliser in accordance with EN 13060.

Directions: Wrap every single part to be treated with a sterile barrier system or packaging in accordance with Norm EN 11607. Place the wrapped parts into the steam steriliser, ensuring that the canister (9) is in an upright position.

Following the device instruction manual, start the sterilisation cycle and select a temperature of:

- 134°C for 10 minutes for a maximum of 50 times for SILICONE TUBING (6,8,14) and the GASKET (11a)
- 121°C and a duration of 15 minutes for a maximum of 50 times for the POLYCARBONATE CANISTER (9).

Storage:

After sanitising, disinfecting or sterilising, reassemble the canister and the connecting tubing following the directions provided in the "CONNECTION DIAGRAM".

At the end of each use store the medical device complete with accessories in a dry place away from dust.

- The manual vacuum flow control (15) is a sterile disposable product and must be replaced after each use.
- ATTENTION: The hydrophobic/antiviral and antibacterial filter (7) supplied with Flaem aspirators must strictly be replaced for each new patient, or if said filter is saturated.

Should the aspirator be used on the same patient, the filter approved for Single-patient use must be replaced every TWO MONTHS at most or if the filter is saturated.

The filter cannot be sanitised, disinfected, or sterilised.

It is in any case good practice to sanitise and disinfect the device and the accessories, with the exception of the Filter, always after each use, whether it is for the same patient or for a new patient.

SINGLE PATIENT CANISTER

SANITISATION

Before and after each use, sanitise the collection canister and the connecting tubing choosing one of the methods described below.

(method A): rinse each part (**6, 8, 9, 10, 10a, 11, 11a, 14**) with hot potable water (approximately 40°C) with gentle dish washing soap (non-abrasive) or in the dishwasher set to a hot cycle

(method B): immerse each part (**6, 8, 9, 10, 10a, 11, 11a, 14**) in a solution of 50% water and 50% white vinegar. Lastly, rinse abundantly with hot potable water (approximately 40°C).

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

Reassemble the canister and the connecting tubing following the directions provided in the "CONNECTION DIAGRAM".

At the end of each use store the medical device complete with accessories in a dry place away from dust.

- The manual vacuum flow control (15) is a sterile disposable product and must be replaced after each use.
- ATTENTION: The hydrophobic/antiviral and antibacterial filter (7) supplied with Flaem aspirators must strictly be replaced for each new patient, or if said filter is saturated.

Should the aspirator be used on the same patient, the filter approved for Single-patient use must be replaced every TWO MONTHS at most or if the filter is saturated.

The filter cannot be sanitised, disinfected, or sterilised.

It is in any case good practice to sanitise and disinfect the device and the accessories, with the exception of the Filter, always after each use, whether it is for the same patient or for a new patient.

STERILISABLE CANISTER	SINGLE PATIENT CANISTER
ACCESSORY SPECIFICATIONS <ul style="list-style-type: none"> - hydrophobic antibacterial/antiviral filter (single-use) - PC collection canister w/ compl. lid - Ø 13 x 7.5 mm L 370 mm silicone tubing - Ø 13 x 7.5 mm L 1300 mm silicone tubing - Ø 13 x 7.5 mm L 80 mm silicone tubing - manual vacuum flow control (single use) 	ACCESSORY SPECIFICATIONS <ul style="list-style-type: none"> - hydrophobic antibacterial/antiviral filter (single use) - PC collection canister w/ compl. lid - Ø 13 x 7.5 mm L 250 mm silicone tubing - Ø 10.5 x 9 mm PVC tubing w/fitt L 2000 mm - Ø 13 x 7.5 mm L 40 mm silicone tubing - manual vacuum flow control (single use)

MICROBIAL CONTAMINATION

In the presence of pathologies with the risk of infection and microbial contamination, we recommend personal use of the accessories, collection canister and connecting tubing (consult your doctor).


Interactions:

The materials used for contact with the secretions are highly stable and chemical resistant thermoplastic polymers (PP, PC, SI). We cannot, however, exclude interactions. Therefore, it is suggested: a) to always avoid prolonged contact of liquid with the canister or tubing and sanitise immediately after use. b) Should anomalous situations occur, i.e. softening or cracking of the accessories, quickly terminate the procedure and substitute the used materials. Contact the authorised service centre and specify how the product was used.

The materials used in the device are biocompatible in accordance with the provisions of Directive 93/42 EC and subsequent amendments. However, the possibility of occurrence of allergic reactions cannot be entirely excluded.

CANISTER HOLDER

For convenience, the canister holder can be detached and reattached on the most suitable side (left or right) as shown in the figure:

Using a tool, unfasten the canister holder from the bottom of the device and rotate it in order to release it. Refasten it, using the inverse procedure, on the selected side. Ensure that the canister holder is properly fastened to the bottom of the device. This system allows up to two canister holders to be connected.

TROUBLE-SHOOTING

Switch off the device before any procedure and unplug the power cable from the socket.

PROBLEM	CAUSE	SOLUTION
The device does not work	The power cable has not been correctly inserted into the socket of the device or the power socket	- Correctly insert the power cable into the sockets
Lack of suction	The collection cannister lid has not been correctly positioned onto the cannister	- Correctly position the collection cannister lid
	Lid gasket not in place	- Correctly position the gasket on the lid
Lack of suction caused by fluid leakage	Blocked filter	- Replace the filter
Blocked float	Incrustation on the float	- Remove the cannister lid and safety device and extract the float. Continue by carrying out the cleaning procedures as described in the paragraph, "CLEANING, SANITISATION, DISINFECTION, STERILISATION"
Vacuum power poor and/or non-existent	Vacuum adjustment knob completely open	- Fully close the adjustment knob and check the vacuum power
	Blocked protection filter	- Replace the filter
	Connecting tubing to the filter and device clogged, bent or disconnected	- Check the condition of the tubing, replace it if blocked and correctly connect it as per the "ASSEMBLY DIAGRAM" of this manual
	Cannister lid overflow valve closed or blocked	- Unblock the overflow valve, keeping the device upright
	Dirty, blocked or damaged pump	- Take the device to your trusted dealer or an authorised FLAEM service centre

If, after checking the above-described conditions the device is still not working properly, we recommend taking it to your nearest trusted dealer or authorised FLAEM service centre.

SPARE PARTS

Description

DURABLE CANISTER (STERILISABLE) complete with tubing and filter
SINGLE PATIENT USE CANISTER complete with tubing and filter

Code

ACO728P
ACO727P

DEVICE DISPOSAL

In compliance with the Directive 2012/19/CE, the symbol printed on the device shows that the device to be disposed of is considered waste and must therefore be an item of "differentiated collection". Consequently, the user must take it (or have it taken) to the differentiated collection sites provided by the local authorities, or turn it in to the dealer when purchasing an equivalent new device. Differentiated waste collection and the subsequent treatment, recycling and disposal procedures promote the production of devices made with recycled materials and limit the negative effects on the environment and on health caused by potential improper waste management. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

ELECTROMAGNETIC COMPATIBILITY

This device was designed to satisfy the current requirements for electromagnetic compatibility (EN 60 601-1-2). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices, in particular with other analysis and treatment devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) might interfere with the operation of electro-medical devices. For further information visit our website www.flaemnuova.it.

The Medical Device may be subject to electromagnetic interference if other devices are used for specific diagnosis or treatments.

Flaem reserves the right to make technical and functional modifications to the product with no prior warning.

SYMBOLS

	Class II device		Device serial number		"ON"	When the device is switched off, the on/off switch stops compressor function in only one phase of the two-phase power supply.
	Type BF applied part		TÜV Certification		"OFF"	
	Attention: check the instructions for use		Minimum and maximum room temperature		EC Marking medical ref. Dir. 93/42 EEC and subsequent updates	
	Single use		Minimum and maximum atmospheric pressure		Production year	
	Alternating current		Minimum and maximum air moisture		Manufacturer	
	Ethylene oxide sterilisation		Without latex		Risk: electrocution. Consequence: Death. Do not use the device while taking a bath or a shower	
	Expiry date		Protection rating of the casing: IP21. (Protected against solid foreign objects larger than 12mm. Protected against access with a finger; Protected against vertically dripping water).			

TECHNICAL SPECIFICATIONS

Non-lubricated rotary piston suction pump with thermal protection.

Mod.	P1211EM/11 I	P1211EM/13.5 I	P1211EM/20 I
Supply/power:	230V~ 50Hz 140VA	230V~ 50Hz 210VA	230V~ 50Hz 210VA
Approvals:			
Fuse:	1 x T 2A 250V	1 x T 2A 250V	1 x T 2A 250V
Device size:	32.5 (W) x 12 (D) x 23.5 (H) cm	32.5 (W) x 12 (D) x 23.5 (H) cm	32.5 (W) x 12 (D) x 23.5 (H) cm
Weight:	2.3 Kg	2.4 Kg	2.4 Kg
Noise (at 1 m):	56 dB (A) (approx.)	56 dB (A) (approx.)	56 dB (A) (approx.)
Use:	Continuous	Continuous	Continuous
Compliant with Directive:	93/42 EEC  0051	93/42 EEC  0051	93/42 EEC  0051
Performance:	high vacuum / low flow	high vacuum / low flow	high vacuum / high flow
Adjustable vacuum level:	from -0.05 to -0.80 approx. (accuracy class 2.5)	from -0.05 to -0.86 approx. (accuracy class 2.5)	from -0.05 to -0.82 approx. (accuracy class 2.5)
Max air flow:	approx 11 l/min	approx 13.5 l/min	approx 20 l/min
Applied parts: Type BF applied parts are:	Possible sterile tube that will be applied to the fitting of the "MANUAL VACUUM FLOW CONTROL".		

Operating conditions: Temperature min. 5°C; max. 40°C
Air humidity min. 10%; max. 95%

Storage conditions: Temperature: min -5°C; max 45°C
Air humidity min. 10%; max. 95%

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© Copyright by FLAEM NUOVA S.p.A., 2019 cod.15481E0-0 rev.03/2019 (TÜV A1)
FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.