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ENGLISH

CONTENTS

0	General Information -----	31
1	Safety Precaution -----	33
2	Caution before use -----	36
3	Product Over View -----	37
4	Operating Instructions -----	39
5	How to Storage and Manage -----	42
6	Trouble Shooting -----	43
7	Warranty -----	44
8	Statements and tables for EMC -----	45
9	Warranty Card -----	54

※ This device is a limb-pressure circulating device, which is a that inflates the sleeve periodically to prevent the blood from accumulating or to circulate blood.

Distributor

Vima s.r.l.
Via Pigafetta 65, 10129, Torino, Italy
P.Iva. 11824140013
R.E.A. Torino 1243754
Registro AEE: iscrizione numero IT18040000010309
Mail: info@mesis.eu - www.mesis.eu

Manufacturer

WelbuTech Co., Ltd.
14, Beakbeom-ro 677 beon-gil, Incheon, republic of Korea

0 General Information

1. Specification

Product Specifications		
Product Name	Sequential Compression Pump	
Model Name	Zam-03	
Dimensions	235mm(L)x245mm(W)x155mm(H)	
Weight (Main Body)	1.8 KG	
Rated Voltage	AC100V~240V,50/60Hz	
Power Consumption	80VA	
Type of protection against electric shock	Class II	
Degree of protection against electric shock	BF-Type Applied Part () to Cuff	
Time (Selectable)	10/20/30 Min±10%	
Maximum Pressure	Less than 297mmHg	
Environment	Operating	<ul style="list-style-type: none"> - a temperature range of +5°C to +40°C; - a relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapour partial pressure greater than 50hPa; - an atmospheric pressure range of 700hPa to 1060hPa.
	Storage and Transport	<ul style="list-style-type: none"> - -25 °C to + 5 °C, and - +5 °C to + 35°C at a relative humidity up to 90%, non- condensing; - >35°C to 70°C at a water vapor pressure up to 50hPa.

2. Intended patient population

- f) Age: 18 years old or greater
- g) Weight: 40 kg or greater
- h) Health: limitation to the people mentioned in contraindication of user manual
- i) Nationality: multiple
- j) PATIENT state:
 - PATIENT is USER: alert, mentally competent (The patient is an intended operator.)
 - PATIENT is not USER: not relevant, unless PATIENT is agitated

3. Part of the body or type of tissue applied to or interacted with

- c) Treatment site: Leg, Arm, and Waist
- d) Condition: No direct contact with the patient, Patient should wear clothes before wearing cuff.

4. Intended user

- f) Education:
 - at least 18 years old
 - have a license of doctor or nurse (professional)
 - no maximum
- g) Knowledge:
 - minimum:
 - read and understand how to operate the device
 - can distinguish: leg, arm, and waist
 - understands hygiene
 - no maximum
- h) Language understanding:
 - understand user manual which is described in English
- i) Experience:
 - No minimum
 - have experience to operate the air compressor for intended use (ex, blood circulation)
- j) Permissible impairments:
 - mild reading vision impairment or vision corrected to log MAR 0.2 (6/10 or 20/32)
 - impaired by 40% resulting in 60% of normal hearing at 500Hz to 2kHz

5. Frequently used functions

- f) connecting/removing the cuffs and air hose to/from main body
- g) wearing/taking off cuffs
- h) switching on/off
- i) setting mode, working time
- j) starting/stopping operation

6. Intended use

This device is intended for use by medical professionals at hospital in treating many conditions, such as:

- ① Primary lymphedema
- ② Edema following trauma and sport injuries
- ③ Post-immobilization edema
- ④ Venous insufficiencies
- ⑤ Lymphedema

The device that increase blood circulation, sequentially pumping air in the cuffs wearred on the limbs, prevent swelling.

I Safety Precaution

This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.

Please comply with "Safety Warnings" since it is aimed to prevent any safety accidents or dangers by the proper and safe use of the product.

Prohibition = Warning of risk of injury or health hazards

Warning = Safety information about possible damage

Caution = Important information

Seller and manufacturer are not responsible for any accidents caused by failure to observe warnings and precautions.

1. Precaution before use

Warning: Please consult the doctor if you have any of the following;

- Person being treated by a doctor,
- Especially the person who feels abnormal in the body
- Person with a malignant tumor
- Person with heart failure
- Person in need of stability
- Person with spinal problems due to accidents or illnesses in the past
- Person with frequent swelling or fever
- The cuff must be cleaned in the specified method and can be reused by other patients after cleaning
- People with diseases such as hip
- Person with a broken calf
- Person with skin-related diseases
- Person with frequent muscle spasms in their calves and thighs
- Person with mental illness
- Person with kidney diseases
- Person with osteoporosis and varicose veins

Prohibition

- Do not use any other method than this manual
- Do not use other devices or electric blankets with this unit
- Do not place heavy objects on the power cord or lay the power cord under the unit. (The power cord may be damaged, causing fire or electric shock.)
- Do not use other than 230V AC. (It may cause electric shock.)
- Do not drop the product from a high place
- Do not use in a bathroom or in a place with high humidity
- Never immerse the appliance in water or other liquids.
- Do not use when the power cord or plug is damaged or the outlet is loose. (It may cause electric shock, shock, or fire.)
- If the product is used at room temperature or below, the operation may not be smooth
- Please limit use within 3 times a day

2. Precaution for use

Warning

- Please keep children and pets away from this unit while using this unit.
- Do not use for anything other than human body.
- Do not use near explosive liquids or flammable liquids or gases.
- Do not use it in an environment with too much electromagnetic wave. It may cause damage to the device and/or cause malfunctioning.
- Do not use the product where it is difficult to disconnect the power plug.
- The cuff must be cleaned in the specified method and can be reused by other patients after cleaning.
- Be careful when using it because it is reused by other patients, and dress should be worn when using cuffs.

Prohibition

- Do not modify, disassemble or repair the product yourself.
- It may cause fire or electric shock. (Do not disassemble or assemble the product.)
- Do not use excessive force on the legs while using the product. It may cause personal injury or product malfunction.
- Do not apply excessive impact while using the product.
- If water has gotten on the unit accidentally, or if you get hurt, turn off the power switch immediately and disconnect the power plug from the outlet.

- If the cuff is worn out, stop using immediately and contact A/S center.
- If the unit is damaged and the inside is exposed, discontinue use and unplug the power cord.
- Do not pull on the power cord strongly.
- Pulling the cord may cause damage to the cord, fire or malfunction.
- Do not unplug the power plug with wet hands. It may cause electric shock or malfunction.
- When unplugging the power plug, be sure to hold the power plug handle and pull it out.
- Please turn off the power switch when you stop using it.
- Do not open the velcro or cuff during operation.

3. Other Precautions

Warning

- When transferring or leasing to another person, attach the instruction manual together with this unit.
- This product is not intended for use in public by many people.
- It is the responsibility of the user to use this product in a public place. In case of product malfunction, A/S will be paid for, and it cannot be exchanged or returned.
- Electromagnetic interference or other interferences may occur between ME devices and other devices, and it is recommended that these devices be used as far away from electronic devices as possible.

Prohibition

- If there is a risk of abrasion, bruising, pain, etc. due to external stimuli or excessive use of cloth, please use after replacing it at the Distributor.
- If you have weak skin or sensitive skin, please put a towel on the site.
- Elderly people, children, and pregnant women may be subject to excessive force when using.
- Be careful not to touch sharp or sharp objects during storage after use.
- Remove dirt and dust, and store in a place with low humidity.
- Do not use when the power cord or power plug is damaged or the plug of the outlet is loose. In this case, please use it after replacing it at the Distributor.
- Please wear thin clothes before using the cuff. There may be damage to the skin.
- Please dispose the product according to the waste procedure. It may cause the environmental pollution.

2 Caution before use

1. Open Box

- Carefully remove the product from its packaging and check the following:
- Check the product carefully for damage.
- Check for missing accessories.
- If the product is damaged or missing, contact your place of purchase or your dealer.

2. Precaution before use

- Be sure to read the instruction manual before using it, and never use it for any purpose other than its intended use.
- Please consult with your doctor if you are already receiving medical treatment.
- The operating pressure can be adjusted to 30 ~ 200mmHg. The cuffs are consumable subject to wear.

Caution

- When using this product for the first time, it is recommended to use low pressure intensities and adjust the operating time.
- Before use, remove any accessories such as rings, bracelets, watches, etc. (Include accessories such as buttons on clothing)
- If you have weak skin or sensitive skin, it may damage your skin. Do not use it directly on the skin.
- This product is equipped with a pressure sensor that allows you to benefit from a constant pressure massage regardless of body measurements.


3 Product Over View

Main Unit and Accessories


Optional items

The below items are not included Cuffs are expendable and can be repurchased when damaged.


Abdominal Cuff Set	Arm Cuff Set	Extender
Abdominal Cuff + Air Hose	Arm Cuff + Air Hose	Extender
	 <p>37</p>	 <p>(For Legs)</p>

1. Control Panel

2. Air Hose Plug Connection


- Connect the air hose plug to the plug-in area before turning on the unit. System alert sound will be played when it is not connected.
- To turn off the sound, connect the plug or turn off the unit.

4 Operating Instructions

1 Putting on Air Massage Cuffs


- Choose a cuff and put it on your body, (legs, arms, or abdomen)
- Connect the air hose plug to the plug-in area.
- ❖ Each massage mode works properly only when the hose plug is connected to the connector in the matching color.

Leg Cuffs

Included Items

- Your legs should be unfolded during the operation, and the acupressure soles can be used together for your foot massaging.
- Zip both cuffs all the way and tight them up with the velcro straps.
- Cuffs can be expanded by applying the extenders (not included).
- Using the rubber band, the acupressure soles can be applied.

(Inner)
Rubber band for holding the soles

Arm Cuff

Not included

Abdominal Cuff

Not included

- Using the band attached to the cuff, fasten the arm cuff enough.
- Cuffs will fall off if they are not fastened up tightly.
- Having the hose located on your left, wear the cuff around your right.
- The cuff should cover your lower back as well as your upper hips.

Wait!

- Remove all jewelry such as necklaces, anklets, rings, and watches before using. (Any buttoned clothing should be avoided also)
- If you have sensitive skin, make sure that your bare skin is not in direct contact with the cuffs.
- Wearing a thin piece of cloth or applying a towel is recommended.


- Arm, abdominal and extenders are not included in the package. There are optional and can be purchased separately. (Contact the seller)

2 Mode Setup

Select Function

- Select your desired mode.
- You can change mode by “Auto” button.
- Default mode is “Auto”.


	ZAM-03
	3 Types:
Mode	Sequence
	Circulatory
	Auto

❖ Sequence of Order for Compression

Leg Cuff
Included

Foot → Calf → Knee → Thigh

3 Controlling Intensity

- Turn clockwise to increase and turn counter clockwise to decrease the intensity.
- Begin at its lowest level and gradually increase the intensity.

Caution

- Make sure that all air hoses don't get folded or twisted during the entire operation and check the air hose plug if you feel any low compression.
- The intensity level should be checked before use.
- If the intensity control knob is pointing **Stop**, the unit will not begin.
 - Turn clockwise to increase to set up your desired intensity.
- To stop the operation, turn the control knob and have it pointed at **Stop**.
 - The unit will stop, and cuffs will be deflated. (Time light blinks.)
- You may feel pain if the intensity level is too high or the duration time is set too long.
- Right level of intensity and time for you should be set.

4 After the use

- Deflation of the cuff will begin when the timed mode is over.
- Remove the air hose plug from the unit.
- Take off the cuff and unplug the power.

5 How to Storage and Manage

1. How to Store

- When not using for a long time, unplug the power plug and store at room temperature.
- Remove dirt and dust, and store in a place with low humidity.
- When not in use for a long time, please cover the cover to keep dust away.
- Detach the accessories from the device.
- Do not store in direct sunlight for a long time every day, or near a hot air balloon (stove) in a hot place.

2. Caring

It is recommended to do after every use

- Before cleaning, be sure to turn off the power switch and disconnect the power plug from the outlet.
- Please wipe the plastic part with neutral detergent by wiping it with a woven cloth and wipe it off with a dry cloth to leave no detergent.
- Do not wash other accessories with a washing machine.
- Do not use water to clean parts made of electronic components.

3. How to clean after use

3.3. How to clean the device

- Disconnect the power cord of the device.
- Clean the outside of the device with a dry cloth. Do not use cleaning solutions

3.4. How to clean the cuff

- Remove the cuff from the appliance.
- Clean the outside of the cuff with a dry cloth.
- Use a cloth soaked in water and connect the cuff to the device after the water has evaporated completely.

4. Precautions

- Do not wipe the main unit or accessories with a damp cloth or wash it with water
- Do not use this product as it may cause discoloration, scratches and cracks. If it is unavoidable, wipe it with a cloth moistened with a little alcohol, then wipe it clean with a dry cloth.
- If you use other cleaning agents or disinfectants, please follow the supplier's instructions.
- If the power cord is damaged or replaced, contact the customer service center for repair.

6 Trouble Shooting

1. When product does not work

- Unplug the power plug once and reconnect it.
 - Make sure that the hose plug is completely connected to the main unit.
 - If the hose plug is not inserted correctly, it is disconnected or it is not connected, it will not work.
 - Insert the hose plug all the way into the body.
- * If the operation does not work, please contact the service center e-mail: assistenza@mesis.eu

2. When the product works but is not smooth

- Press the POWER button to turn on the device.
- Some functions may not work properly if the device has been kept at an unsuitable temperature for a long time.
 - Please use after 2 hours at room temperature
- Check if the air hose is bent.
- Check for air leaks.
 - 1) If the hose is missing from the cuff Insert the missing hose into the hose inlet of the cuff.
 - 2) If the device and hose plug are connected incorrectly -> Insert the hose plug firmly into the device.
 - 3) If air leaks from the cuff or air hose, etc. – Please contact the customer service center for repair and use.
- It may cause static electricity during use, but it is not a malfunction of the product.
- There may be a wind noise while using the product, but this is to prevent the safety of the user and damage to the product. It is a safety device and it is not a malfunction.

3. If the plug or power cord is damaged

- Please contact Customer Service Center for repair.

7 Warranty

1. Warranty Card (Attached to this manual)

- Please read through instructions and warranty details described here
- The warranty card must be presented for any necessary services.
- The manufacturer offers two-years standard warranty from the date of purchase.
- All the cuffs are expendable. (No exchange or replace when damaged)

2. When submitting A Warranty Claim

- Please read through this manual before making your warranty claim
- Following information must be provided for your warranty.
1) Your name, address, and phone number, 2) Model name and serial number, 3) Date of your purchase, 4) Product condition (as much detail as possible).
- The warranty card must be presented for any services within your warranty period. The manufacturer warrants this product against defects under normal use, and you may pay for the service in other cases. See this manual for details.
- If your warranty is expired, please contact the seller.
Before carrying out the repair service, an investigation will be made on the status of the product by telephone or e-mail, to verify the usefulness of any repairs.

3. Other Information

- Please contact the seller or our customer care center for more information.
- The unit comes with two-years standard warranty. (Not for cuffs)
- (Fill out the warranty card upon opening the package for your warranty)
- Date of purchase is valid only with your legible receipt or the seller's sales record.
- This product has been manufactured and distributed after a thorough quality control and inspection. The manufacturer warrant the quality of this product.
- All repair, maintenance, and regular technical inspections must be performed by only service representatives authorized by the manufacturer*.
- * The authorized service center is Vima s.r.l.

8 Statements and tables for EMC

Table 1 – ELECTROMAGNETIC EMISSIONS – for ZAM-01, 02, 03

Guidance and manufacturer's declaration – electromagnetic emissions		
<p>ZAM-01,02,03 is intended for use in the electromagnetic environment specified below. The customer or the user of the ZAM-01,02,03 should assure that it is used in such an environment</p>		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR II	Group I	The ZAM-01,02,03 uses RF energy only for its internal Function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR II	Class A	<p>The ZAM-01,02,03 is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded.</p> <p>Warning : This equipment/system is intended for use by professionals only. It may be necessary to take mitigation measures, such as re-orienting or relocating the ZAM-01, 02,03 or shielding the location.</p>
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

Table 2 – electromagnetic IMMUNITY – for ZAM-01,02,03

Guidance and manufacturer's declaration – electromagnetic immunity			
The ZAM-01,02,03 is intended for use in the electromagnetic environment specified below. The customer or the user of the ZAM-01,02,03 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD) EN 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should ne wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EN 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environmental
Surge EN 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environmental
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	< 5% UT (> 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) (> 95% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (> 95% dip in UT) for 5 s	< 5% UT (> 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) (> 95% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (> 95% dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ZAM-01,02,03 requires continued operation during power mains interruptions, it is recommended that the ZAM-01,02,03 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field EN 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3 - electromagnetic IMMUNITY – for ZAM-01,02,03 that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity			
The ZAM-01,02,03 is intended for use in the electromagnetic environment specified below. The customer or the user of the ZAM-01,02,03 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF EN 61000-4-6	3 Vrms 150kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the [ME EQUIPMENT or ME SYSTEM], including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter Recommended separation distance $d = 1.17 \sqrt{P}$ $d = 1.17 \sqrt{P}$ 80MHz to 800MHz $d = 2.33 \sqrt{P}$ 80MHz to 800MHz where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in meters(M)
Radiated RF EN 61000-4-3	3 Vrms 150kHz to 80 MHz	3 Vrms	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range b. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80MHz and 800MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption reflection from structures, objects and people.			
<p>a Field strengths from fixed transmitters, such as base stations for radio(cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ZAM-01,02,03 is used exceeds the applicable RF compliance level above, the ZAM-01,02,03 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ZAM-01,02,03.</p> <p>b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.</p>			

Table 4 – Recommended separation distances between portable and mobile RF communications

equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEM that are not LIFE-SUPPORTING

Recommended separation distance between portable and mobile RF communications equipment and the ZAM-01,02,03			
The ZAM-01,02,03 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ZAM-01,02,03 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ZAM-01,02,03 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency transmitter m		
	150kHz to 80MHz \sqrt{P} $d=1.17$	80MHz to 800MHz \sqrt{P} $d=1.17$	800Hz to 2.5GHz \sqrt{P} $d=1.17$
0.01	0.117	0.117	0.233
0.1	0.370	0.370	0.736
1	1.17	1.17	2.33
10	3.70	3.70	7.36
100	11.7	11.7	23.3
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

Table 5 – Enclosure port

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS
		Professional healthcare facility environment
ELECTROSTATIC DISCHARG	IEC 61000-4-2	±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM fields a)	IEC 61000-4-3	3V/m f) 80MHz – 2.7GHz b) 80% AM at 1kHz c)
RATED power frequency magnetic field d) e)	IEC 61000-4-8	30A/m g) 50Hz or 60Hz

a) The interface between the PATIENT physiological signal simulation, if used, and the ME EQUIPMENT or ME SYSTEM shall be located within 0.1m of the vertical plane of the uniform field area in one orientation of the ME EQUIPMENT or ME SYSTEM.

b) ME EQUIPMENT and ME SYSTEM that intentionally receive RF electromagnetic energy for the purpose of their operation shall be tested at frequency of reception. Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test assesses the BASIC SAFETY and ESSENTIAL PERFORMANCE of an international receiver might not achieve normal reception during the test.

c) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.

d) Applies only to ME EQUIPMENT and ME SYSTEMS with magnetically sensitive components or circuitry.

e) During the test, the ME EQUIPMENT or ME SYSTEM may be powered at any NOMAL input voltage, but with the same frequency as the test signal (see Table 1).

f) Before modulation is applied.

g) This test level assumes a minimum distance between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic field of at least 15cm. If the RISK ANALYSIS shows that the ME EQUIPMENT or ME SYSTEM will be used closer than 15cm to sources of power frequency magnetic field, the IMMUNITY TEST LEVEL shall be adjusted as appropriate minimum expected distance.

Table 6 – Input a.c power PORT

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS
		Professional healthcare facility environment
Electrical fast transients/bursts a	IEC 61000-4-4	Not applicable ±2kV 100kHz repetition frequency
Surges a) b) j) o) Line – to - line	IEC 61000-4-5	Not applicable ±0,5kV, ±1kV
Surges a) b) j) k) o) Line - to - ground	IEC 61000-4-5	Not applicable ±0,5kV, ±1kV, ±2kV
Conducted disturbances induced by RF fields c) d) o)	IEC 61000-4-6	Not applicable 3V m) 0,15MHz – 80MHz 6V m) in ISM bands between 0,15MHz and 80MHz n) 80% AM at 1kHz e)
Voltage dips f) p) r)	IEC 61000-4-11	Not applicable 0% UT; 1 cycle and 70% UT; 25/30cycles h) Single phase : at 0°
		Not applicable 0% UT; 1 cycle and 70% UT; 25/30cycles h) Single phase : at 0°
Voltage interruptions f) i) o) r)	IEC 61000-4-11	Not applicable 0% UT; 250/300 cycle h)
<p>a) The test may be performed at any one power input voltage within the ME EQUIPMENT or ME SYSTEM RATED voltage range. If the ME EQUIPMENT or ME SYSTEM is tested at one power input voltage, it is not necessary to re-test at additional voltages.</p> <p>b) All ME EQUIPMENT and ME SYSTEM cables are attached during the test.</p> <p>c) Calibration for current injection clamps shall be performed in a 150Ω system.</p> <p>d) If the frequency stepping skips over an ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.</p>		

- f) ME EQUIPMENT and ME SYSTEMS with a d,c power input intended for use a,c – – d.c converters shall be tested using a converter.
- g) Applicable only to ME EQUIPMENT and ME SYSTEMS connected to single-phase a.c mains.
- h) E.g 10/20 means 10 periods at 50Hz or 12 periods at 60Hz.
- i) ME EQUIPMENT and ME SYSTEM with RATED input current greater than 16A/phase shall be interrupted once for 250/300 cycles at any angle and at all phases at the same time (if applicable). ME EQUIPMENT and ME SYSTEM with battery backup shall resume line power operation after the test. For ME EQUIPMENT and ME SYSTEMS with RATED input current not exceeding 16 A, all phase shall be interrupted simultaneously.
- j) ME EQUIPMENT and ME SYSTEMS that do not have a surge protection device in the primary power circuit may be tested only at ± 2 kV line(s) to earth and 1kV line(s) to line(s)
- k) Not applicable to CLASS II ME EQUIPMENT and ME SYSTEMS.
- l) Direct coupling shall be used.
- m) r.m.s., before modulation is applied.
- n) The ISM (industrial, scientific and medical) bands between 0,15MHz to 80MHz are 6,765MHz to 6,795MHz; 13,553MHz to 13,567MHz; 26,957MHz to 27,283MHz; and 40,66MHz to 40,70MHz. The amateur radio bands between 0,15MHz and 80MHz are 1,8MHz to 2,0MHz, 3,5MHz to 4,0MHz, 5,3MHz to 5,4MHz, 7MHz to 7,3MHz, 10,1MHz to 10,15MHz, 14MHz to 14,2MHz 18,07MHz to 18,17MHz, 21,0MHz to 21,4MHz, 24,89MHz to 24,99MHz, 28,0MHz to 29,7MHz and 50,0MHz to 54,0MHz.
- o) Applicable to ME EQUIPMENT and ME SYSTEM with RATED input current less than or equal to 16A/phase and ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16A/phase.
- p) Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16A/phase.
- q) At some phase angles, applying this test to ME EQUIPMENT with transformer mains power input might cause an overcurrent protection device to open. This can occur due to magnetic flux saturation of the transformer shall provide BASIC SAFETY during and after the test.
- r) For ME EQUIPMENT and ME SYSTEMS that have multiple voltage settings or auto ranging voltage capability, the test shall be performed at the minimum and maximum RATED input voltage. ME EQUIPMENT and ME SYSTEMS with a RATED input voltage range of less than 25% of the highest RATED input voltage shall be tested at one RATED input voltage within the range. See table I Note c) for examples calculations.

Table 7 – Input d.c power PORT

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS
		Professional healthcare facility environment
Electrical fast transients/bursts a) g)	IEC 61000-4-4	±2kV 100kHz repetition frequency
Surges a) b) g) Line – to - line	IEC 61000-4-5	±0,5kV, ±1kV
Surges a) b) g) Line – to - ground	IEC 61000-4-5	±0,5kV, ±1kV, ±2kV
Conducted disturbance induced by RF fields a) c) d) i)	IEC 61000-4-6	3V m) 0,15MHz – 80MHz 6V m) in ISM bands between 0.15MHz and 80MHz n) 80% AM at 1kHz e)
Electrical transient conducting along supply lines f)	ISO 7367-2	Not applicable

a) The test is applicable to all d.c power PORTS intended to be connected permanently to cables longer than 3m.
b) All ME EQUIPMENT and ME SYSTEMS cables shall be attached during the test.
c) INTERNALLY POWERED ME EQUIPMENT is exempt from this test if it cannot be used during battery charging, is of less than 0.4m maximum dimension including the maximum length of all cables specified and has no connection to earth, telecommunications systems, any other equipment or a PATIENT.
d) The test may be performed with the ME EQUIPMENT or ME SYSTEM powered at any one of its NOMINAL input voltage.
e) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
f) For ME EQUIPMENT and ME SYSTEMS intended to be installed in passenger cars and light commercial vehicles including ambulances fitted with 24V electrical systems.
g) Direct coupling shall be used.
h) r.m.s., before modulation is applied.
i) If the frequency stepping skips over an ISM or amateur radio band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. this applies to each ISM and amateur radio band within the specified frequency range.

j) The ISM (industrial, scientific and medical) bands between 0,15MHz and 80MHz are 6,765MHz to 6,795MHz; 13,553MHz to 13,567MHz; 26,957MHz to 27,283MHz; and 40,66MHz to 40,70MHz. The amateur radio bands between 0,15MHz and 80MHz are 1,8MHz to 2,0MHz, 3,5MHz to 4,0MHz, 5,3MHz to 5,4 MHz, 7MHz to 7,3MHz, 10,1MHz to 10,15MHz, 14MHz to 14,2MHz, 18,07MHz to 18,17MHz, 21,0MHz to 21,4MHz, 24,89MHz to 24,99MHz, 28,0MHz to 29,7MHz and 50,0MHz to 54,0MHz.

Table 8 – Patient coupling PORT

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS
		Professional healthcare facility environment
ELECTROSTATIC DISCHARGE c)	IEC 61000-4-2	Not applicable ±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Conducted disturbances induced by RF fields a)	IEC 61000-4-6	Not applicable 3V b) 0,15MHz – 80MHz 6V b) in ISM bands between 0,15MHz and 80MHz 805 AM at 1kHz
<p>a) The following apply:</p> <ul style="list-style-type: none"> - All PATIENT-COUPLED cables shall be tested, either individually or bundled - PATIENT-COUPLED cables shall be tested using a current clamp unless a current clamp is not suitable. In cases where a current clamp is not suitable, an EM clamp shall be used. - No international decoupling device shall be used between the injection point and the PATIENT COUPLING POINT IN any case. - Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS. - Tubes that are intentionally filled with conductive liquids and intended to be connected to a PATIENT shall be considered to be PATIENT-COUPLED cables. - If the frequency stepping skips over an ISM or amateur radio band, This applies to each ISM and amateur radio band within the specified frequency range. - The ISM (industrial, scientific and medical) bands between 0,15MHz and 80MHz are 6,765MHz to 6,795MHz; 13,553MHz to 13,567MHz; 26,957MHz to 27,283MHz; and 40,66MHz to 40,70MHz. The amateur radio bands between 0,15MHz and 80MHz are 1,8MHz to 2,0MHz, 3,5MHz to 4,0MHz, 5,3MHz to 5,4MHz, 7MHz to 7,3MHz, 10,1MHz to 10,15MHz, 14MHz to 14,2MHz, 18,07MHz to 18,17MHz, 21,0MHz to 21,4MHz, 24,89MHz to 24,99MHz, 28,0MHz to 29,7MHz and 50,0MHz to 54,0MHz. <p>b) r.m.s., before modulation is applied</p> <p>c) Discharges shall be applied with no connection to an artificial hand and no connection to PATIENT simulation. PATIENT simulation may be connected after the test as needed in order to verify BASIC SAFETY and ESSENTIAL PERFORMANCE.</p>		

9 Warranty Card

· Product Name :

· Type Designation:

· Date of Purchase:

· Customer :

· Type:

· Term of Warranty:

· Retailer:

This product has been manufactured and distributed after a thorough quality control and inspection. We warrant the quality of this product.

This product has been manufactured and distributed after a thorough quality control and inspection. We warrant the quality of this product.

No cost will be incurred and free repair will be available for this product in case of any troubles or breakdowns originating from this product despite normal and accurate usage during the warranty term (two years) calculated from the date of purchase according to the Act on Protection of Customers.

However, if the warranty term has expired or in case of any of the following during the term, the actual costs of repair will be requested.

- 1) Breakdown or damage coming from willful or gross negligence or poor/misdirected handling
- 2) Damages coming from natural disasters and accidents (such as fire, smoke pollution, gas pollution, earthquake, storm, lightning, etc.)
- 3) Breakdown coming from inappropriate repair or remodeling
- 4) No provision of the warranty or false information on the warranty
- 5) Any breakdown coming from customer's negligence or failure to keep the safety rules provided on the manual

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