
GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 - fax 199 400 403
Export: tel. +39 02 953854209/221/225 fax +39 02 95380056
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Sfigmomanometro digitale senza mercurio su carrello
Digital sphygmomanometer without mercury on cart
Tensiomètre numérique sans mercure avec pied à roulettes
Digitales, quecksilberfreies Blutdruckmessgerät mit Ständer
Esfigmomanómetro digital sin mercurio sobre carro
Esfigmomanômetro digital sem mercúrio sobre rodas
Ψηφιακό πιεσόμετρο χωρίς υδράργυρο σε καροτσάκι
جهاز رقمي لقياس الضغط الدموي بدون زئبق على عجلة

MANUALE D'USO E MANUTENZIONE
USE AND MAINTENANCE BOOK
INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN
BETRIEBS UND WARTUNGS ANWEISUNGEN
MANUAL DE USO Y MANTENIMIENTO
MANUAL DE USO E MANUTENÇÃO
ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ
دليل الإستعمال والرعاية

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TECHNICAL FEATURES

The electronic sphygmomanometer without mercury is a control instrument for the measurement of blood pressure. Made out of ABS plastic and with an ergonomic design it assures precise measurements.

The reading of the values is made easier thanks to the double display (graduated and digital scale). An LCD column replaces the mercury column thus avoiding possible leaks of mercury into the environment that could be dangerous for people also.

Display method: high definition LCD / digital display

Measurement range: 0-300mmHg (0-40kPa)

Precision: +/-3mmHg (+/-0.4kPa)

Power supply: two "AA" alkaline batteries

Environmental conditions of use: 10-40° C, 30-85% RH

Environmental conditions for storage: -20 +55°C. 10-95% RH

Dimensions: Body: 360x96x66 mm

Cuff: from 22 to 33 cm

Weight : 1500g (batteries included)

PRESCRIPTIONS

Do not use the equipment in case it is damaged. Apply to your retailer.

Avoid precarious repairs. Repairs shall be carried out with original spare parts only, which shall be installed according to the intended use.

Since the product is made of corrosion-proof materials suitable for the environmental conditions foreseen for its normal use, does not require special care, however it is necessary to store it in a closed place making sure that is protected from dust and dirt to assure its hygienic conditions. Moreover, it is recommended to store the product in a place which can be reached easily by the personnel in case of necessity.

USE

Always follow the indications of your doctor as concerns modalities and frequency of measurements.

Do not position the cuff on parts of the body other than the arm and do not use if not properly fastened.

Before the measurement

Perform measurement in a comfortable environment. Very hot or very cold temperatures could affect the measurement.

If drinks containing caffeine, such as coffee or cola, have been taken wait approx. 30-45 minutes. Do not smoke right before the measurement. Lay down and relax for approx. 10 minutes before starting the measurement.


During the measurement

Sit down and stay still as much as possible during the measurement. Talking or moving may increase the value of the measurement. Do not cross legs and do not touch the cuff during the measurement.


To obtain proper results to compare, always measure the pressure on the same arm and possibly at the same time. We suggest measuring the blood pressure on the left arm.

To perform multiple measurements wait at least 10 minutes between each measurement.

Fasten the cuff around the upper part of the arm. The cuff must be tightly fastened but not too tight. Remove anything that may block circulation (rolled-up sleeves or tight clothes). Make sure the cuff is at the same height of the heart.



Auscultation Mode

1. Set the mode selector to the high position (AUSC) to select the auscultation method of measurement.


2. Press the 'ON/OFF' button. After the loading phase, '000' will appear on the display and the LCD scale will be set to zero. The device has an auto-diagnosis function: if residual air


is detected in the cuff, a flashing 'P' will appear. Once the residual air has been eliminated, '000' will appear on the display.


3. Insert the stethoscope's bell into the cuff, positioning it over the brachial artery in the crook of the elbow.

4. Close the air vent valve and start to inflate the cuff using the pump. When the indicator reaches a value higher than around 2.5~4.0 kPa (18.75~30 mmHg) compared to the normal systolic pressure, inflation can be stopped.


5. At the end of the inflation, the air vent valve will open automatically. Air is slowly released and the pressure decreases. The first pulsation detected by the stethoscope corresponds to the value of the 'systolic or maximum pressure' on the LCD column display. During decompression, the pulsations will continue and then decrease until they suddenly disappear or are so faint that they are imperceptible. The pressure shown when the pulsations disappear corresponds to the 'diastolic or minimum pressure'.


Oscillographic Mode


1. Set the mode selector to the low position (OSCI) to select the oscillographic method of measurement.

2. Press the 'ON/OFF' button. After the loading phase, '000' will appear on the display and the LCD scale will be set to zero. The device has an auto-diagnosis function: if residual air is detected in the cuff, a flashing 'P' will appear. Once the residual air has been eliminated, '000' will appear on the display.



2. Press the 'ON/OFF' button. After the loading phase, '000' will appear on the display and the LCD scale will be set to zero. The device has an auto-diagnosis function: if residual air is detected in the cuff, a flashing 'P' will appear. Once the residual air has been eliminated, '000' will appear on the display.

3. Check that the symbol on the cuff is positioned in correspondence with the brachial artery, in the crook of the elbow.

4. Close the air vent valve and start to inflate the cuff using the pump. When the indicator reaches a value higher than around 2.5~4.0 kPa (18.75~30 mmHg) compared to the normal systolic pressure, inflation can be stopped.


5. At the end of the inflation, the air vent valve will open automatically. Air is slowly released and the pressure decreases. The device will begin to measure the blood pressure and heart rate automatically.

SYS: 112 mmHg
DIA: 76 mmHg
PUL: 72 n/min

6. On the LCD column display 2 bars will appear. The bar on top is the 'systolic or maximum pressure' value, the bar at the bottom shows the 'diastolic or minimum pressure'. The heart rate is shown on the number display, underneath the LCD column display.

* When the LCD display shows the symbol in the illustration it means the batteries are low and the sphygmomanometer cannot operate. Replace the "AA" alkaline batteries with two new ones.

After the measurement, completely open the valve to let the residual air out, remove the cuff, press the ON/OFF button to switch device off. (This device automatically switches off after five minutes).

Clean the cuff and place it in its case along with the tube which, to avoid damages, must not be tightly coiled up.

HOW TO STORE AND MAINTAIN THE DEVICE

Keep the product far from direct sunlight, humidity and rapid temperature changes.

Avoid crashes or falls.

Do not keep body detached from cuff and do not fix it in a place different from its stand. If you want to clean the device, use a soft dry cloth or a slightly damp cloth. Do not use alcohol, benzene solvents or other aggressive chemical products to clean the device or cuff.

Before using the device, wash hands. Do not wash or wet the cuff or the pump.

Note: When the device is not used for a long period of time, remove the batteries because they could damage the device.

What is blood pressure?

Blood is sent to the arteries with the action of a pump called heart (contraction and dilation). The pressure of blood leaving the heart is called "blood pressure", Blood Pressure pulsate with each beat of the heart. The high blood pressure when the heart contracts are called "Systolic pressure" and the low blood pressure when the heart dilates is called "diastolic pressure", The threshold value for hypertension in adults is defined by the World Health Organization (WHO) as 140/90 mmHg.

Health and blood pressure!

When people reach middle age, the risk of hypertension markedly increase. With aging, the blood vessels age rapidly, Furthermore, because of the obesity and lack of exercise, cholesterols stick to blood vessels, causing them to lose elasticity, Therefore, watching daily blood pressure help to evaluate our health condition.

Why do we need to monitor blood pressure at home?

By recording the blood pressure values and the measuring conditions such as the measuring time or living state every day, you can know the fluctuation tendency of your blood pressure, which helps control your health. Furthermore, recording of daily blood pressure values is very helpful for your doctor to diagnose.

How to manage your Blood Pressure if it works unsuitable?

If correct measurement is impossible even after checking the above-mentioned points, consult at the store where you have purchased the unit or the nearest dealer without touching the internal mechanism.

In some very rare cases, there may be error due to the physical condition of the person. In such cases, please consult your doctor.

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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