GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 - fax 199 400 403
Export: tel. +39 02 953854209/221/225 fax +39 02 95380056
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Termometro auricolare a raggi infrarossi
Infra-red ear thermometer
Thermomètre auriculaire à rayons infrarouges
Termómetro infrarrojo de oído

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
The ear thermometer is a device capable of achieving infra-red temperature measurement when placed in the auditory canal of a subject. It is a safe means of measuring human body temperature through the ear.

Specification

The standard we adopted is EN12470-5:2003 Clinical thermometers -- Part 5: Performance of infrared ear thermometers (with maximum device).

This device complies with the requirements of IEC 60601-1, IEC 60601-1-11 and IEC 60601-1-2.

Range of displayed temperature: 34.0° C ~ 44.0° C (93.2° F ~111.2° F)

Display L °C (L °F), when the temperature is under 34,0° C (93,2° F)

Display H °C (H °F) when the temperature is over 44,0° C (111,2° F)

Operating ambient:

Temperature: 16°C ~35°C (60.8°F ~95.0°F)

Humidity: ≤80% UR

Storage and transport condition:

Temperature: da -25°C a 55°C (da -13°F a 131°F)

Humidity: 85% UR

Display resolution temperature range: 0.1°C (0.1°F)

Accuracy: ±0.2°C (da 35.5°C a 42.0°C) / ± 0.4°F (da 95.9°F a 107.6°F)

Display: liquid crystal display, 3,5 digits

temperature value: display the maximum temperature in measuring process

temperature unit: centigrade or fahrenheit

display of memory: last ten memories

low voltage warning: the LCD display  and then a beep sound is heard

Power consumption: 9 milliwatt in measurement mode

Battery: one 3V Lithium battery (CR 2032)

Battery life: 2000 takes

Dimension: 110 mm x 35 mm x 35 mm

Net weight: Approx.37g

Beeper sign: on/off, measuring finish and low voltage warning, etc.

Self-testing sequence: Press the 'on/recall' button to turn on the thermometer and all of symbol (See Fig A) should be displayed on the LCD in one second.

Manufacturing and Calibration date: See the label in the battery box.

Cautions

1. Before measuring, the thermometer shall be stabilized at least 30min under operating condition.
2. Please keep your ear canal clear, if not, the measuring result shall be inaccurate.
3. The symbol ' ? ' on the LCD shows you can measure in the ear;

4. Please don't scratch the probe membrane, if not the thermometer shall lose the efficacy;
5. If the LCD no changes under pressing any button, please take the batteries out and put in again;
6. Please don't use the thermometer if your ear canal has been inflammation.
7. : Type BF equipment.

How to use

Measuring human body temperature in the Ear Canal

1. Press the 'on/recall' button to turn on the thermometer, a beep sound is heard and the LCD displays as Fig A in one second.


Fig. A

Fig. B

When the LCD displays as Fig B, now ready to measure.

2. Straighten the ear canal by pulling the outer ear up and back to give a clear view of the eardrum.

For children under 1 year, Pull the ear straight back.

For children ages 1 year to adult, Pull the ear up and back. (Caution: Many pediatricians suggest use of ear thermometers for children older than 6 months.)

Infants

Adults

3. While pulling the outer ear, insert the probe snugly into the ear canal (best to insert deep), press the 'scan' button, release it. When you hear a beep sound, the measurement is complete. Then remove the thermometer from the ear. The display show the measure temperature.

4. Measuring again: If measuring again, please unclinch the 'scan' button and wait for the '?' symbol to display, then press 'scan' button to measure again. Advice: The time between each reading should be not less than 1 mins.


Fig. C


Fig. D

5. Memory search: Press the 'on/recall' button, the last ten memories (NO.9-NO.0) will take turns to display on the LCD as Fig C and Fig D.


Fig. E

6. The LCD will displays 'Er 1' as Fig E and cannot displays temperature when the ambient temperature is out of the range of 10 Cto 35 C.

Note: Before measuring, the thermometer shall be stabilized at the operating ambient condition for a minimum of 30 min.

7. Shut off: The thermometer shall automatically shut off without manipulating in one minute.

Changing from Celsius to Fahrenheit (only for switchable type)

To switch the display between C and F, turn the unit on. Press and hold the on/recall button for 3 seconds, the display will show as Fig. F, the unit will prepare for a measurement.



Fig. F

Cleaning and Storage

1. Store thermometer in a dry location free from dust and contamination and away from direct sunlight. The ambient temperature at the storage location should remain fairly constant and within the range of -25 C to 55 C.
2. Use an alcohol swab or cotton swab moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer.
3. Never use abrasive cleaning agents, thinners or gasoline for cleaning and never immerse the instrument in water or other cleaning liquids. Take care not to scratch the surface of the probe lens or display.

Replacing the batteries

1. When voltage of the battery is low, the LCD will display “” symbol, please replace a new battery in the unit. The thermometer cannot work accurately under the condition of low voltage.
2. The thermometer is supplied with a 3 V LITHIUM BATTERY (CR 2032), Insert a new battery when the low voltage symbol appears on the LCD.
3. Remove the battery cover and pull down the metal button(see Fig),the battery will be bounced out, then you can take out the old battery.
4. Place a new battery, the positive pole “+ “ up and negative pole “ - “ down.
5. Please take out the battery to avoid battery leaking if unit not used for over six months.
6. The disposal of the battery and the device shall comply with the local environment requirements.

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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