GIMA

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
Conservare a 2-35°C

INTENDED USE

The pH test strip is a fast pH value test of vaginal fluid specimens.

Vaginal pH is a very useful screening parameter as an indication for most clinical and sub-clinical infections in the vagina, a common cause of free-term delivery, low weight birth and reduced fertility.

The pH test has been designed to give the physician a clinical reference and increase chances of a successful pregnancy and normal development of patient's babies. It can be used as periodical check even for woman who isn't pregnant to avoid vaginal infections.

This product is intended for professional use only. The assay should not be used without proper supervision and is not intended for over the counter sale to lay persons.

SUMMARY

pH test strip is effective risk free and totally painless test. It consists of two parts: a swab and the test strip are packed separately to assure sterile conditions for swab and a dry place for the test strip. Thus, the test can be performed with totally hygienically conditions. The test strip is ready-to-use and disposable.

It is dry chemical reagent strip for testing pH value of vaginal fluid specimens, and is intended for in vitro diagnostics use. It is a visual test.

TEST PRINCIPLE

The pH test strip is a fast dry chemical method to test pH value of vaginal fluid specimens. It is made of one reagent pad which is affixed to a strip of plastic.

When the strip is wetted with a sterile swab, the pH value of the specimen will define its color. By comparing the color of reagent pad with the standard color chart on the package foil, the detail result will be obtained.

BASICS

Clinical or not clinical infection in the vagina, can provoke a pre-term delivery or pre-term water breaking. Most of these infections can be previously determined by modifications of vaginal acidity level .

Sub-clinical infections can reduce fertility or provoke sterility in no-pregnant woman. A routine-test like vaginal pH test can help you to previously determine infections like these. Normal vaginal fluid contains a variety of micro organisms. Among these lactobacilli predominate. Lactobacilli produce lactic acid causing a protective acid environment. A healthy vagina has a pH lower than 4.5. The acidity of the vaginal environment will make the growth of infections germs difficult and prevent their migration to the uterus where they could cause early contractions or rupture of the membranes. The growth of the germs in the vagina mostly goes unnoticed. The earlier this disorder is detected by a regular testing of the vaginal pH, the earlier appropriate action can be taken to reduce the risk of premature birth.

REAGENTS AND MATERIALS SUPPLIED


- One individually wrapped test strip
- One individually wrapped swab
- One instruction sheet

STORAGE

- Store at room temperature or between 2°C – 35°C.
- All test strips should be stored in the sealed pouch.
- Do not expose to direct sunlight.
- Do not use after expiration date.

PRECAUTIONS

- For *In Vitro* use only
- For professional use only
- Every test (swab and test strip) is for single use only. Use only once
- Do not touch the swab or the test strip at the reaction pad to avoid contamination
- Keep the swab sterile before use
- Do not spill the samples into the reagent pad
- Do not use urine as liquid
- Vaginal specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Do not use test devices after expiration date listed on pouch
- Do not use the test after damage of the packaging foil
- Use test right after unwrapping
- Please be aware of the developing time of the test (1 min) before evaluation
- It is recommended that all laboratory personnel interpreting the results of these strips be tested for not color blindness.
- Please take the specificity and the cross reactivity into account for evaluation
- Store and transport the test device always at 2-35°C

SPECIMEN COLLECTION AND HANDLING

- Wait 15 minutes after last urine flow to perform the test. Patient should make sure (by wiping) that no urine is left at the vaginal entrance, if she just has used the toilet, as urine may change the acidity and the test result.
- It is important not to perform the test within 12 hours of having a sexual intercourse or 12 hours after introducing the patient's vagina any substances like spray, cream, gel and so on.
- The pH test strip can't be used during the patient's period.
- pH test strip can be performed 3-4 days after the end of monthly period.
- The pH test strip can't be used when the patient is taking antibiotics. These factors may change the vaginal pH and give invalid results.
- Before using the test the patient should not have any bath or shower.
- Patient should not wash, rinse vagina before test because it may change her vagina pH value.

TEST PROCEDURE


Review "Specimen collection" instructions. Do not open package foils until ready to perform the assay.

1. Open a swab pouch and remove a swab.
2. Insert the swab into the vagina no deeper than 3 cm (1.5 inch). It is necessary to leave the swab in the vaginal tract to soak the swab
3. Withdraw the swab.
4. Do not touch the swab with fingers after withdrawing it.
5. Touch the reagent pad with the swab and transfer vaginal fluid to the reagent pad. It should be totally soaked.
6. Read the result in one minute by comparing the color of the reagent pad with the color chart printed on the pouch. Changes in color after 2 minute are of no diagnostic value.

If the patient's pH value is outside the normal pH range on one occasion perform a second test a few hours later or the next morning. If the second measurement still displays an elevated pH value a further investigation by a physician should be done. Please keep in mind, that stress or certain food may temporarily change vaginal pH value.

INTERPRETATION OF RESULTS

Semi-quantitative results are obtained by visually comparing the color of pad with the corresponding test color block displayed on the container label. No equipment is required.

A healthy vagina has a pH lower than 4.5. Thus, patients with a pH higher 4.5 should be further checked by a physician.

QUALITY CONTROL

For best results, performance of reagent strips should be confirmed by testing known negative and positive specimens or controls whenever a new test is performed. Each laboratory should establish its own goals for adequate standard of performance, and should question handling and testing procedures if these standards are not met.

LIMITATION

Comparison to the color chart is dependent on the interpretation of the individual. It is therefore, recommended that all laboratory personnel interpreting the results of these strips be tested for color blindness. Like with all laboratory tests, definitive diagnostic or therapeutic decisions should not be based on any single test result or method.

Performance Characteristics

The test device can be evaluated 1 minute after addition of specimen to the test device.

Accuracy

The accuracy of the pH test was compared and checked against a commercially available test. 30 samples were examined under both tests. The results were 100% in agreement.

Sensitivity

The analytical sensitivity of pH test was established by repetitive testing of buffers solutions containing pH 3.8 / 4.1 / 4.4 / 4.6 / 4.9 / 5.1.

The test was in 100% agreement with a pH meter.

Specifity

Following actions and substances are possible interfering:

- Urine is left at the vaginal entrance, as urine may change the acidity and the test result.
- Performing the test within 12 hours of having a sexual intercourse or 12 hours after introducing the patient's vagina any substances like spray, cream, gel and so on.
- Performing the test during the patient's period.
- Taking antibiotics. This factor may change the bacterial environment of the vaginal tract. Thus, a change in the pH might occur and give invalid results.
- Having any bath or shower before using the test, because soaps and detergents might change the pH of the vaginal fluid.
- Wash, rinse vagina before test because it will change the vagina pH value.

The test is influenced by many detergences, soaks, disinfectious agents and alcohols if these substances are used at the vagina. But in addition, antibiotics and some pharmaceuticals might influence the test indirectly by changing the bacterial environment.

Food or drinks have little influence on the vaginal pH test. Eating food or taking drink will not have a big influence on body fluid, particularly on vaginal pH.

SIMBOLS

For in-vitro
diagnostic use only

For single use only

Content

Expiry date

Lot number

Store at 2-35°C