

ERKA. KOBOLD
Ø48 MM

RELIABLE AND ACCURATE
IN EVERYDAY MEDICAL PRACTICE



ERKA. KOBOLD silver anodised 48 mm with one GREEN CUFF Superb rapid cuff grey

Size 1	227.28493
Size 2	227.27493
Size 3	227.26493
Size 4	227.20493
Size 5	227.24493
Size 6	227.29493

The following variations are also available upon request
via made@erka.org:

- Surface in 5 different colors
- Custom printed cuffs and scales,
e.g. with your company name and company logo

The ERKA. Kobold is a robust and reliable partner for every professional user.

The ergonomically balanced design combined with a weight-optimized design and the extra-sensitive screw valve guarantee maximum user comfort in the physician's daily routine – with each individual measurement.

- 1_ Special tempered copper-beryllium membrane ensures durability with consistently high precision.
- 2_ With a diameter of 48 mm, the robust scale is optimally readable.
- 3_ The anodized aluminum housing is easy to clean and ensures optimal surface hygiene.
- 4_ The selection of ERKA. GREEN CUFF cuffs is comprehensive. They are durable, skin-friendly, biocompatible and represent optimal wearer comfort.

REPLACEMENT PARTS AND CUFFS

BALLS



PUMP BALL black
000.20400

CUFFS

Single tube-Cuffs



ERKA. GREEN CUFF
Superb rapid cuff
grey



ERKA. GREEN CUFF
Superb hook cuff
grey

Size 1 infant	10-15 cm
Size 2 small child	14-21.5 cm
Size 3 adult/small	20.5-28 cm
Size 4 adult	27-35 cm
Size 5 adult/large	34-43 cm
Size 6 thigh	42-54 cm

024.28493
024.27493
024.26493
024.20493
024.24493
024.29493

025.26293
025.20293
025.24293

**ERKA. KOBOLD
SMART Ø48 MM**

100% HYGIENE IN INTERACTION
WITH PATIENTS



With the innovative Smart Concept, ERKA. provides optimal infection prevention along with absolutely precise measured values. The use of an ERKA. Kobold Smart in conjunction with an ERKA GREEN CUFF Smart Single patient soft cuff – specially designed for single use – significantly reduces the risk of cross-contamination.

In addition, the ERKA. hygiene concept is continuous. The reusable ERKA. GREEN CUFF Smart cuffs can be superbly cleaned, as well as quickly, simply and absolutely reliably disinfected.

- 1_ The anodized aluminum housing of the ERKA. Kobold Smart is easy to clean and guarantees maximum surface hygiene.
- 2_ With a diameter of 48 mm, the robust scale is optimally readable.
- 3_ The appropriate cuff size for every patient - the selection of the ERKA. set of cuffs is comprehensive and ranges from Size 1 (infant) to 6 (thigh).
- 4_ The Smart concept is simple and effective: the size is selected on the basis of color-standardized cuff tubing; cuffs can be changed in seconds with the ERKA. quick closure.
- 5_ For maximum hygiene, the device is also available with disposable cuffs.

REPLACEMENT PARTS AND CUFFS

BALLS / CONNECTORS



PUMP BALL black
000.20400



QUICK-CONNECTOR
995.00001



CUFFS

Single tube-Cuffs

Size 1 infant	10-15 cm
Size 2 small child	14-21.5 cm
Size 3 adult/small	20.5-28 cm
Size 4 adult	27-35 cm
Size 5 adult/large	34-43 cm
Size 6 thigh	42-54 cm

ERKA. GREEN CUFF
Smart rapid cuff
grey

024.28493
024.27493
024.26493
024.20493
024.24493
024.29493

ERKA. GREEN CUFF
Smart single patient soft
white, Box of 10

024.58473
024.57473
024.56473
024.50473
024.54473

TUBES

TUBE grey
995.00053

ERKA. KOBOLD SMART 48 mm
silver anodised with GREEN CUFF
Smart rapid cuff grey

Set	
Size 1-6	227.71193
Size 3-5	227.71293
Size 1-3	227.71393

ERKA. KOBOLD SMART 48 mm
silver anodised with one GREEN
CUFF Smart single patient soft white

Box of 10	
Size 1	227.91873
Size 2	227.91773
Size 3	227.91673
Size 4	227.91073
Size 5	227.91473

The following variations are also available
upon request via made@erka.org:

- Surface in 5 different colors
- Custom printed cuffs and scales,
e.g. with your company name
and company logo