

## Column and floor scales



Eye-level beam scale, readable from both sides, equipped with a telescopic height measuring rod as standard and two transport castors for easy movability. This solid and robust scale is made of steel and comes with a heavy-duty cast-iron lever system, guaranteeing a long and trouble-free life time. The large base is covered by a removable non slip cover for easy cleaning and secure use.

### Technical specification:

Capacity:	200 kg
Graduation:	100 g
Dimensions:	280 x 560 x 1470 mm
Weight:	15 kg
Platform size:	275 x 375 mm
Measuring range:	60 – 210 cm by 1 mm