

**ERKA. SWITCH 2.0
COMFORT Ø56 MM**

FUNCTIONAL AND RELIABLE
IN EVERYDAY USE



48

Based on ERKA's specifications for a modern blood pressure monitor and in close co-operation with physicians worldwide, ERKA. – a pioneer in blood pressure monitoring – has developed the ERKA. SWITCH 2.0 Comfort.

Outstanding in quality and precision, convenient in operation and handling, impressive in efficiency and advantages, and responsible in managing environmental resources – the modern standard, defined by ERKA.

- 1_ The ability of the user to change the measuring vessel in seconds, directly and in situ, saves time and resources.
- 2_ The device meets the requirements for use in an MR environment and is shock resistant.
- 3_ The use of a special tempered copper-beryllium membrane ensures durability along with the consistent high precision of the measuring device.
- 4_ The sensitive screw valve enables the user to take perfect measurements.
- 5_ The innovative design allows the device to be switched from a right-handed to a left-handed device in a single step.
- 6_ The selection of ERKA. GREEN CUFF cuffs is comprehensive. They are durable, skin-friendly, biocompatible and represent optimal wearer comfort.

49

REPLACEMENT PARTS AND CUFFS

BALLS/MEASURING VESSELS



PUMP BALL grey
000.20900



MEASURING VESSEL SWITCH 2.0 COMFORT
000.30200



CUFFS
Single tube-Cuffs

ERKA. GREEN CUFF
Superb rapid cuff
grey



ERKA. GREEN CUFF
Superb hook cuff
grey

ERKA. SWITCH 2.0 COMFORT 56 mm with one GREEN CUFF Superb rapid cuff grey

Size 1	293.28493
Size 2	293.27493
Size 3	293.26493
Size 4	293.20493
Size 5	293.24493
Size 6	293.29493

The following variations are also available upon request
via made@erka.org:

_ Custom printed cuffs and scales,
e.g. with your company name and company logo

Size 1 infant	10-15 cm	025.28493	
Size 2 small child	14-21.5 cm	025.27493	
Size 3 adult/small	20.5-28 cm	025.26493	025.26293
Size 4 adult	27-35 cm	025.20493	025.20293
Size 5 adult/large	34-43 cm	025.24493	025.24293
Size 6 thigh	42-54 cm	025.29493	