

ERKA Finesse

The professional standard.

The stethoscope for every medical professional!

With a switchable double bell, this stethoscope covers every eventuality of daily life in a medical practice. The open funnel side is intended especially for low frequencies. The membrane side covers the full frequency range. A chrome-plated, high-gloss polished surface is free of pores, and all of the materials used can be cleaned perfectly and quickly for maximum surface hygiene.



ERKA standard membrane with convex tension for outstanding resonance and sound transmission. g für eine hervorragende Resonanz- und Schallübertragung.



The brass bell, polished or with a matt chrome finish – both versions available in adult and child sizes. A ring for protection against cold on both sides maximises patient comfort.



Two parallel interior channels run inside the preformed tube – available in numerous different colours. This avoids friction noise which occurs when two separate tubes are used.



Thanks to the spring integrated in the tube, acoustic signals are conducted from the bell to the ear with virtually no loss.



The ergonomically shaped, super soft ear buds offer a good fit for enhanced wearer comfort and prevent acoustic loss. Naturally, the soft ear buds shown are also included as an alternative.

	Material · Surface	black	grey	blue	green	burgundy	red
Finesse	Brass · bright chrome	550.00000	550.00005	550.00020	550.00055	550.00060	550.00030
	Brass · satin chrome	550.01000	550.01005	550.01020	550.01055	550.01060	550.01030
Finesse Child	Brass · bright chrome	549.00000	549.00005	549.00020	549.00055	549.00060	549.00030
	Brass · satin chrome	549.01000	549.01005	549.01020	549.01055	549.01060	549.01030
Finesse ²	Brass · bright chrome	535.00000	535.00005	535.00020	535.00055	535.00060	535.00030
	Brass · satin chrome	535.01000	535.01005	535.01020	535.01055	535.01060	535.01030
Finesse ² Child	Brass · bright chrome	536.00000	536.00053	536.00020	536.00055	536.00060	536.00030
	Brass · satin chrome	536.01000	536.01005	536.01020	536.01055	536.01060	536.01030

Special colors as pink  or bright blue  are available by request.

Spare-part overview

	Type	Sensitive	Precise	Classic	Finesse	Finesse ²	ERKAphon duo	ERKAphon
Membrane	Dual-Membrane	000.40200	000.40200	000.40200	-	-	-	-
	Standard Membrane baby	-	-	-	-	000.40102	-	-
	Standard Membrane chid	-	-	-	000.40100	000.40100	000.40100	000.40100
	Standard Membrane adult	-	-	-	000.40000	000.40000	000.40000	000.40000
non-chill ring	Dual-Membrane	000.53104	000.53104	000.53104	-	-	-	-
	Size baby	-	-	-	-	000.53103	-	-
	*optional for ERKAphon Stethoscops	-	-	-	000.53142	000.53142	000.53142 *	000.53142 *
	Size adult	-	-	-	000.53141	000.53141	000.53141 *	000.53141 *
Non-chill ring for funnel		-	000.53304	-	000.53304	-	000.53304	-
Eartips	Soft	-	-	000.45004	-	-	000.45004	000.45004
	*optional for ERKAphon Stethoscops	000.45104	000.45104	-	000.45104	000.45104	-	-
	Super Soft	000.44900	000.44900	000.44900	000.44900	000.44900	000.44900 *	000.44900 *

**Device type
Variations**

Stethoscope with headpiece integrated in dual-duct tube with membrane and funnel side, switchable, with ring as protection against cold

UMDNS code:13-755 Medical device according to class 1

Stethoscope with dual bell in decorative carton with instructions for use

Variations. Sizes (bell): Ø 47mm or Ø 33mm, 6 different colours (tube), electroplating: gloss or matt finish

Dimensions

Overall length adult/child: 730 mm/716 mm

Tube length: 555 mm

Earpiece length/Ø: 150 mm/6 mm

Bell dimensions

Adult: 47 mm/44 mm/30 mm/22 mm

Child: 33 mm/30 mm/30 mm/22 mm

(Ø bell/Ø membrane/Ø funnel/height)

Material

Bell: Brass, polished and chrome-plated Earpiece: Brass, chrome-plated Yoke spring: Spring steel, chrome-plated

Ear buds: TPE/TPS, colour: grey Ring as protection against cold/funnel protection ring: NBR

Membrane: boPET (Hostaphan RN 250), 0.25 mm thick Moulded tube: Soft PVC, latex-free, resistant against common cleaners and disinfectants

Weight

Adult (brass): 209 g

Child (brass): 183 g

Shipping data

Gross weight, packaging & device

Adult (brass): 419 g

Child (brass): 393 g

Box (L x W x H): 377 mm x 196 mm x 43 mm Volume packaging/device: 3177 cm³

Guarantee

5 years

Manufacturer

ERKA. Kallmeyer Medizintechnik GmbH & Co. KG

Im Farchet 15 D 83646 Bad Tölz

Phone: +49 (0)8041 8009-0 Fax: +49 (0)8041 8009-39 E-mail: info@erka.org

Country of origin

Germany