



# Thermal patient cover

TD-57-I.2.c-1.1

edition I

date: 2021-05-26

amendment of: N/A

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland; T: +48 32 376 07 00-02,06 F: +48 32 370 38 08; www.zarys.com

## CLASSIFICATION

Class I medical device compliant with Medical Devices Regulation (EU) 2017/745 (MDR)

## INTENDED PURPOSE

- disposable medical device intended to cover the patient to prevent perioperative hypothermia, used as an insulator to protect against hypothermia caused by impaired thermoregulatory mechanism and to achieve normothermia in the patient
- single use for one patient minimizes the risk of cross infection
- used in perioperative care, intensive care units, hospital emergency departments and other medical departments where it is necessary to ensure the patient's thermal comfort

## FEATURES

- made of polypropylene non-woven with polyester filling ensuring thermal comfort
- ultrasonic seams over the entire surface
- can be heated up in the heating device
- latex-free
- controlled level of microbiological purity
- single-use
- non-sterile
- shelf life: 5 years



CATALOGUE NUMBER	DIMENSIONS	OUTER MATERIAL	KIND OF INDIVIDUAL PACKAGING	SALES PACKAGING	BULK PACKAGING
BT-TOP-11P	110 cm x 120 cm	polypropylene non-woven	foil	1 pc.	50 x 1 pc.
BT-TOP-12P	110 cm x 210 cm	polypropylene non-woven	foil	1 pc.	25 x 1 pc.
BT-TOP-13P	150 cm x 210 cm	polypropylene non-woven	foil	1 pc.	20 x 1 pc.