



UltraMax O₂

Oxygen Analyzer

The UltraMax O₂ oxygen analyzer is designed to measure oxygen concentration, flow and pressure using ultrasound technology that uses a piezoresistive silicon pressure sensor. All relevant measurement data is displayed on the unit's LCD screen.

TECHNICAL SPECIFICATIONS

Oxygen Measurement Range	20.95 – 96% ± 1.5%
Flow Measurement Range	0 – 10 LPM ± 0.2%
Pressure Measurement Range	3.4 – 344 kPa ± 0.5%
Battery Life (2 x AA)	1,100 hours (16,500 read cycles)
Dimensions	8.6cm x 13cm x 2.5cm
Weight	181 gram

Product Features:

- No oxygen sensor required
- No calibration necessary
- Self-diagnostics
- Silicone protective case
- Accurate oxygen measuring
- Durable, compact design
- Convenient pressure unit measurement (PSI or kPa)
- Large, easy-to-read, liquid crystal display (LCD)
- Shielded, reinforced sample gas inlet port
- Auto-off function (4 minutes)
- Low battery indication
- Easy to clean
- Maintenance-free
- Made in USA



DRIVE MEDICAL GMBH & CO. KG

Leutkircher Str. 44
88316 Isny im Allgäu

DE +49 (0) 7562 9724 0
info@drivedevilbiss.com

www.drivedevilbiss.com

drive
DeVilbiss
HEALTHCARE

V2_20200211_EN_0325