



PILL SPLITTER



Easily divides large or small pills

- 1) Set on firm surface
- 2) Place pill inside “V” so pill touches both sides of “V”
- 3) Lower the lid and press lightly
- 4) Store unused portions in base



CAUTION: *sharp blade, keep away from children*



GIMA

**TAGLIAPILLOLE
PILL SPLITTER
COUPE-PILUL
CORTAPÍLDORAS
CORTA-COMPRIMIDOS**

Κομματιαστής
Χαπιών

قاطع لأقراص الدواء

ATTENZIONE: Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto.

ATTENTION: The operators must carefully read and completely understand the present manual before using the product.

AVIS: Les opérateurs doivent lire et bien comprendre ce manuel avant d'utiliser le produit.

ATENCIÓN: Los operadores tienen que leer y entender completamente este manual antes de utilizar el producto.

ATENÇÃO: Os operadores devem ler e entender completamente este manual antes de usar o produto.

ΠΡΟΣΟΧΗ: Οι χειριστές αυτού του προϊόντος πρέπει να διαβάσουν και να καταλάβουν πλήρως τις οδηγίες του εγχειριδίου πριν από την χρήση του.

الحذر: على العمال قراءة وفهم هذا الدليل بكامله قبل البدء باستعمال المنتج.



M25730-M-Rev.3-12.22

Gima S.p.A.

Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com