

ERKA. CLASSIC

EXPERIENCE THE ESSENTIAL



ERKA. CLASSIC

Black	570.00000
Light grey	570.00045
Navy blue	570.00020
Burgundy	570.00060
Dark green	570.00055

The following variations are also available upon request  
via [made@erka.org](mailto:made@erka.org):

- \_ Custom-printed membrane
- \_ Matte chrome surface
- \_ Personalized laser engraving

The best of the ERKA. tradition, reduced to the essentials.

The geometry of the acoustic chamber – perfected over decades – in combination with the modern ERKA. dual membrane provides optimal auscultation over the entire frequency range. The auscultation is optimized by the high specific gravity of the brass body, as this reliably prevents noise interference.

USER GROUPS:

Cardiologists, emergency physicians, paramedics, family physicians,  
assistant physicians, medical students, specialists in respiratory diseases

- 1\_ The construction of the patented ERKA. dual-membrane allows for an additional resonance chamber which optimises the sound.
- 2\_ The non-chill ring ensures patient comfort and protects the membrane.
- 3\_ Highly polished and chromium-plated for a non-porous surface, the brass chest-piece ensures perfect surface hygiene.
- 4\_ Through the separately integrated spring in the tube, there is no acoustic loss from the chest-piece to the ear.
- 5\_ The ergonomically designed, super soft ear tips create an excellent acoustic seal against external noise interference and maximize wearer comfort and fit.
- 6\_ The standard 15° inclined binurals can be individually adapted to the ear canal characteristics of the user.

REPLACEMENT PARTS



EAR TIPS

EAR TIP SOFT  
000.45004

EAR TIP SUPERSOFT  
000.44904

NON-CHILL RINGS

NON-CHILL RING  
ADULT  
000.53104