

## Shiley® Tracheostomy Tube Holder

**A comfortable,  
convenient alternative  
to twill tape.**

- One size adjusts to fit most patients from pediatric to adults, with neck size up to 18".
- Made of soft, foam-padded cotton for patient comfort.
- Velcro® hooks attach easily to tracheostomy tube neck flange.
- Wide band tapers at hooks to reduce excess secretion contact.
- Soft materials and wide band promote better skin care.
- Latex free and designed for single-patient-use.



## Shiley® Stronghold® Anti-Disconnect Devices

An innovative solution to minimize tracheostomy tube disconnects.

- Easy to apply without removing tracheostomy ties.
- Adjusts to accommodate various adaptors.
- Does not restrict movement of circuit or swivel-type adaptors.
- Fastens above the tracheostomy tube neck flange to remain clean.
- Velcro fasteners allow one-handed removal.
- Comfortable for the patient.



**tyco**

Healthcare

**Nellcor**

4280 Hacienda Drive  
Pleasanton, CA 94588  
Tel 925.463.4000  
Toll Free 1.800.635.5267  
www.nellcor.com

Tyco Healthcare UK LTD.  
154 Fareham Road  
Gosport, UK PO13 0AS  
Tel +44.1329.224000

© 2003 Nellcor Puritan  
Bennett, Inc.  
All rights reserved.

B,p 0866-0403  
TR01998

### Ordering Information

Description	Units / Box	Catalog Number
<b>Shiley Tracheostomy Tube Holder</b> For pediatrics and adults, neck size to 18"	10	TTH
<b>Stronghold Anti-Disconnect Device</b> For <i>Shiley</i> XLT tubes and all <i>Shiley</i> adult tracheostomy tubes with rigid neck flange	20	ADD
<b>Stronghold Anti-Disconnect Device</b> For <i>Shiley</i> adult tracheostomy tubes with soft swivel flange	20	SSDD
<b>Stronghold Anti-Disconnect Device</b> For <i>Shiley</i> pediatric and neonatal tracheostomy tubes	20	PDD

For more information contact your local Nellcor sales representative, call 1-800-NELLCOR or visit [www.nellcor.com](http://www.nellcor.com).