Shiley Phonate

SPEAKING VALVE

Provide easier vocalization for tracheostomy patients.

The Shiley® Phonate® Speaking Valve is designed for use on alert, awake tracheostomy patients who can breathe independently without assisted

mechanical ventilation. The

Phonate Speaking Valve connects to the tracheostomy tube

and directs airflow past the vocal

cords to give patients the ability to speak. It permits easier vocalization without finger occlusion, which may cause hygiene and infection concerns. The airflow through the oral and nasal chambers also decreases secretions and improves swallowing and sense of smell.







The Shiley product family from Nellcor includes a complete line of tracheostomy tubes and accessories.

Valve design minimizes work of breathing.

Lightweight valve minimizes additional pressure on the stoma site.

Connector adapts to any tracheostomy tube with a 15-mm hub.

Convenient, hinged cap allows easy cleaning.

Option with port and cap for supplemental oxygen.



Shiley Phonate

SPEAKING VALVE

Specifications

• Airflow Resistance (valve only) at 50 L/min

• Work of Breathing (with size 10)

Connector

2.0 cm H₂0 0.95 joules/L 15 mm

Ordering Information

Description	Quantity	Catalog Number
Phonate Speaking Valve	1 per box	SSV
Phonate Speaking Valve with O ₂ Port	1 per box	SSV0

Nellcor Offers the Complete Family of Shiley Tracheostomy Tube Products

Trusted for exceptional quality and design, *Shiley* has been the leading brand of disposable tracheostomy tubes for over three decades. The comprehensive line of *Shiley* products includes:

- Dual-Cannula Cuffed/Cuffless Tracheostomy Tubes with disposable inner cannula
- Low Pressure Cuffed Tracheostomy Tubes
- Cuffless Tracheostomy Tubes
- TracheoSoft® XLT Extended-Length Tracheostomy Tubes with disposable inner cannula
- Single-Cannula Tracheostomy Tubes
- Pediatric and Neonatal Tracheostomy Tubes
- Tracheostomy Tube Accessories



tyco

Healthcare

Nellcor

4280 Hacienda Drive Pleasanton, CA 94588 Tel 925.463.4000 Toll Free 1.800.635.5267 www.nellcor.com

Tyco Healthcare UK LTD. 154 Fareham Road Gosport, UK PO13 0AS Tel +44.1329.224000

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Specifications subject to change without notice.

© 2003 Nellcor Puritan Bennett Inc. All rights reserved.

B.b 1322-0403 TR00698