

ERKA Sensitive

**The auscultation standard.**

The combination of the sound chamber geometry and ERKA dual membrane technology permits optimum auscultation of both low and high tones. Second to none for blood pressure measurement and cardiology. The perfection of inner values is consummately complemented by the exterior design of this high-end stethoscope. High-gloss polished with chrome plating that is free of pores, the brass bell offers perfect surface hygiene.



Composition and effect of the patented ERKA Dual Membrane:



The floating membrane prevents sound transmission at deformation. A ring as protection against cold improves patient comfort and protects the membrane.



Two parallel interior channels run inside the preformed tube – available in numerous different colours. This avoids friction noise which occurs when two separate tubes are used.



Thanks to the spring integrated in the tube, acoustic signals are conducted from the bell to the ear with virtually no loss.



The ergonomically shaped, super soft ear buds offer a good fit for enhanced wearer comfort and prevent acoustic loss. Naturally, the soft ear buds shown are also included as an alternative.

**ERKA Stethoscope with ERKA Dual Membrane**

	Material · Surface	black	blue	burgundy	grey	green	red
Sensitive	Brass · bright chrome	525.00000	525.00020	525.00060	525.00005	525.00055	525.00030
Precise	Brass · bright chrome	531.00000	531.00020	531.00060	531.00005	531.00055	531.00030
	Brass · bright chrome	531.01000	531.01020	531.01060	531.01005	531.01055	531.01030
Classic	Brass · bright chrome	570.00000	570.00020	570.00060	570.00005	570.00055	570.00030
	Brass · satin chrome	570.01000	570.01020	570.01060	570.01005	570.01055	570.01030

Special colors as pink  or bright blue  are available by request.

Spare-part overview

	Type	Sensitive	Precise	Classic	Finesse	Finesse <sup>2</sup>	ERKAphon duo	ERKAphon
Membrane	Dual-Membrane	000.40200	000.40200	000.40200	-	-	-	-
	Standard Membrane baby	-	-	-	-	000.40102	-	-
	Standard Membrane child	-	-	-	000.40100	000.40100	000.40100	000.40100
	Standard Membrane adult	-	-	-	000.40000	000.40000	000.40000	000.40000
non-chill ring	Dual-Membrane	000.53104	000.53104	000.53104	-	-	-	-
	Size baby	-	-	-	-	000.53103	-	-
	* optional for ERKAphon Stethoscops	-	-	-	000.53142	000.53142	000.53142 *	000.53142 *
	Size adult	-	-	-	000.53141	000.53141	000.53141 *	000.53141 *
Non-chill ring for funnel		-	000.53304	-	000.53304	-	000.53304	-
Eartips	Soft	-	-	000.45004	-	-	000.45004	000.45004
	* optional for ERKAphon Stethoscops	000.45104	000.45104	-	000.45104	000.45104	-	-
	Super Soft	000.44900	000.44900	000.44900	000.44900	000.44900	000.44900 *	000.44900 *

<b>Device type</b> <b>Variations</b>	<b>Stethoscope with dual-duct tube and integrated headpiece with dual membrane and ring as protection against cold</b> UMDNS code: 13-755      Medical device according to class 1 Stethoscope with dual bell in decorative carton with instructions for use Variations: 6 different colours (tube)
<b>Dimensions</b>	Overall length: 715 mm Tube length: 555 mm Earpiece length/Ø: 150 mm/6 mm
<b>Bell dimensions</b>	Bell Ø: 52 mm Membrane Ø: 44 mm Height: 25 mm
<b>Material</b>	Bell: Brass, polished and chrome-plated      Earpiece: Brass, chrome-plated      Yoke spring: Spring steel, chrome-plated Ear buds: TPE/TPS, soft, colour: grey      Ring as protection against cold: NBR      Membrane: boPET (Hostaphan RN 250), 0.25 mm thick CuBE Moulded tube: Soft PVC, latex-free, resistant against common cleaners and disinfectants
<b>Weight</b>	Device: 217 g
<b>Shipping data</b>	Gross weight, packaging & device: 427 g Box (L x W x H): 377 mm x 196 mm x 43 mm Volume packaging/device: 3,177 cm <sup>3</sup>
<b>Guarantee</b>	8 years
<b>Manufacturer</b>	ERKA. Kallmeyer Medizintechnik GmbH & Co. KG Im Farchet 15      D 83646 Bad Tölz Tel.+49 (0)8041 8009-0      Fax:+49 (0)8041 8009-39      E-Mail: made@erka.org
<b>Country of origin</b>	Germany