



ELECTRODE PLACEMENT

T4

LONGISSIMUS DORSI, ILIOCOSTAL AND LUMBAR TRACT




Working position



Working position


T5

ABDOMINIS RECTUS AND EXTERNAL OBLIQUE



T6

ABDOMINIS RECTUS



Working position



Working position



I1

GASTROCNEMIUS



I2

FEMORAL BICEPS, SEMIMEMBRANEUS, SEMITENDINEUS



Working position



Working position

I3

Tensor Fascia Lata

