



Batteria / Battery



Motore / Motor



Ruota / Wheel

**T5011**

36V 8A

350W

16"

- Grande comfort di marcia grazie ad un'apposita progettazione
- Comandi ottici e sonori direttamente dal manubrio
- Cruise control e retromarcia di serie
- Soft start
- Settaggio personalizzabile della macchina
- Bracci di aggancio richiudibili
- Attacchi in lega leggera anodizzata con morsetto universale anti-graffio adatto ad ogni tipo di tubolare della carrozzina
- Inserimento rapido ed intuitivo dei bracci di aggancio con battuta di fine corsa pre-registrabile e sicurezza elettronica di serie per l'arresto macchina in caso di emergenza
- Sicurezza meccanica opzionale anti-sgancio accidentale
- Cover in ABS personalizzabile

- Great ride comfort thanks to a special design
- Optical and sound controls directly from the handlebars
- Cruise control and reverse gear as standard
- Soft start
- Customizable machine setting
- Closable hooking arms
- Anodized light alloy attachments with universal clamp anti-scratch suitable for any type of wheelchair tubular
- Quick and intuitive insertion of the hooking arms with pre-adjustable end run stop and electronic safety standard for machine stop in case of emergency
- Optional mechanical safety against accidental release
- Customizable ABS cover

 Cod. **T5011** 

COD.ISO 12.24.09.015


**Caratteristiche tecniche: / Technical features:**

Struttura / Structure	acciaio verniciato / painted steel
Motorizzazione / Motorization	350 W
Pneumatici / Tires	a scelta in base all'utilizzo / selectable tires based on use
Freni / Brakes	doppio meccanico / mechanical double
Manubri / Handlebars	a scelta tra 4 modelli / to be chosen among 4 models
Display / Display	multifunzione retroilluminato / backlit multi-function
Velocità massima / Maximum speed	programmabile / programmable
Cover / Cover	personalizzabile (opzionale) / customizable cover (optional)



# Caratteristiche tecniche

## Technical features



Dettaglio faro  
Light detail



Dettaglio display  
Display detail



Dettaglio disco  
Disk brake detail



Dettaglio morsetto  
Clamp detail

4 tipologie di manubri  
4 handlebar types

HM1



HM2



HM3



HM4



Configurazioni Configurations	Light T5001	Light T5002	T5011	T5012	T5021	T5022	Cloud T5031
Manubrio 1 / Handlebar 1	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Manubrio 2 / Handlebar 2	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Manubrio 3 / Handlebar 3	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Manubrio 4 / Handlebar 4	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Pneumatico 16x2,05 / Tyre 16x2,5 slike	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Pneumatico 16x2,5 slike / Tyre 16x2,05	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
Pneumatico 16x2,5 tassell / Tyre 16x2,5 tessell	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Acceleratore dx a manopola / Right hand throttle accelerator	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Acceleratore dx a pollice / Right throttle with thumb	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Acceleratore sx a pollice / Left hand throttle accelerator	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Acceleratore sx a manopola / Left throttle with thumb	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Maniglia serraggio large / Large clamping handle	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Cover nero standard / Black standard cover	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
Cover personalizzata / Customized cover	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

● di serie /  
standard  
○ optional /  
optional  
— non disponibile /  
not available

Accessori (a richiesta) Accessories (on demand)	Light T5001	Light T5002	T5011	T5012	T5021	T5022	Cloud T5031
Specchio retrov. Sx / Left rear-view mirror	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Specchio retrov. Dx / Right rear-view mirror	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Catadiottri posteriori / Rear retro-reflectors	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Cavalletto / Easel	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Sicurezza mec. / Mechanical security	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Specchio retrov. compact Sx. / Compact left rearview	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Specchio retrov. compact Dx. / Compact right rearview	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>