

Agitadores magnéticos con calefacción | **Magnetic stirrers with heating**

Marca|Trademark: Nahita

Pequeños agitadores magnéticos con calefacción de regulación analógica e ideales para aquellos lugares en los que se disponga de poco espacio debido a su reducido tamaño. Permiten agitar todo tipo de soluciones en las que se pueda introducir un imán de agitación, hasta un volumen máximo de 2 ó 5 litros dependiendo del modelo. Su estructura metálica recubierta de resina epoxi es resistente a las agresiones químicas mientras que la placa calefactora está fabricada en aluminio anodizado resistente a la corrosión. Presentan un interruptor lateral de encendido/apagado y el panel de mandos frontal incluye un mando para la regulación de la velocidad, un mando para la regulación de la temperatura y dos pilotos de señalización de funcionamiento y de encendido de la placa calefactora.

Los agitadores con calefacción vienen además equipados con una práctica varilla para poder sujetar termómetros o una sonda de pH y se suministran también con un imán revestido de teflón de 9 mm ø y 28 mm de largo.

These small analogical magnetic stirrers with heating are ideal for those places with little space due to their compact design. For mixing up to 2 or 5 L (depending on the model) of all types of solutions in which introduce a stirring magnet is possible. Their metallic structure is covered with epoxy assuring a high resistance to chemicals while the heating plate is made of corrosion resistant anodized aluminium. They present a lateral on/off switch and the front panel includes a speed regulation knob, a temperature regulation knob and two pilot (power and heating) lights.

These magnetic stirrers with heating are also equipped with a useful rod to hold a termometer or a pH probe and are also supplied with a 9x28 mm teflon covered magnet.



Referencia / Code: 50690010



Referencia / Code: 50690050

Referencia / Code	50690010	50690050
Rango velocidad / Speed range	100-2000 r.p.m.	100-2000 r.p.m.
Temperatura máxima (placa) / Max. temperature (plate)	180° C ± 10%	350° C ± 10%
Volumen agitable / Stirrable volume	50-2000 mL	50-5000 mL
Dimensiones plato / Pan dimensions	120x120 mm	170x170 mm
Alimentación / Power	220 V / 50 Hz	220 V / 50 Hz
Dimensiones (LxAxH) / Dimensions (LxWxH)	200x120x90 mm	230x180x120 mm
Peso / Weight	2 Kg	3 Kg

\*Nota: el fabricante se reserva el derecho de modificar las características y fotografías de los productos.

\*Note: Manufacturer reserves the right to modify features and photographs of the goods.