

Medical Anti-fog Disposable face Shield (Isolation face Mask)



protective
mask

Blocking
virus

ented foam headband for increased
rflow and comfort. Anti-fog, anti-glare
ylar polyester film shield. Not made
ith natural rubber latex.



Specific Foam band can absorb sweat and provide enough room for eye glasses or safety goggles
Foam design for increase air flow and comfort.

Excellent protection against potential contamination from blood pathogens, body fluids, or harmful chemical splash

Double side Anti-fog, anti-static, No-glare lens;
Lightweight and comfortable to wear , quick and easy to don.

Our various face shields offer different types of protection and the same anti-fog, anti-glare performance you've come to expect. Bacterial Filtration Efficiency (BFE) is the effectiveness of a material to filter bacteria of a specified particle size, $3.0 \mu\text{m}$ (micron) $\pm \mu\text{m}$. This expressed as a percentage of a quantity that does not pass through the material Particle Filtration Efficiency (PFE) is the effectiveness of a material to filter aerosol particles in the size range $0.1 \mu\text{m}$ to $5.0 \mu\text{m}$.

LENGTH	22 cm
WIDTH	33 cm
SINGLE USE	YES



WUXI EMSRUN TECHNOLOGY CO.,LTD
No.100 Fengbin Road, Wuxi city,
Jiangsu Province, 214000, P.R. Of China
Tel:+86(510)82418838 Fax:+86(510)824199
Email: info@emsr.com
Web: www.emsr.com