


DATI TECNICI / technical sheet	BENDA EMOCONTROL
	compressive bandage

2007/47/CE

93/42/CE

MATERIALI / components

A	TAMPONE ASSORBENTE WP253005(T) <i>absorbent pad</i>	PE bassa densità ; TNT 100% viscosa <i>perforated LDPE ; wnw 100% viscose</i>
B	COMPRESSE DI GARZA 10 X 10 CM <i>cotton gauze compresses</i>	100% cotone <i>100% cotton</i>
C	BENDA ELASTICA m 4,5 x 10 cm <i>elastic bandage</i>	65% cotone 35% poliammide <i>65% cotton 35% polyamide</i>
D	FERMABENDA (2) <i>clamps (2)</i>	gomma e poliester ; alluminio <i>rubber and polyester ; aluminium</i>
E	FOGLIO ISTRUZIONI <i>instruction sheet</i>	carta patinata <i>paper</i>



dim. 120 mm x 55mm

cod. BEN239

PVS S.p.A. - Cassina de Pecchi - Italy	20/02/10
---	-----------------