



The family **Respireo** Primo







We are all different. The **Respireo Primo** series is a range of adjustable masks, specially designed to give a personalized fit to your patients and adapt to each facial shape.

References

Respireo
Primo P



<p>KM219900 > RESPIREO PRIMO P VENTED KM220000 > RESPIREO PRIMO P NON VENTED 3 sizes of cushions S-M-L are included inside the packaging</p>	
<p>KM220700 > CUSHION RESPIREO PRIMO P - Size S (x10) KM220800 > CUSHION RESPIREO PRIMO P - Size M (x10) KM220900 > CUSHION RESPIREO PRIMO P - Size L (x10)</p>	

Respireo
Primo F

<p>KM217900 > RESPIREO PRIMO F - Size S KM218000 > RESPIREO PRIMO F - Size M KM218100 > RESPIREO PRIMO F - Size L</p>	
<p>KM218500 > CUSHION RESPIREO PRIMO F - Size S (x10) KM218600 > CUSHION RESPIREO PRIMO F - Size M (x10) KM218700 > CUSHION RESPIREO PRIMO F - Size L (x10)</p>	
<p>KM220500 > SINGLE PIECE SILICON PAD (x10)</p>	
<p>KM224000 > HEADGEAR RESPIREO Size Standard KM224100 > HEADGEAR RESPIREO Size Large</p>	
<p>KM218200 > RESPIREO PRIMO F Non Vented - Size S KM218300 > RESPIREO PRIMO F Non Vented - Size M KM218400 > RESPIREO PRIMO F Non Vented - Size L</p>	
<p>KM220400 > SIZING GAUGE RESPIREO PRIMO F (x20)</p>	

For other references, please contact your sales manager.

Respireo
Primo B

<p>KM209500 > RESPIREO PRIMO B - Size S KM209400 > RESPIREO PRIMO B - Size M KM209600 > RESPIREO PRIMO B - Size L</p>	
<p>KM210800 > HEADGEAR RESPIREO PRIMO B - Size S KM209700 > HEADGEAR RESPIREO PRIMO B - Size M KM210900 > HEADGEAR RESPIREO PRIMO B - Size L</p> <p><i>Vented and Non-Vented bend are both included inside the packaging</i></p>	

Respireo
Primo N

<p>KM223200 > RESPIREO PRIMO N INT - Size S KM223300 > RESPIREO PRIMO N INT - Size M KM223400 > RESPIREO PRIMO N INT - Size L KM223500 > RESPIREO PRIMO N INT - Size XL</p>	
<p>KM217700 > NON VENTED BLUE BEND (x10)</p>	
<p>KM221900 > DIFFUSE LEAKAGE (x10)</p>	
<p>KM219100 > SILICON PAD (x20)</p>	
<p>KM220500 > SINGLE PIECE SILICON PAD (x10)</p>	
<p>KM224000 > HEADGEAR RESPIREO - Size Standard KM224100 > HEADGEAR RESPIREO - Size Large</p>	
<p>KM223600 > CUSHIONS RESPIREO PRIMO N Size S (x10) KM223700 > CUSHIONS RESPIREO PRIMO N Size M (x10) KM223800 > CUSHIONS RESPIREO PRIMO N Size L (x10) KM221800 > CUSHIONS RESPIREO PRIMO N Size XL (x10)</p>	
<p>KM219800 > SIZING GAUGE RESPIREO PRIMO N (x20)</p>	

Contact :

Air Liquide Medical Systems
Parc de Haute Technologie
6, rue Georges Besse
92182 Antony cedex - France
Tél : +33 (0)1.40.96.66.00
Fax : +33(0)1.40.96.67.00

Customer service

0820 146 359
0,12 € TTC / min
depuis un poste fixe



Air Liquide Healthcare is a world leader in medical gases, home healthcare, hygiene products and healthcare specialty ingredients. It aims to provide customers in the continuum of care from hospital to home with medical products, specialty ingredients and services that contribute to protecting vulnerable lives.

Please read user manual before use. Medical devices class IIIa CE 0051. Manufactured by Air Liquide Medical Systems Spa. Distributed by Air Liquide Medical Systems SA.



Respireo Primo

Simply your well-being



Document non contractuel - produit modifiable à tout moment / YN091511 - V2.2 - 02/05/2014 - Ref:fablo Fischer / Crédit photo : Le Square, Thinkstock.

www.airliquidemedicalsystems.com

Respireo Primo **N**


First-line mask for obstructive sleep apnea



Simple and comfortable, Respireo Primo N is an **adjustable nasal mask** that will successfully assist your patients equipped with CPAP (Continuous Positive Airway Pressure) or NIV (Non Invasive Ventilation) ventilators.

Respireo Primo **F**

From the start of non invasive ventilation



Simple and convenient, Respireo Primo F is a **full face mask with adjustment dial** and superior ease of use.

Respireo Primo **B**

To support situations such as nasal bridge injury, nasal congestion or nasal septal deviation



Simple and flexible, Respireo Primo B is **the first oral mask with no intra-oral mouthpiece** for ultimate patient satisfaction.

Respireo Primo **P**

Recommended for patients with injuries to the bridge of their nose, who suffer from claustrophobia and who may have sensitive eye areas.



Simple and ergonomic, **Respireo Primo P is an adjustable nasal pillow mask** whose silence and ease of use will quickly win your patient's compliance.