


DATI TECNICI / technical sheet	BARELLE PIEGHEVOLI
	<i>foldable stretchers</i>

93/42/CE 2007/47/CE

MATERIALI / components

1	STRUTTURA ESTERNA <i>external frame</i>	ALLUMINIO ANODIZZATO <i>anodised aluminium</i>
2	SNODI <i>joint</i>	FERRO ZINCATO <i>zinc-plated iron</i>
3	TELO PER BARELLA <i>stretcher sheet</i>	PVC <i>pvc</i>



BAR016 pieghevole in 2 parti
 BAR021 pieghevole in 4 parti
 BAR022 pieghevole

foldable in 2 parts
foldable in 4 parts
foldable