

<b>DATI TECNICI / <i>technical sheet</i></b>	<b>COPERTA DI EMERGENZA</b>
	<b><i>emergency blanket</i></b>

## SIKUR TREK

### MATERIALI / *components*

1	BORSA <i>bag</i>	POLIETILENE <i>polyethylene</i>
2	COPERTA LATO SUPERIORE <i>upper side blanket</i>	100% TESSUTO ACRILICO DRALON <i>100 Dralon acrylic tissue</i>
3	COPERTA LATO INFERIORE PLASTIFICATO <i>plastic-coated down blaket side</i>	100% POLIETILENE <i>100% metallized polyethylene</i>



dimensioni /size : 150 x 200 cm

cod : COP210