

<b>DATI TECNICI / technical sheet</b>	<b>SET 4 STECCHE PNEUMATICHE</b>
	<i>inflatable 4 splints set</i>

**QUICK SPLINT**  
 custodia/bag dim : 290 x 380 mm

93/42/CE  
 2007/47/CE

**MATERIALI / components**

1	CUSTODIA PER STECCHE GONFIABILI <i>inflatable splints bag</i>	NYLON <i>nylon</i>
2	STECCHIE GONFIABILI <i>inflatable splints</i>	PVC GONFIABILE <i>inflatable pvc</i>
3	POMPA PER GONFIAGGIO STECCHE <i>inflating pump</i>	POLYPROPILENE <i>polypropilene</i>



**composizione set**  
 pompa di gonfiaggio  
 stecca gamba intera  
 stecca piede/caviglia  
 stecca braccio intero  
 stecca mano/polso

**set list**  
*inflating pump*  
*full leg splint*  
*foot/ankle splint*  
*full arm splint*  
*hand/wrist splint*

2