

**ERKA. SWITCH 2.0
SIMPLEX Ø56 MM**

MAXIMUM INNOVATION
FOR 125 YEARS.



52

ERKA. SWITCH 2.0 SIMPLEX 56 mm with one GREEN CUFF Superb rapid cuff grey

Size 3	294.46894
Size 4	294.40894
Size 5	294.44894

The following variations are also available upon request
via made@erka.org:

— Custom printed cuffs and scales,
e.g. with your company name and company logo

Our aim – Economic efficiency and environmental responsibility.

Our implementation – The new ERKA. with lightning fast, intuitive operation and the ability for users to change the measuring vessel. SWITCH 2.0.Simplex.

Our result - The ERKA. Switch 2.0 Simplex, because

his measuring instrument saves the medical technician time; as part of the metrological control (MTC), the old measuring vessel can simply be replaced by a brand-new, precalibrated measuring vessel.

Saves the user time, because the device is ready for patient use once more after a few seconds.

Saves the administration time, because complex processes in procurement, logistics and storage are significantly simplified.

Conserves resources, because regardless of whether MTC or replacement is involved – only the measuring vessel needs to be changed.

- 1_ The use of a special tempered copper-beryllium membrane ensures durability along with the consistently high precision of the measuring device.
- 2_ The device is shock resistant in accordance with the international standards ANSI/AAMI SP 10.
- 3_ The design of the SWITCH 2.0 is ergonomically well-balanced and the materials used can be cleaned easily and thoroughly.
- 4_ The ERKA. selection of GREEN CUFF cuffs is comprehensive. They are skin-friendly, biocompatible and represent optimal wearer comfort.

53

REPLACEMENT PARTS AND CUFFS

BALLS/MEASURING VESSELS



PUMP BALL grey
000.20900



PUMP BALL grey
with reducing valve grey
000.21100



MEASURING VESSEL SWITCH 2.0 SIMPLEX
with clip/hook
000.30202



CUFFS
Double tube-Cuffs

ERKA. GREEN CUFF
Superb D-Ring cuff
grey

Size 1 infant	10-15 cm	
Size 2 small child	14-21.5 cm	
Size 3 adult/small	20.5-28 cm	025.46894
Size 4 adult	27-35 cm	025.40894
Size 5 adult/large	34-43 cm	025.44894
Size 6 thigh	42-54 cm	