

# COMEN

Share With The World.

In collaborazione con:



## C30 Ambulance Transport Patient Monitor

New generation of specialized ambulance transport patient monitor

# Applications



1 First-aid; on site/ambulance rescue



2 Co-work with C90 for continuous monitoring



3 Patient transferring



4 Bedside monitoring

## Ward round

It's convenient operation for medical staff carrying C30 to measure patients.  
Medical staff can store patient's measure information through C90A.

# C30 Ambulance Transport Patient Monitor

## Features



built-in  
Lithium battery



4.3"  
touch screen



3-level  
alarm



data  
saving and review



light weight  
1.2kg



Linux  
operating system



1.2m drop  
test proof

- Built-in 4400mAh lithium battery in specialized first aid monitoring package.
- With double lithium batteries: Built-in 2,600mAh lithium battery in C30 and built-in 4,400mAh lithium battery in first aid monitoring package. over 8 working hours, support vehicle power.
- 4.3-inch, 16:9 color LCD display, touch screen operation
- Three-level acoustic and optic intelligent alarm and intelligent alarm lamp switching for red, yellow and blue color; Alarm volume also can be adjusted.
- Trend storage & review for 48 hours.
- Weight: 1.2kg; light and Portable.
- With Linux operation system
- 1.2m drop test proof



- Waterproof, shock & drop test proof
- The specialized waterproof and anti-shock C30 first aid bag is portable and easy carrying with enough space for first-aid kits.
- Multi-carry options: carry on shoulder, hung on the stretcher, hung on the ambulance.
- To place first aid articles in left storage bag and all instrument fittings in right one.



- Automatic data exchange with C90
- C30 is multi-parameter module of modular monitor C90



- Multi operation interfaces.
- C30 inherits the concise operation style of COMEN monitor with 9 common shortcut keys, suitable for emergency treatment and first aid.
- Short trend display for 1 hour
- 100 ARR alarm record review.
- II channel ECG waveform review for 5 minutes.

# C30 Ambulance Transport Patient Monitor

## New generation of specialized ambulance transport patient monitor

**COMEN C30** is a new generation of specialized ambulance&transport monitor. It is designed according rescuing site and ambulance transportation monitoring requirements for practicability, convenience and reliability of emergency treatment, first aid and transportation. It adopts PC/ABS high-strength project plastics, high-reliability Linux operation system and high-brightness scratching-prevention LCD. It is equipped with specialized first-aid monitor bag with fitting management pack and first aid device collection pack able for multi-carry . The entire machine is portable, solid, reliable, stable, waterproof, fireproof and anti-falling. It also can be displayed clearly under strong sunshine in the field. It also can be used under various kinds of severe environment to meet demand of medical agents such as troops and hospital, etc during emergency treatment, first aid and transportation in the field.

### Standard Configuration:

ECG, HR, NIBP, SpO<sub>2</sub>,RESP, Dual TEMP, PR, Masimo SpO<sub>2</sub>

### Optional Configuration:

EtCO<sub>2</sub>, Nellcor SpO<sub>2</sub>

First aid monitoring kit





Numero Verde  
**800.777.534**

