


# HEIGHT ROD



## Folding height rod 5003



- > Ultrasound measuring technology, temperature-compensated
- > Tilt sensor to correct if measuring instrument is out of horizontal and error message
- > Small space requirement due to fold function
- > Powered by two 1.5V AA batteries
- > Plastic housing
- > LCD, digit height: 9mm
- > Acoustic acknowledgement after measurement is completed
- > Unit switchover to "inches"

| Order Number | Measuring Range | Divisions | Dimensions | Weight |
|--------------|-----------------|-----------|------------|--------|
|              | cm              | cm        | mm         | g      |